

### 3.6 Service Mode for the Zent2 Version Centrifuge Electronics

To activate the service mode:

Press Cursor Up key

Press Enter three times, holding down the key the third time until the Service Mode appears

The time display = 'Ser'

The program and speed fields display the following parameters:

program	speed	description
0	ver	software version and centrifuge type
1	sen	sensor status
2	imb	N/A
3	error	error table
4	eep—	erase EEPROM or Door Lockout Timer
5	imbs	N/A
6	imbf	N/A
7	led888	LED test
8	zeit	accumulated run time
9	zykl	cycle mode

Select the program using the ↑ and ↓ keys.

To exit any service program press STOP.

Zent2 hardware overspeed jumper settings:

Centrifuge Model	Overspeed (Y Mode)	Overspeed (Δ Mode)
Spinchron DLX	4424rpm (equivalent to 7413 rpm for 10ppr tach) S5, S3, S1	N/A

## 0 Software Version and Centrifuge Type

Press ENTER to activate.

The time and speed fields display the following parameters:

time = Software Version

speed = Centrifuge Type

## 1 Sensor Status

Press ENTER to activate.

The program, time and speed fields display the following parameters:

ACC/DEC =

digit 1, i.e. left = 1 inverter bus charged, 0 inverter bus discharged

digit 2, = 1 inverter enabled, 0 inverter not enabled

TIME =

digit 1, i.e. left = 1 door closed, 0 door not closed

digit 2 = n/a

digit 3 = 1 chamber overtemperature (O/T), 0 no chamber O/T

SPEED =

digit 1, i.e. left = 1 no inverter O/T, 0 inverter O/T

digit 2 = 1 imbalance, 0 no imbalance

Note: for the DLX this signal is inverted

digit 3 = 1/0 rotor ID sensor output

digit 4 = 1/0 tach sensor output

digit 5 = 1 rotor turning, 0 rotor stopped at least 2 seconds

## 2 Imbalance Test N/A

## 3 Error Table

Press ENTER to activate.

The 14 most recent errors are stored in the error table. Pressing the ↑ and ↓ keys scrolls through the list of errors. Pressing RPM deletes the error from the list.

The time and speed fields display the following parameters:

time = sequence number

speed = error number

#### 4 Erase EEPROM or Door Lockout Timer

The speed field displays 'ERASE'. Two options are available by pressing the following keys:

ENTER = deletes the error table; exits Service Mode

not effected: elapsed run time, cycle counter,

RPM = deletes errors and overrides Door Lockout Timer

not effected: elapsed run time, cycle counter; exits Service Mode

The accumulated run time and cycle counter values can be deleted in programs 8 and 9, respectively.

#### 5 Static Imbalance Mode N/A

#### 6 Dynamic Imbalance Mode N/A

#### 7 LED test

Press ENTER to activate.

Every discrete LED and seven segment LED digit is turned ON.

#### 8 Accumulated Run Time

Press ENTER to activate.

The speed and time fields display the following parameters:

speed =      time =  
    hhh    &    hh-mm

A maximum of 99999 hours and 59 minutes can be displayed, about 11 years.

Pressing RPM clears the stored elapsed time in the EEPROM.

Pressing STOP or ENTER exits the program.

#### 9 Cycle Mode

Press ENTER to activate.

The centrifuge can be programmed to cycle automatically. The time and speed fields display the following parameters:

time = wait time (minutes)

speed = cycle counter

The following keys are used to program the cycle mode parameters:

↑ and ↓ = wait time (1 minute increments)

RPM = clear the cycle counter

START = begin cycling and start the centrifuge

FSTOP = end cycling and stop the centrifuge

STOP or ENTER = exit the program

### 3.7 Error Code listing for the Zent2 Version Centrifuge Electronics

No.	Type	Description	Service Procedure
		System CPU:	
1	fail	CPU test after power up failure	replace EPROM, replace Control Board
2	stack	stack overflow	same as error #1
3	CPU RAM	CPU RAM failure	same as error #1
4	RAM external	RAM failure	same as error #1
5	EPROM	checksum invalid	same as error #1
6	comm	communication fault	same as error #1
7	comm	Keypad to Display communication error	connection, replace Keypad and/or Display Board
8	watchdog	frequent tripping	+5volt power supply on Control Board, replace EPROM, replace Control Board
9	status	disallowed SW traps	replace EPROM, replace Control Board
		Motor CPU:	
11	fail	CPU test after power up failure	replace EPROM, replace Control Board
12	stack	stack overflow	same as error #11
13	CPU RAM	CPU RAM failure	same as error #11
14	RAM external	RAM failure	same as error #11
15	EPROM	checksum invalid	same as error #11
16	comm	communication fault	same as error #11
17	reset	Motor CPU reset	same as error #11
18	watchdog	frequent tripping	5 volt power supply on Control board, replace EPROM, replace Control Board
19	status	disallowed SW traps	replace EPROM, replace Control Board
		Tach, System CPU:	
20	unstable	interrupt faulty	tach sensor to Control Board connection tach sensor cable grounding, tach sensor / magnet ring gap, replace tach sensor
21	overspeed	RPM exceeds limit	same as error #20
22	unstable	noisy or jittery	same as error #20
23	compare	inconsistent reading between motor and system CPUs	same as error #20
24	control	final speed not reached	AC input voltage, how easily the rotor turns, easily the rotor turns, replace Power board, replace drive
25	drive	inverter problem	AC input voltage, replace Power board

No.	Type	Description	Service Procedure
27	PWM check	check PWM signal or feedback signal absent	same as error #26
28	no tach	start without tach	same as error #20 lid lockout requires power to be left on for 12-13 minutes, or can be cleared with service mode #4, replace Power board, replace latch
29	no tach	lose tach during run	same as error #28
30	compare	inconsistent overspeed limit between motor and system CPU's	same as error #26
		Tach, Motor CPU:	
31	overspeed	RPM exceeds limit	same as error #20
32	unstable	noisy or jittery	same as error #20
33	compare	inconsistent tach reading between motor and system CPU's	same as error #20
34	control	final speed not reached	same as error #24
35	drive	inverter problem	same as error #25
36	monoflop	no indication of rotor turning even though tach O.K.	replace EPROM, replace Control board
37	PWM	check PWM signal or feedback signal absent	same as error #26
38	no rotor	start without rotor	same as error #20
39	no tach	no tach indication	same as error #20
40	overspeed	hardware overspeed test failure	overspeed jumper settings, replace Control Board
41	overspeed	hardware overspeed detection, at max rotor speed	same as error #40, Check Tach cable routing
42	imbalance	hardware imbalance test failure	imbalance sensor and Control board connection, imbalance sensor position, replace imbalance sensor
43	imbalance offset	DC offset not in range	same as error #42
50	status	floating point operation error, overflow	replace EPROM, replace Control board
51	status	floating point operation error, overflow	same as error #50
52	status	floating point operation error, underflow	same as error #50

No.	Type	Description	Service Procedure
		EEPROM:	
69	type	not compatible	EEPROM loose,
70	type	version not compatible with SW version	same as error #69 same as error #69
71	hardware	not responding	same as error #69
72	verify	not responding	same as error #69
73	checksum	program data not valid	same as error #69
74	checksum	free area not valid	same as error #69
75	checksum	imbalance not valid	same as error #69
76	checksum	offset not valid	same as error #69
77	checksum	program variables not valid	same as error #69
		Door:	
80	transistor	transistor failure	replace Power Board,
81	monoflop	rotor turning with door open	turning rotor by hand with the door open, replace door latch, tach sensor
82	test	door lock test signal failure	door latch to Power Board connection, replace Power, replace Control board
83	open failure	door won't open	replace latch, same as error #82.
98	rotor ID	invalid rotor ID	rotor ID sensor to rotor gap, rotor missing magnets, rotor ID sensor and cable
99	rotor ID	invalid rotor ID	wrong rotor ID keyed in by user, same as error #98