

ROTINA 380 BENCHTOP CENTRIFUGE

FOR **CLINICAL & RESEARCH** APPLICATIONS

FEATURES:

- Robust & steady, temperature control system with pre-cool function from -20 to +40 °C (Optional)
- Versa-Rotor™ for multiple applications (no need to exchange rotors)
- Run-Silent™ Operation (≤ 65 dB, as soft as normal conversation)
- Easy-Lift™ / Safety-Lock™ Lid
- Five (5) Year Warranty
- FDA registered centrifuge manufacturer

Hettich

CENTRIFUGES

CENTRIFUGE PACKAGES

NEW

- BLOOD
- CELL CULTURE
- PLATES
- BOTTLES

(See reverse for more info)



↓
Find out more by visiting:
www.hettweb.com/rotina-380

MAX. CAPACITY:
4 x 290 mL


DIMENSIONS (HxWxD) ROTINA 380: 16.4" x 18" x 23.6" | ROTINA 380R: 16.4" x 18" x 29.5"

TUBE TYPES MICROLITER, ROUND BOTTOM, BLOOD, URINE, KOVA, CONICAL, SKIRTED, SCREW TOP, BOTTLES, LARGE CONICAL, PLATES, ZENTRIFORM


ROTINA 380 CENTRIFUGE PACKAGES FOR CLINICAL & RESEARCH APPLICATIONS

(BELOW IS A LIST OF THE **MOST POPULAR** CENTRIFUGE PACKAGES AVAILABLE FOR THIS UNIT. **ADDITIONAL CONFIGURATIONS** ARE AVAILABLE DEPENDING ON YOUR APPLICATION. PLEASE INQUIRE FOR MORE DETAILS)

ROTINA 380 BLOOD TUBE PACKAGE - (4 or 2)

	Tube Qty	Volume (mL)	Size (mm)	RPM	G-FORCE	PACKAGE CONTENTS
	64 (32)	1.6 - 7	13 x 75/ 100	5,000	3,913	(1) Centrifuge (1) 90° Swing-out Rotor (4-place) (4) Buckets (4) Adapters for 13 x 75/100 mm (4) Adapters for 16 x 75/100/125 mm or (2) Adapters for 13 x 75/100 mm (2) Adapters for 16 x 75/100/125 mm
	52 (26)	4 - 10	16 x 75/ 100			
	12 (6)	8	16 x 125			
Pkg No.:				380BLOOD (4 or 2)		
Pkg No. w/Temperature Control:				380RBLOOD (4 or 2)		
Add Bio-containment Lids				(-BC)		


ROTINA 380 CELL CULTURE PACKAGE - (4 or 2)

	Tube Qty	Volume (mL)	Size (mm)	RPM	G-FORCE	PACKAGE CONTENTS
	36 (18)	15 Conical	17 x 120	5,000	4,863	(1) Centrifuge (1) 90° Swing-out Rotor (4-place) (4) Buckets (4) Adapters for 15 mL (4) Adapters for 50 mL or (2) Adapters for 15 mL (2) Adapters for 50 mL
	16 (8)	50 Conical	29 x 115			
	Pkg No.:					
Pkg No. w/Temp. Control:				380RCCELLCULTURE (4 or 2)		
Add Bio-containment Lids				(-BC)		

ROTINA 380 PLATE PACKAGE

	Tube Qty	Volume (mL)	Size (mm)	RPM	G-FORCE	PACKAGE CONTENTS	
	10	Microtiter	--	4,000	2,397	1) Centrifuge (1) 90° Swing-out Rotor (4-place) (2) Buckets	
	Pkg No.:				380PLATE		
	Pkg No. w/Temperature Control:				380RPLATE		

ROTINA 380 BOTTLE PACKAGE

	Tube Qty	Volume (mL)	Size (mm)	RPM	G-FORCE	PACKAGE CONTENTS	
	4	250	61 x 122	5,000	4,863	1) Centrifuge (1) 90° Swing-out Rotor (4-place) (4) Buckets (4) Adapters	
	Pkg No.:				380BOTTLE		
	Pkg No. w/Temperature Control:				380RBOTTLE		



Find out more by visiting:
www.hettweb.com/rotina-380

DISTRIBUTED BY:

Hettich
CENTRIFUGES

HETTICH NORTH AMERICA
Toll Free: 1-866-370-4388
Email: info@hettweb.com

→ Discover more solutions online at: www.hettweb.com/centrifuges

Copyright © 2012 Hettich Instruments, LP. All rights reserved. For the most up to date information on Hettich please visit hettweb.com.
Printed in the USA using soy based inks on FSC certified recycled paper (25% pre-consumer waste & 25% post consumer waste).

PS-ROTINA-380-E RevA 01-2013