

REPLACEMENT PARTS

NOTE: Parts listed are common to both 115V and 230V models, except where indicated otherwise.

Item Part No.		Schematic Code
Motor Control PCB Assy,	88225	
115V Motor Control PCB	88226	
Assy, 230V. Motor, 115V	SPN38218	MTR
Motor, 230V	SPN38871	MTR
Brushes (pair) Line Cord,	SPN106166	
115V	SPN04538	P1
Line Cord, 230V	20865	P1
Thumbscrew, Rotor	38816	
Washer	05222	
Bushing, Rotor Adapter	38869	
Timer, 0-30 Min.	SPN38822	TMR
Knob, Speed Control	SPN20845	
Knob Assy. TIME Control	36491	
POWER Lamp, 115V	SPN09712	PL1
POWER Lamp, 230V	04589	PL1
Resistor, 50ohms, 50W,	38890	R5
115V	38988	R5
Resistor, 200 ohms, 50W,		
230V		
BRAKE Switch, DPDT	38823	S1
Brush Caps	SPN106550	