

8 General arrangement of the components

Item	Designation
1	Lid
2	Housing
3	Leg spring
4	Sealing ring
5	Centrifuge chamber
6	Motor cover
7	Motor
8	Speed sensor
9	Rubber-metal bearing
10	Upper anti-twist device
11	Lower anti-twist device
12	Imbalance switch
13	Lid lock complete
14	Appliance plug (without fuse holder)
15	Fuse holder
16	Fuse
17	Electronics (230V version)
18	Electronics (120V version)
19	Control panel
20	Flat ribbon cable (4-pole)
21	Cover foil
22	Rubber foot

spare part list 01.12.2010

Material: 104.01.00.00 EBA 270, 100-127V, 50-60Hz Alternative: 1

Lev	Item	Componet	Description	Quantity	Un	Old material no.	Suggested stock levels for 100 field installed units
2	30	103.00.10.05	Lower anti twist device	3,000	PC.	E604	
2	40	104.01.16.00	Lid Lock assy complete	1,000	PC.	E3368	3
2	50	104.01.91.00	Main Electronics board	1,000	PC.	E3306	3
2	60	104.00.10.09	Rubber foot	4,000	PC.	E3352	
2	70	205.200.00.06	Upper anti twist device	3,000	PC.	E2648	
2	100	902.20.26	Rubber-metal-bearing M6 (20x30) 45Sh	3,000	PC.	E2111	
2	110	905.10.132	Motor	1,000	PC.	E3303	3
3	9000	905.10.132.02	Speed sensor	1,000	PC.	E3307	3
2	120	905.19.97.01	Connector plug, combination	1,000	PC.	E3329	3
2	130	905.19.97.02	Fuse holder for E3329	1,000	PC.	E3330	
2	140	905.07.16.05	Fuse T3,15AH (5x20)	2,000	PC.	E997	
3	10	902.32.45	Self-adhesive label - Read	1,000	PC.	L037	
2	20	103.10.70.05	Outer sight-glass and gluing ring	1,000	PC.	E1323	
2	9000	VP.06.300.158	Labels "scratched out waste container"	1,000	PC.	E3010	
2	150	104.00.15.00	Imbalance switch kpl. EBA 270	1,000	ST	E3358	3
2	20	102.85.92.01	Control panel - LED	1,000	ST	E1422	3
2	60	905.55.06.56	Ribbon cable 170mm/4-conductor	1,000	ST	E3357	
2	20	103.10.70.05	View window in lid	1,000	ST	E1323	
1	70	104.00.00.05	bowl seal	1,000	ST	?	3

