

Component	Material Description	Old material number	Component quantity
112.01.00.00	ROTINA 35, 110-127V 50-60Hz	1705-01	1
902.19.06.01	Rubber foot	2251	1.000
902.20.21.03	Rubber-metal-bearing M8 (25x40) 45Sh	E1889	3.000
905.10.22	Frequency motor 115V 60Hz	OCM-33G	1.000
D7343 2.5X24A2	Spiral spring pin ISO 8750-2,5x24-A2	E1126	1.000
205.26.50.01	Spacer	E3294	1.000
111.01.21.00.1	Speed sensor on top of the motor	E730	1.000
905.59.03.04	Cable tie 102 mm	E205	2.000
205.17.00.05	Anti-twist-device	E732	6.000
112.00.13.00	Ventilator housing	E1240	1.000
905.11.29.02	Radio suppression filter 115V	E1284	1.000
112.01.91.00	Frequency converter	E1184	1.000
112.01.10.07	Braking resistor cpl.	E1469	1.000
103.10.10.10	Imbalance switch sub-assembly	E2590	1.000
112.00.12.03.1	Imbalance switch with angle	E2531	1.000
905.19.49.01	Connector plug with line filter	E507	1.000
905.15.79.10	Power switch ON/OFF, 2-pole, 9A	E1639	1.000
112.02.10.02	Motor housing	E1395	1.000
112.02.00.01	Hinge	E1316	2.000
103.05.93.00.B	Supply board 115V 60Hz	E1888	1.000
905.55.06.14	Control lead 170/10 pole	E1227	1.000
905.55.03.13	Control lead 230/16 pole	E1311	1.000
111.51.93.00.3	Locking shackle 16-pole	E776	1.000
112.02.20.05.1	Upper part	E1697	1.000
103.12.20.14	Lock	E1477	1.000
103.11.20.01	Lid lock with stroke magnet	E1478	1.000
905.15.71	Micro switch 10A 250V	E1012	1.000
112.00.00.04	Sealing ring	E1242	1.000
112.02.70.00.1	Lid complete	E1699	1.000
112.02.70.01	Shackle	E2170	1.000
103.12.20.12	Roller	E1562	2.000
103.10.70.05	Outer sight-glass and gluing ring	E1323	1.000
103.12.20.13	Handle rail	E1481	1.000
112.02.00.02.1	Bolzen (Montagebaugruppe)		2.000
112.02.00.02	Bolt	E1317	1.000
D0471 8X0.8	Washer	E1132	1.000
D0988 10X16X0.5	Shim ring	E1318	4.000
D0471 8X0.8	Washer	E1132	2.000
902.13.05	Gas spring 130N with bolt	52004	1.000
902.13.05.03.2	U-head M5x10	E2697	1.000
D0471 6X0.7	Washer	E509	1.000
111.51.93.00.3	Locking shackle 16-pole	E776	1.000
112.01.00.03.1	Front panel	E1682-1	1.000
111.20.92.02	Cover foil	E2671	1.000
112.01.92.07.1	E-Prom ROTINA 35 -3.12	E1690	1.000
112.02.23.00	Back wall	E1389	1.000
205.26.00.05	Motor cover	E907	1.000
902.32.47	Self-adhesive label - arrow black -	E2286	1.000

103.05.93.00.A	Supply board 120V 60Hz	E949	1.000
103.05.93.00.B	Supply board 115V 60Hz	E1888	1.000
103.05.93.00.A.9	Supply board reconditioned	E1369	1.000
103.10.00.05	Adjusting knob complete	E1315	1.000
111.00.20.05.2	Flat spring	E727	2.000
103.03.20.02	Lid lock with with stroke	E1385	1.000
902.10.03.02	Traction spring	E2152	1.000
111.00.22.00.1	Lock without stroke magnet	E426	1.000
902.10.03.02	Traction spring	E2152	1.000
112.00.20.03	Push rod	E1229	1.000
112.01.92.00	Control panel complete, 115V 60Hz	E1500	1.000
103.10.92.03	Control panel	E3504	1.000
112.01.92.07	E-Prom ROTINA 35 -1.12	E1501	1.000
112.01.92.00.1	Control panel complete 115V 60Hz	E1692	1.000
111.20.92.01	Control panel	E3502	1.000
111.20.92.02	Cover foil	E2671	1.000
112.01.92.07.1	E-Prom ROTINA 35 -3.12	E1690	1.000
112.01.92.07	E-Prom ROTINA 35 -1.12	E1501	1.000
112.02.00.03.1	Front panel	E1321	1.000
112.02.20.03.1	Upper part	E1391	1.000
112.02.70.00	Lid complete	E1387	1.000
103.10.70.05	Outer sight-glass and gluing ring	E1323	1.000
112.02.71.00.1	Lining plate complete	E1319	1.000
112.02.71.00.2	Lining plate complete	E1688	1.000
905.15.79.04	Power switch ON/OFF	E1008	1.000