



The world leader in serving science

GP3 Spareparts - Ersatzteile

Heraeus: Multifuge X3, Megafuge 40
Sorvall: Legend XT, ST 40

Overview - Übersicht

- [Outer Casing – Aussengehäuse](#) [Parts List](#)
- [Lid Cover – Deckel](#) [Parts List](#)
- [Electronics – Elektronik](#) [Parts List](#)
- [Service Kit](#) [Parts List](#)

Power cord P/N 20190346

Centrifuges Overview – Zentrifugen Übersicht

Heraeus Megafuge 40 / Sorvall ST 40



Heraeus Multifuge X3 / Sorvall Legend XT



Outer Casing Overview – Aussengehäuse Übersicht



Key Foils – Tasterfolien



Pos. 1120

Pos. 1140




Pos. 1110

Pos. 1130

Spare Parts Outer Casing – Ersatzteile Aussengehäuse

Pos	Component	Outer casing	Aussengehäuse	Multi X3	Mega 40	Legend XT	ST 40
1000	20058548	Front Cover	Frontblende	X	X	X	X
1110	20058565	Overlay Multifuge X3	Tasterfolie Multifuge X3	X			
1120	20058566	Overlay Megafuge 40	Tasterfolie Megafuge 40		X		
1130	20058567	Overlay Legend XT	Tasterfolie Legend XT			X	
1140	20058568	Overlay ST 40	Tasterfolie ST 40				X
1200	20058539	Display frame high end	Display Rahmen high end	X		X	
1250	20058538	Display frame routine	Display Rahmen routine		X		X
1300	20150266	Display to Overlay adaptor	Display Folienadapter	X	X	X	X
1400	20150281	Display 16+24 high end	Display 16+24 high end	X		X	
1450	20150282	Display 26x1 routine	Display 26x1 routine		X		X


Lid air cooled – Deckel Luftgekühlt



Spare Parts Lid Cover – Ersatzteile Deckel

Pos	Component	Lid Cover	Deckel
2000	20058531	Lid Cover	Deckel Oberteil
2010	20058525	Lid Hinge	Deckelscharnier
2020	20050543	Window	Schauglas
2030	20058558	Cover Foil	Abdecker Deckelwulst
2040	20056650	Holder	Lager bock
2050	20058480	Hook	Kloben GP
2060	20058553	Insert	Schaumeinleger Deckel
2070	20058524	Insert	Deckelblech

Electronics – Elektronik



Pos. 3000


Pos. 3600

Pos. 3110

Pos. 3120

Pos. 3130


Pos. 3000 umfasst die gesamte Elektronikassette inkl. Pos. 3110, 3120 und 3130, aber nicht 3600!
Pos. 3000 includes all parts 3110, 3120 and 3130, but not 3600!



Pos. 3300

Pos. 3400

Pos. 3200



Pos. 3500


Spare Parts Electronics – Ersatzteile Elektronik

Pos	Component	Electronics	Elektronik	230V	120V	100V	230V US
3000	20150258	HP General Purpose 120V	HP General Purpose 120V		X		
3000	20150259	HP General Purpose 230V	HP General Purpose 230V	X			X
						X	
3110	20210443	Fan	Elektronik Lüfter	X	X	X	X
3120	70904621	Receptacle	Gerätestecker TYP 711-2-PW30		X	X	
3120	20180748	Receptacle	Gerätestecker TYP 722-1-PW30	X			X
3130	20180747	Overcurrent switch	Überstromschutzschalter	X	X	X	X
3200	20210429	Motor GP	Motor General Purpose	X	X	X	X
3300	70902026	Autolock Assy SH65y	Autolock komplett	X	X	X	X
3400	20150270	PCB speed imbalance rotor	LP Drehzahl Unwucht	X	X	X	X
3500	70902692	Lid Latch	Verschluss vollständig	X	X	X	X

NV-RAM

Pos	Component	Version: Air-cooled / Luftgekühlt	Equipment	230V	120V	100V	230V US
3600	70904639	Heraeus Multifuge X3	75004500	X			
	70904639	Sorvall Legend XT	75004505	X			
	70904635	Heraeus Megafuge 40	75004503	X			
	70904635	Sorvall ST 40	75004509	X			
	70904638	Heraeus Multifuge X3	75004501		X		
	70904638	Sorvall Legend XT	75004506		X		
	70904634	Heraeus Megafuge 40	75004504		X		
	70904634	Sorvall ST 40	75004510		X		
	70904650	Sorvall Legend XT	75004507			X	
	70904648	Sorvall ST 40	75004511			X	
	70904652	Sorvall Legend XT	75004508				X
Pos	Component	Version: Refrigerated / Gekühlt	Equipment	230V	120V	100V	230V US
3600	70904641	Heraeus Multifuge X3R	75004515	X			
	70904641	Sorvall Legend XTR	75004520	X			
	70904637	Heraeus Megafuge 40R	75004518	X			
	70904637	Sorvall ST 40R	75004524	X			
	70904640	Heraeus Multifuge X3R	75004516		X		
	70904640	Sorvall Legend XTR	75004521		X		
	70904636	Heraeus Megafuge 40R	75004519		X		
	70904636	Sorvall ST 40R	75004525		X		
	70904651	Sorvall Legend XTR	75004522			X	
	70904649	Sorvall ST 40R	75004526			X	
	70904653	Sorvall Legend XTR	75004523				X

Service Kit



[Pos. 4000](#)

Service Kit

Pos	Component	Mechanics	Mechanik	air-cooled	refrigerated
4000	70904694	Service Kit anti-vibration mounts	Service Kit Gummipuffer	X	X