



6 SPARE PARTS - DIAGRAMS

6.1. GENERAL




Item	B4i Cat. N°	BR4i Cat. N°	Description	Qty
1	89000838	86001434	Power switch	1
2	85220838	85220761	Keyboard	1
3	89000840	89000840	Gas spring	1
4	89000839	89000839	Motor chamber/bowl seal	1
5	89000841	89000943	Lid/bowl seal	1
-	85280460	85280460	Main p.c.b.	1
-	85280514	85280514	Display p.c.b.	1

6.2. MOTOR



Item	B4i Cat. N°	BR4i Cat. N°	Description	Qty
1	89000878	89000878	Motor	1
2	89000843	89000843	Load imbalance sensor	1
3	89000853	89000853	Elastic suspension	1
4	86001483	86001483	Spring washers	3
5	86001484	86001484	Washers	2
6	85240971	85240971	Motor shaft clamp	1
7	86001449	86001449	O-Ring	1
-	89000856	89000856	Tachometer sensor assy	1


6.3. AUTO-LOCK SYSTEM




Item	B4i Cat. N°	BR4i Cat. N°	Description	Qty
1	89000914	89000914	Guide collar	1
2	85230694	85230694	Latch	1
3	86001482	86001482	Spring	1
4	86001481	86001481	Circlip	1
5	89000915	89000915	Latch guide	1
6	85241222	85241222	Rotating equipment unlocking device	1

6.4. LOCK

B4i lid lock




BR4i lid lock




Item	B4i Cat. N°	BR4i Cat. N°	Description	Qty
1	89000842	89000842	Lid microswitch	2
2	89000844	89000844	Lock	1
3	89000845	89000845	Solenoid switch	1
4	89000846	89000846	Locking system spring	1
5	89000847	89000847	Lid solenoid	1
6	89000848	89000901	Locking system assembly	1

6.5. FEEDER GROUP

B4i feeder group




BR4i feeder group




Item	B4i Cat. N°	BR4i Cat. N°	Description	Qty
1	89000900	89000900	Mains connector	1
2	-	85280383	Compressor p.c.b	1
3	89000854	89000854	Braking resistor	2
4	89000837	89000837	Transformer	1
5	89000855	89000882	Feeder group	1

6.6. REFRIGERATION GROUP (BR4i)



Item	B4i Cat. N°	BR4i Cat. N°	Description	Qty
1	-	89000902	Bowl assembly	1
2	-	89000903	Condenser group assembly 230 V - 50 Hz	1
2	-	89000904	Condenser group assembly 120 V - 60 Hz	1
3	-	89000905	Compressor 230 V - 50 Hz	1
3	-	89000906	Compressor 120 V - 60 Hz	1

6.7. CONDENSER GROUP (BR4i)



Item	B4i Cat. N°	BR4i Cat. N°	Description	Qty
1	-	89000907	Condenser	1
2	-	89000908	Screw	1
3	-	89000909	Fan blades	1
4	-	89000910	Motor nuts	1
5	-	89000911	Carter fan blades and motor	1
6	-	89000912	Motor 230 V - 50 Hz	1
6	-	89000913	Motor 120 V - 60 Hz	1