







Biofuge 13 # 75003613/02			120v 60Hz		
Index	Part Number	Description	Index	Part Number	Description
000	70901018				
101	20460133	Screw M3x14	403	20230149	Fuse Holder
102	20022598	Front Panel	404	20230190	Fuse 10a Slow
103	20022611	Key Foil	405	20460100	Screw M3x10
110	20150072	Main PCB	406	20420076	Hex. Nut M4
111	20460137	Screw 3x8	407	20480108	Fan Washer A4,3
112	70901049	NV.Ram	408	20220289	Terminal Strip
113	20230189	Fuse 6.25a slow	409	20460087	Allen Screw M3x16
114	20230189	Fuse 6.25a slow			
116	20022615	Wire netting		Revised 2/22/02	
117	20022616	Wire netting			
201	70901006	Armoured bowl			
202	20049932	Foot			
203	20310294	Securing plug			
204	20460088	Allen screw M5x8			
205	20020706	Housing			
206	20020712	Seal			
207	20020713	Seal			
210	20020743	Lid			
211	20050543	Window			
212	20350211	Collar			
213	20480250	Spring washer			
214	20460074	Screw M6x16			
215	20290419	O-Ring 4.47x1.78			
220	20180581	Mains Switch			
300	20003620	Lid lock assy.			
301	20110297	Solenoid			
302	20180557	Microswitch			
303	20460135	Screw 1.8x10			
304	20460138	Screw 4x10			
305	20460094	Screw M4x8			
306	20020745	Shim			
310	20210290	Motor			
311	20290520	O-Ring 15x2.5			
313	20420070	Nut M6			
314	20480168	Fan Washer			
315	20480301	Emery paper washer			
316	20020732	Antivibration mount			
317	20480168	Fan Washer			
318	20430216	Screw M6x12			
321	20460094	Screw M4x8			
322	20480108	Fan Washer A4,3			
323	20901059	Mains Cable			
324	20220286	Grip 6.5			
325	20480108	Fan Washer A4,3			
326	20460094	Screw M4x8			
327	20540048	Capacitor			
328	20480180	Fan Washer			
401	20170291	Transformer			

Bold = SPL YES

402	20220472	Insulation			
-----	----------	------------	--	--	--

Bold = SPL YES

Biofuge 13 #75003637/02			120v 60Hz		
Index	Part Number	Description	Index	Part Number	Description
000	70901018	Front panel complete	403	20230149	Fuse Holder
101	20460133	Screw M3x14	404	20230190	Fuse 10a Slow
102	20022598	Front Panel	405	20460100	Screw M3x10
103	20022601	Key Foil	406	20420076	Hex. Nut M4
110	20150072	Main PCB	407	20480108	Fan Washer A4,3
111	20460137	Screw 3x8	408	20220289	Terminal Strip
112	70901048	NV.Ram	409	20460087	Allen Screw M3x16
113	20230189	Fuse 6.25a slow			
114	20230189	Fuse 6.25a slow			
201	70901006	Armoured bowl		Revised 4/9//02	
202	20049932	Foot			
203	20310294	Securing plug			
204	20460088	Allen screw M5x8			
205	20020706	Housing			
206	20020712	Seal			
207	20020713	Seal			
210	20020743	Lid			
211	20050543	Window			
212	20350211	Collar			
213	20480250	Spring washer			
214	20460074	Screw M6x16			
215	20290419	O-Ring 4.47x1.78			
220	20180581	Mains Switch			
300	20003620	Lid lock assy.			
301	20110297	Solenoid			
302	20180557	Microswitch			
303	20460135	Screw 1.8x10			
304	20460138	Screw 4x10			
305	20460094	Screw M4x8			
306	20020745	Shim			
310	20210290	Motor			
311	20290520	O-Ring 15x2.5			
313	20420070	Nut M6			
314	20480251	Securing Washer			
315	20480301	Emery paper washer			
316	20020732	Antivibration mount			
317	20480251	Securing Washer			
318	20430216	Screw M6x12			
321	20460094	Screw M4x8			
322	20480108	Fan Washer A4,3			
323	20901059	Mains Cable			
324	20220286	Grip 6.5			
325	20480108	Fan Washer A4,3			
326	20460094	Screw M4x8			
327	20540048	Capacitor			
328	20480180	Fan Washer			
401	20170291	Transformer			
402	20220472	Insulation			

Bold = SPL YES

Medifuge #75003632/03**120V 60Hz.**

Index	Part Number	Description	
000	70901013	Front Panel Complete	
001	70901064	Service Kit	
002	20053358	Rotor Nut	
101	20460133	Screw 3x14	
102	20022598	Front Panel	
103	20022900	Key Foil	
110	20150178	Main Board	20150062 original pcb. obsolete
111	20460137	Screw 3x8	
112	70901066	NV-Ram	
113	20230189	Fuse 6.25 slow	
114	20230189	Fuse 6.25 slow	
201	70003705	Bowl, complete	
202	20049932	Foot	
203	20310248	Securing Plug	
204	20460088	Allen screw M5x8	
205	20020726	Housing	
206	20020717	Seal	
210	70020803	Lid	
211	20050543	Window	
212	20350211	Collar	
213	20480250	Spring Washer	
214	20460074	Screw M6x16	
215	20290419	O-Ring 4.47x1.78	
220	20180581	Mains Switch	
300	70003621	Lid Lock Assy.	
301	20110298	Solenoid	
302	20180557	Microswitch	
303	20460135	Screw 1.8x10	
304	20460138	Screw 4x10	
305	20460094	Screw M4x8	
306	20020745	Shim	
310	20210276	Motor	Motor clinifuge: 20210275
311	20020738	Motor Cover	
313	20420070	Nut M6	
314	20480251	Securing Washer	
315	20480301	Emery paper washer	
316	20020744	Anti-vibration mount	
317	20480251	Securing Washer	
318	20430216	Screw M6x12	
321	20460094	Screw M4x8	
322	20480108	Fan washer A 4.3	
323	20901059	Mains cable	
324	20220286	Grip 6.5	
325	20480108	Fan washer A 4.3	
326	20460094	Screw M4x8	
327	20540344	Capacitor	
328	20480180	Fan Washer	
	20490376	Motor shaft pin	