



5. SPARE PARTS

5.1. C4i GENERAL SCHEMATIC




Mark	Reference	Description	Qty
1	86000287	Breaker for C4i 230 V	1
1	86001434	Breaker for C4i 120 V	1
2	89002885	Blue wing + keyboard	1
3	26962008	Gas spring	2
4	25439072	Motor / bowl seal	1
5	25437049	Lid / bowl seal	1
6	85241270	Lid hinge	2
7	26566011	Rubber foot	4
8	89002886	Complete lid	1

5.2. CR4i GENERAL SCHEMATIC




Mark	Reference	Description	Qty
1	86001434	Breaker	1
2	89002885	Blue wing + keyboard	1
3	26962008	Gas spring	2
4	25439072	Motor / bowl seal	1
5	25437049	Lid / bowl seal	1
6	85241270	Lid hinge	2
7	26566011	Rubber foot	4
8	89002887	Complete lid	1

5.3. MOTOR




Mark	Reference C4i/CR4i	Description	Qty
	89002888	Motor kit (includes shock absorbers)	1
1	26007026	Shock absorber (1 pce)	3
2	85480878	Imbalance detector	1

5.4. LOCK



Mark	Reference C4i/CR4i	Description	Qty
1	89002889	Complete lock	1
	89000309	Lid solenoid (without switch)	1
2	26487030	Micro switch (without lever)	1
	26461053	Flat lever	1
3	86001406	Safety micro switch	2

5.5. ELECTRONIC RACK




5.5.1. C4i RACK

Mark	Reference 230 V	Reference 120 V	Description	Qty
1	89002891	89002892	Main PCB	1
2	86004585	86004586	Relay	1
3	26387034	26387045	Fuse	1
4	86003242	86003243	Transformer	1
5	86004598	86004598	Diodes bridge	1
	89002895	89002896	Complete electronic rack	1

5.5.2. CR4i RACK


Mark	Reference 230 V	Reference 120 V	Description	Qty
1	89002893	89002894	Main PCB	1
2	86004585	86004586	Relay	1
3	26387034	26387045	Fuse	1
4	86003242	86003243	Transformer	1
5	86004598	86004598	Diodes bridge	1
6	89002899	89002899	Cooling unit PCB	1
	89002897	89002898	Complete electronic rack	1

5.6. FRONT PANEL



Mark	Reference	Description	Qty
1	89002885	Blue wing + keyboard	1
2	89002890	Display PCB	1

5.7. COOLING UNIT (CR4i)



Mark	Reference 230 V	Reference 120 V	Description	Qty
1	86001138	86001142	Compressor	1
2	26213006	26213006	Condensor	1
3	26408017	26408017	Fan	1
4	26960023	26960024	Motor + stand (without fan)	1