

11. PARTS LIST

11.1 Model CL

			IEC PART NO.	
UNIT	DESCRIPTION	QTY	115V	230V
(1)	Base		42985	42985
(2)	Top Cap		42470	42470
(3)	Assembly, Cover and Hinge		43933	-43933
(4)	Motor Armature		42469	43007
(5)	Motor Field		4173	9629
(6)	Brush and Spring (2)		1780	1780
(7)	Brush Holder (2)		9037	9037
(8)	Brush Holder Cap (2)		9036	9036
(9)	Nut		1729	1729
(14)	Rubber Cushion (3)		1719	1719
(15)	Base Screen		36401	36401
(16)	Switch and Transformer		36536	36537
(17)	Nameplate		42813	42814
(18)	Switch Knob, Raytheon #409D2K3		42812	42812
(19)	Cord Set		38805	38806
(20)	Cap Plug, #BP1½		44241	44241
(21)	Rnd Head Phillips Screw, 8-32 UNC - 2A, x 5/16" long (2)		COM'L	COM'L
(22)	Hex Nut, 8-32 UNC - 2B (2)		COM'L	COM'L
(23)	Lock Washer #8 Internal (2)		COM'L	COM'L
(24)	Rnd Head Phillips Screw, 10-32 UNC - 2A, x 2" long (2)		COM'L	COM'L
(25)	Screw, Special (3)		45250	45250
(26)	Square Nut, 10-32 UNC - 2B (3)		COM'L	COM'L
(27)	Rnd Hd. Screw, ¼-20 UNC - 2A, x 3/8" long (3)		COM'L	COM'L
(28)	Lock Washer, ¼" Split (3)		COM'L	COM'L
(29)	Drive Screw, #0 x ¼" long (4)		COM'L	COM'L
(30)	Lock Washer, #10 Split (3)		COM'L	COM'L
(31)	Set Cup Pt. Screw, 10-32 UNC - 2A, x 1/8" long (2)		COM'L	COM'L
(32)	Rivet, Dome Hd. Avex #1601-0613 (2)		COM'L	COM'L
(34)	Rd. Hd. Phillips Screw, #10-32 UNC-2A x 3/8" long		COM'L	COM'L

<u>UNIT</u>	<u>DESCRIPTION</u>	<u>QTY</u>	<u>IEC PART NO.</u>	
			<u>115V</u>	<u>230V</u>
(36)	Strain Relief, Heyco SR-6L-1		COM'L	COM'L
(38)	Guard and Hinge Assembly		4154	4154
(39)	Button		1779	1779
(41)	Knob, #244 DIMCO		38595	38595
(42)	Washer, #8 Shakeproof		COM'L	COM'L
(43)	Pan Hd. Slotted Screw, #8-32 UNC - 2A x .31" long		COM'L	COM'L
(46)	Spring Finger Washer Spec. #S-10		41778	41778
(47)	Bearing-Upper and Lower, Fafnir 200-PP or Equiv.		9939	9939
	Adapter Hub for fastening Hematocrit to Clinical		1116	1116

11.2 Model CH

UNIT	DESCRIPTION	QTY	IEC PART NO.	
			115V	230V
(1)	Guard Bowl Assembly		1399	1399
(7)	Rubber Cushion (3)		1719	1719
(8)	Cord Set		38805	38806
(10)	Special Screw (3)		45250	45250
(11)	Square Nut, Stl. #10-32 UNF - 2B (3)		COM'L	COM'L
(12)	Rd. Hd., Phillips Screw, Stl., #8-32 UNC - 2A x 5/16" lg. (3)		COM'L	COM'L
(13)	Hex Nut, Stl., #8-32-2B (2)		COM'L	COM'L
(14)	Lock Washer, Stl., #8 Internal (2)		COM'L	COM'L
(15)	Switch and Transformer		36536	36537
(16)	Base Screen		36401	36401
(17)	Assembly Base		42985	42985
(20)	Rd. Hd. Phillips Screw, Stl., #10-32 UNF-2A x 3/8" lg.		COM'L	COM'L
(21)	Armature		42469	43007
(22)	Field		4173	9629
(24)	Brush and Spring (2)		1780	1780
(25)	Brush and Holder (2)		9037	9037
(26)	Brush Holder Cap (2)		9036	9036
(27)	Nut Keps, Stl. Shakeproof #511-101800-00 (3)		COM'L	COM'L
(29)	CH Nameplate		42815	42816
(30)	Switch Knob, Raytheon - 409D2K3		42812	42812
(31)	Drive Screw, Steel (4)			
(32)	Set Cup Point Screw, Stl., #10-32 UNF-2A x 1/8" (2)		COM'L	COM'L
(33)	Rnd. Hd. Phillips Screw, Stl., #10-32 UNF - 2A x 2" (2)		COM'L	COM'L
(34)	Lock Washer, Stl., #10 - Split		COM'L	COM'L
(35)	Nut		1729	1729
(36)	Top-Cap		42470	42470

<u>UNIT</u>	<u>DESCRIPTION</u>	<u>QTY</u>	<u>IEC PART NO.</u>	
			<u>115V</u>	<u>230V</u>
(38)	Strain Relief Heyco SR-6L-1		COM'L	COM'L
(40)	Spring Washer Finger Spec. #S-10		41778	41778
(41)	Bearing, Fafnir 200-PP or Equivalent	(4)	COM'L	COM'L
(42)	Caplug #BP-1½		44241	44241

NOTE: When ordering replacement parts, the serial number of the centrifuge is necessary.