


PARTS LIST

See exploded view with reference numbers opposite.

120V CENTRIFUGE

1. **MOTOR & HARDWARE ASSEMBLY:** 42010101
Nut, Flex Loc. No. 8-32 (3)
Seal, Motor
Spring Pin, 3/32 Dia. × 7/16 Lg.
Head Retainer, (Collar Clamp)
Set Screw, Socket (Allen) (2)
Head Screw Assembly
Wiring Instructions
2. **MOTOR BRUSH KIT (for 120 V & 220 V Motors):** 42000102
Brush Assembly (2)
Outer Caps (2)
3. **TIMER WITH KNOB:** 42010103
Knob
Timer
Lockwasher, Int. Tooth
Wiring Instructions
4. **KNOB ONLY, TIMER:** 42010104
5. **SWITCH, BRAKE:** 42010105
Switch & Locknut
Ring Nut
Wiring Instructions
6. **CORDSET ASSEMBLY:** 42054109
Cord Set
Wiring Instructions
7. **COVER ASSEMBLY:** 42010308
8. **LATCH ASSEMBLY** 42010309
9. **SPEED CONTROLLER ASSEMBLY:** 42010108
Wiring Instructions
10. **KNOB ONLY, SPEED CONTROL:** 42010109
11. **LIGHT ASSEMBLY & CONNECTOR:** 42010110
Light with Spring Fastener
Push-On Connector
Wire Black 18 GA. 4" Lg.
Wiring Instructions
Connector (2)
12. **HEAD SCREW ASSEMBLY & RETAINER:** 42010111
Head Screw Assembly
Head Retainer (Collar Clamp)
Set Screw, Socket (Allen)
No. 6-32 N.C. × 1/4 D.P. (2)
13. **RUBBER BUMPER:** 42010112
3 per Package
14. **RUBBER FOOT & HARDWARE:** 42010113
Rubber Foot with Washer (4)
Screw, Phil. Hd.
No. 10-32 × 5/8 Lg. (4)
15. **LID CUSHION:** 42010310
16. **LATCH PLATE:** 42010312



Exploded View, DYNAC Centrifuge
Models 420101 and 420102

220V CENTRIFUGE

17. **MOTOR & HARDWARE ASSEMBLY:** 42010201
Motor
Nut, Flex. Loc. No. 8-32 (3)
Seal, Motor
Spring Pin, 3/32 Dia. × 7/16 Lg.
Scotch Foam, 11" Lg.
Head Retainer (Collar Clamp)
Set Screw, Socket (Allen) (2)
Head Screw Assembly
Wiring Instructions
18. **SPEED CONTROLLER ASSEMBLY:** 42010202
Wiring Instructions

- *19a. Screw, Pan, Head
No. 6-32 × 3/8 Lg.
- *20a. Screw, Flat Head
No. 6-32 × 5/16 Lg.
- *21a. Nut, Hex No. 6-32
Lockwasher Split No. 6

SPECIFICATIONS (both models)

Height: 10 in. (25.4 cm)
Bowl Diameter: 15 in. (38.1 cm)
Overall Width: 18.5 in. (47 cm)
Net Weight (base only): 27 lbs. (12.25 kg)

*Attaching hardware not included. Purchase locally.