



- 1 Lid.....E1666
- 2 Housing upper.....E1664 part
- 3 Cover foil.....E1665
- 4 Control panel.....E1422R
- 5 Lid lock.....20161
- 6 Motor.....OCM-33B
- 7 Rubber-metal-.....E239 bearing
- 8 Power board.....52009
- 9 Appliance plug.....E1661
- 10 Housing lower.....E1663 part
- 11 Rubber foot.....E1659
- 12 Friction rubber.....E1667
- 13 Hinge bolt.....E1668
- 14. Motor Brush.....CB-203