

EBA 21 Parts Listing		
Item	Description	Part No.
1	Rubber Foot, Each	E605
2	Turning Handle	E1875
3	Sight glass and gluing ring	E1323
4	Control Panel	52013
5	Frictional Rubber	E610
6	Leg Spring	52002
7	Appliance Plug	E1661
8	Fuse Holder	E1671
9	Fuse, 6.3 Amp	09946OR
10	Lid Moulding	E1874
11	Hook	E532
12	Centrifuge Chamber	E627
13	Motor Cover	E839
14	Packing Ring	52001
15	Covering Foil	E1876
16	Speed Sensor	E825
17	Motor Assy.	OCM-33
18	Motor Mount, Each	E1574
19	Anti-twist washer	E604
20.1	Imbalance Switch - Old Style	E2527
20.2	Imbalance Switch - New Style	E2590
21	Supply Board	52011
22	Frequency Converter	52012
23	Flat Ribbon, 16-pin	E835
24	Flat Ribbon, 10-pin	E842
25	Ferrite Ring	N/A
26	Mains choke coil	N/A
27	Brake Resistor	E849
28	Overtemperature fuse	E886
29	RF suppressor	E826
30	Lid Lock	20158

