

01-5

THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION OF BECKMAN INSTRUMENTS, INC. AND IS FURNISHED SUBJECT TO THE CONDITIONS THAT THE INFORMATION IS TO BE KEPT IN CONFIDENTIALITY, NOT BE REPRODUCED OR COPIED IN WHOLE OR IN PART AND NOT BE USED OR INCORPORATED IN ANY PRODUCT, EXCEPT UNDER AN EXPRESS WRITTEN AGREEMENT WITH BECKMAN INSTRUMENTS, INC.				
CHG.	E.O.NO.	DATE	BY	CHK
G		SEE SHEET 1		

							40		
							39		
							38		
							37		
							36		
							35		
							34		
						2A	33	241-961187	FAN, ENCLOSED - H5VAC
							32		
							31		
							30		
							29		
							28		
							27		
						2C	26	449-362336	DRIVE AND MOTOR MOUNT ASSY, 120V
							25		
							24		
			4		4		3B	23	121-823254 RIVET, BLIND - .094 DIA X .222 LG
							22		
							3B	21	261-234832 PLATE, FRENCH COMPLIANCE
(2)							3B	20	261-362113 PLATE, NAME RATING - 100V, 50/60HZ, GS-6
							3B	19	261-362112 PLATE, NAME RATING - 220V, 50HZ, GS-6
							3B	18	261-362111 PLATE, NAME RATING - 120V, 60HZ, GS-6
							3B	17	261-358708 PLATE, NAME RATING - 100V, 50/60HZ, SPINCHRON
(2)							3B	16	261-358707 PLATE, NAME RATING - 220V, 50HZ, SPINCHRON
							3B	15	261-358706 PLATE, NAME RATING - 120V, 60HZ, SPINCHRON
							14		
						2B	13	270-233469 LABEL, SERIAL NO.	
							12	060-927459 LABEL, YWR	
F							11	035-360515 SHIPPING LIST, GS-6 220V	
F							10	035-360514 SHIPPING LIST, GS-6 100/120V	
F							9	035-360506 SHIPPING LIST, SPINCHRON 220V	
F							8	305-360505 SHIPPING LIST, SPINCHRON 100/120V	
						3A	7	270-361976 OVERLAY, GS-6	
						3A	6	270-362250 OVERLAY, SPINCHRON	
							5		
							4		
							3		
							2		
							1	449-358701 CENTRIFUGE ASSY, COMMON	

QTY	ZONE	ITEM	PART NO.	DESCRIPTION
(2)		GS-6 CENTRIFUGE 100V, 50/60 HZ 050-362110		
(2)		GS-6 CENTRIFUGE 220V, 50 HZ 050-362109		
(2)		GS-6 CENTRIFUGE 120V, 60 HZ 050-362108		
(2)		SPINCHRON CENTRIFUGE 100V, 50/60 HZ 050-358704		
(2)		SPINCHRON CENTRIFUGE 220V, 50 HZ 050-358703		
(2)		SPINCHRON CENTRIFUGE 120V, 60 HZ 050-358702		
		BY	DATE	BECKMAN BECKMAN INSTRUMENTS, INC. SPINO DIVISION ONE FIVE FIVE FIVE FIVE SIX FIVE SIX FIVE SIX FIVE SIX
		OR NGAL	11-29-90	
		CHK		
		DSGN		
		ENGR		TITLE
				TABULATION, CENTRIFUGE ASSY
		C SIZE	CODE IDENT NO 07978	DWG NO. 035-358700
		MOD SPINCHRON/GS-6	SCALE	1st USE

C 035-358700
DWG NO.

THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION OF BECKMAN INSTRUMENTS, INC. AND IS LOANED SUBJECT TO THE CONDITIONS THAT THE INFORMATION BE KEPT IN CONFIDENCE, BE NOT BE REPRODUCED OR COPIED IN WHOLE OR IN PART AND BE NOT BE USED OR INCORPORATED IN ANY PRODUCT, EXCEPT UNDER AN EXPRESS WRITTEN AGREEMENT WITH BECKMAN INSTRUMENTS, INC.

CHG.	E.O.NO.	DATE	BY	CHK
G				

SEE SHEET 1

								80		
								79		
(C)							3C	78	449-360549	LATCH ASSY
								77		
								76		
								75		
								74		
								73		
								72		
								71		
								70		
								69		
								68		
								67		
							2B	66	604-927526	RESISTOR, W.W. 13, 100W, 5%
							2B	65	604-927839	RESISTOR, W.W. .5, 25W, 5%
								64		
								63		
							4C	62	701-927905	SWITCH, ROCKER - DPST 220V (POWER)
							4C	61	701-927904	SWITCH, ROCKER - DPST 120V (POWER)
								60		
								59		
							3B	58	664-355810	CORD, POWER-EUROPEAN-STR, PLUG
							3B	57	664-889097	CORD, POWER-60HZ
								56		
								55		
								54		
								53		
								52		
								51		
								50		
(F)							(I)	49	662-358307	HARNESS, AUTO TRANSFORMER 220V/JAPAN
								48		
								47		
							3B	46	792-356029	AUTO TRANSFORMER 100V TO 120V
							1A	45	792-362338	AUTO TRANSFORMER 220V TO 120V
								44		
								43		
							2B	42	795-355823	CHOKE, 10 MH
								41		

QTY		ZONE	ITEM	PART NO.	DESCRIPTION
GS-6 CENTRIFUGE 100V, 50/60 HZ 050-362110	(21)		BY	DATE	BECKMAN BECKMAN INSTRUMENTS, INC. 5100 DOW ST FULLER, CALIF 94604
GS-6 CENTRIFUGE 220V, 50 HZ 050-362109	(21)		DR NGALE	11-29-90	
GS-6 CENTRIFUGE 120V, 60 HZ 050-362108	(21)		CHK		
SPINCHRON CENTRIFUGE 100V, 50/60 HZ 050-358704	(21)		DSGN		
SPINCHRON CENTRIFUGE 220V, 50 HZ 050-358703	(21)		ENGR		TITLE
SPINCHRON CENTRIFUGE 120V, 60 HZ 050-358702	(21)				TABULATION, CENTRIFUGE ASSY
REPLOTTED 4-9-91			C	CODE IDENT NO.	DWG NO.
			SIZE	07978	035-358700
MOD SPINCHRON/GS-6			SCALE	1st USE	SHEET 3 OF 8

		150			
		149			
		148			
		147			
		146			
	2B	145	1	792-358941	TRANSFORMER, POWER S.C.D.
		144			
	2D	143	1	337-358009	BELLOW, MOTOR
		142			
	2A	141	4	313-358798	FOOT, GPII S.C.D.
		140			
		139			
	2A	138	3	297-358774	SUPPORT CAPS, BOWL
	2A	137	3	350-349799	SUPPORT, BOWL
	2D	136	1	175-358730	GASKET, BOWL
	2D	135	1	349-358938	BOWL, PAINTED
		134			
		133			
	3D	132	1	449-360246	SAFETY LOCK ASSY
	3D	131	1	729-362262	ACTUATOR, WELDMENT ASSY
	3D	130	1	405-362264	GUIDING, BLOCK
	4C	129	2	325-362050	KNOB, MOLDED
	4C	128	1	161-355103	SPACER, POTENTIOMETER
	3D	127	2	161-361939	SPACER, LATCH
		126			
		125			
		124			
		123			
5	2B	122	(1)	254-361993	FOAM, BACK PANEL
	2D	121	(1)	254-355826	FOAM, BOWL LINER
	3C	120	(1)	254-355827	FOAM, CONTROL HEAD
		119			
		118			
		117			
	4B	116	(1)	254-361938	PANEL, DAMPING, FRONT
	3C	115	(2)	449-359412	FOAM ASSY, SIDE PANEL
	2A	114	(1)	449-361907	FOAM DUCT ASSY
5		113	1	060-358768	JIT KIT, FOAM
		112			
		111			
5	4C	110	(1)	855-357190	BOARD ASSY, P.W. DISPLAY
	4B	109	(1)	855-358946	BOARD ASSY, CONTROL BOARD
		108			
		107			
	1	106	(1)	662-359329	HARNESS, SWITCH-POWER/DOOR
	1	105	(1)	662-359328	HARNESS, DOOR LOCK
	1	104	(1)	662-359330	HARNESS, LOW VOLTAGE (NON-REFR)
	1	103	(1)	662-358948	HARNESS, CONTROL PANEL
	1	102	(1)	662-358947	HARNESS, POWER
5		101	1	060-358739	JIT KIT, ELECTRICAL
		100			
		99			
	1C	98	1	349-362145	DOOR ASSY, NON-REFRIG.
		97			
	1B	96	2	317-361945	HINGE ASSY
		95			
		94			
		93			
		92			
		91			
		90			
		89			
5	3A	88	(1)	445-358731	CONTROL PANEL
	2C	87	(1)	299-361942	HOUSING
		86			
	2A	85	(1)	430-358745	CHASSIS
5		84	1	060-358741	JIT KIT, MECHANICAL
		83			
10		82	1	284-359631	PACKAGING
		81			
	ZONE	ITEM	QTY	PART NO.	

CENTRIFUGE ASSY, COMMON
449-358701

5-12

DRAWING NO.
035-358700
4 OF 8

REV. C SEE SHIT
RECKMAN INSTR., INC.
DRAWING NUMBER
035-358700
SHEET 4 OF 8

F

	290			
	289			
2D	288	A/R	65-040-04	RELEASE AGENT, EPOXY PARFILM
3B	287	A/R	65-033-03	COMPOUND, HEATSINK
3C	286	A/R	65-013-01	GREASE, MOLYBDENUM DISULFIDE
	285			
	284			
	283			
	282			
	281			
	280			
	279			
	278			
	277			
	276			
	275			
	274			
	273			
	272			
	271			
	270			
	269			
	268			
2C	267	I	165-961099	WASHER, FLT - .75 O.D., .401 I.D. .03 THK CRES
	266			
2C	265	I	165-802184	WASHER, SPLIT - #10 S-Zi
2B	264	I	165-808255	WASHER, LKG, EXT #10 S-C
	263			
2A	262	4	165-810690	WASHER, FLT #6 S-Zi
	261			
	260			
	259			
2A	258	2	165-802163	WASHER, LKG, INTL #4 S-C
	257			
	256			
4B, 4C	255	9	161-897483	SPACER, CYL, NON-THD, .530 LG .25 O.D. 148 I.D.
	254			
3D	253	2	161-818024	SPACER, CYL, NON-THD, .125 LG .375 O.D. 25 I.D.
	252			
	251			
2C	250	I	160-874672	STANDOFF, INT - 1.5 LG, 10-32 THD, 375 HX
	249			
	248			
2C	247	I	105-927746	NUT, MID X 1.25 S-Zi
	246			
2A	245	4	105-878544	NUT, 1/4-20 SLFLKG-S-Zi
2A	244	3	105-848760	NUT, 1/4-20 SLFLK-S-C
	243			
4C	242	2	105-878636	NUT, 10-32 SLFLKG-S-C
	241			
2A	240	4	105-024687	NUT, 6-32 SLFLK-S-Zi
3D, 4B	239	11	105-853432	NUT, 6-32 KEP, EXT S-Zi
	238			
	237			
2A	236	4	101-804522	SCREW, .250-20 X .75 P-S-C-PHIL
	235			
3D	234	I	101-844970	SCREW, 10-32 X .375 BU-S-BLK Zi-HXS0
4D	233	2	101-961517	SCREW, DECOR - 10-32 X 2.5 BU-S-Zi-HXS0
2B	232	2	101-803257	SCREW, 10-32 X .500 FL-S-Zi-PHIL
3A, 3C	231	11	101-927750	SCREW, DECOR - 10-32 X .375 BU-S-Zi-HXS0
2A, 2B, 3B, 3A	230	52	101-878542	SCREW, 10-32 X .375 SLFLKG-S-Zi
	229			
1B	228	12	101-878635	SCREW, 8-32 X .50 SLFLKG-S-C
	227			
2A	226	4	101-804428	SCREW, 6-32 X 2.0 P-S-Zi-PHIL
2B	225	2	101-889105	SCREW, 6-32 X 25 SLFLKG-S-C
3B	224	2	101-803124	SCREW, 6-32 X 25 FH-S-Zi-PHIL
	223			
2A	222	2	101-823155	SCREW, 4-40 X .562 P-S-C-SLTD
	221			
ZONE				
CENTRIFUGE ASSY, COMMON			5-14	035-358700
449-358701				6 OF 8

REV. 6 SEE SHI 1
BECKMAN INSTR. INC.
DRAWING NUMBER
035-358700
SHEET 6 OF 8