



















Spare-Part-List**Kendro Laboratory Products**

Plant Osterode

75008157 -01**LABOFUGE 400 120V+-10% 50/60HZ HI(VWR)**

from Serial-No.

to Serial-No. 258402

Index	Partno.	Text
00101	20360047	BOX SPANNER
00103	20054407	NUT
00105	20054544	SEAL
00108	20420025	NUT
00110	20054524	BOWL
00111	20054408	COVER
00120	20150026	CONTROL BOARD
00121	20460122	SCREW 3X6
00122	20055577	FOIL
00123	20054377	FRONTPANEL
00125	70902501	WIRING - DISPLAY
00127	20310302	BOLT
00128	20310301	CLAMP
00129	20290534	NOZZLE
00130	20180581	MAINS SWITCH
00201	70054369	HOUSING - TOP
00202	70054368	HOUSING
00203	70054367	LID
00204	20050543	WINDOW
00205	20054449	HANDLE
00210	20053064	HOOK
00211	20420151	NUT M10 X 1.25
00212	20480154	WASHER
00214	20052364	HINGE
00217	20420155	ACORN NUT
00220	20310311	GAS LID STAY
00221	20250085	WASHER 6,4
00230	70908552	LID LOCK
00231	20440306	SCREW M5X40
00232	70022477	FIXING CLAMP
00234	20310316	CLIP
00235	70055643	SHEET
00245	20290454	RUBBER BAND
00303	20054530	BASIS
00304	20300059	FOOT
00306	20420022	NUT-WASHER
00310	20054540	RUBBER MOUNT
00311	10104080	WASHER (SAND)

Spare-Part-List**Kendro Laboratory Products**

Plant Osterode

75008157 -01**LABOFUGE 400 120V+-10% 50/60HZ HI(VWR)**

from Serial-No.

to Serial-No. 258402

Index	Partno.	Text
00312	20430220	SCREW M6X25
00313	20480251	SECURING WASHER
00314	20480193	WASHER
00315	20420070	NUT M 6
00316	20220530	ADHESIVE CLIP
00320	20023748	HINGE BRACKET
00324	20490287	PIN
00330	20210255	FAN
00331	20210225	PROTECTIVE GRILL
00340	20160686	RESISTOR 470R
00341	20220518	ELBOW
00342	20420125	CAGE NUT M 4
00344	70054543	BASISPLATE
00350	20901059	MAINS CABEL
00351	20220411	CORD GRIP
00354	20180138	CONNECTOR
00360	70054405	DEFLECTOR
00364	20220414	EDGE PROTECTION
00400	20210330	MOTOR
00410	20150028	SPEED CONTROL BOARD
00420	70054912	PLATE
00424	20140289	INSULATION SOCKET
00425	70054910	IMBALANCESWITCH - BASIS
00432	20250031	WASHER
00433	70900003	IMBALANCESWITCH
00443	20028940	SPRING
00450	70902502	TACH CABLE
00500	20150107	FRONTPANEL COMPL.
00501	20230189	FUSE 6,25 A SLOW
00502	70903004	EE-PROM
00503	70000467	PROCESSOR

Spare-Part-List**Kendro Laboratory Products**

Plant Osterode

75008157 -02**LABOFUGE 400 120V+-10% 50/60HZ HI(VWR)**

from Serial-No. 258403

to Serial-No.

Index	Partno.	Text
00101	20360047	BOX SPANNER
00103	20054407	NUT
00105	20054544	SEAL
00108	20420025	NUT
00110	20054524	BOWL
00111	20054408	COVER
00120	20150026	CONTROL BOARD
00121	20460122	SCREW 3X6
00122	20055577	FOIL
00123	20054377	FRONTPANEL
00125	70902501	WIRING - DISPLAY
00127	20310302	BOLT
00128	20310301	CLAMP
00129	20290534	NOZZLE
00130	20180581	MAINS SWITCH
00201	70054369	HOUSING - TOP
00202	70054368	HOUSING
00203	70054367	LID
00204	20050543	WINDOW
00205	20054449	HANDLE
00210	20053064	HOOK
00211	20420151	NUT M10 X 1.25
00212	20480154	WASHER
00214	20052364	HINGE
00217	20420155	ACORN NUT
00220	20310311	GAS LID STAY
00221	20250085	WASHER 6,4
00230	70908552	LID LOCK
00231	20440306	SCREW M5X40
00232	70022477	FIXING CLAMP
00234	20310316	CLIP
00235	70055643	SHEET
00245	20290454	RUBBER BAND
00303	20054530	BASIS
00304	20300059	FOOT
00306	20420022	NUT-WASHER
00310	20054540	RUBBER MOUNT
00311	10104080	WASHER (SAND)

Spare-Part-List**Kendro Laboratory Products**

Plant Osterode

75008157 - 02**LABOFUGE 400 120V+-10% 50/60HZ HI(VWR)**

from Serial-No. 258403

to Serial-No.

Index	Partno.	Text
00312	20430220	SCREW M6X25
00313	20480251	SECURING WASHER
00314	20480193	WASHER
00315	20420070	NUT M 6
00316	20220530	ADHESIVE CLIP
00320	20023748	HINGE BRACKET
00324	20490287	PIN
00330	20210255	FAN
00331	20210225	PROTECTIVE GRILL
00340	20160686	RESISTOR 470R
00341	20220518	ELBOW
00342	20420125	CAGE NUT M 4
00344	70054543	BASISPLATE
00350	20901059	MAINS CABEL
00351	20220411	CORD GRIP
00354	20180138	CONNECTOR
00360	70054405	DEFLECTOR
00364	20220414	EDGE PROTECTION
00400	20210343	MOTOR
00410	20150028	SPEED CONTROL BOARD
00420	70054912	PLATE
00424	20140289	INSULATION SOCKET
00425	70054910	IMBALANCESWITCH - BASIS
00432	20250031	WASHER
00433	70900003	IMBALANCESWITCH
00443	20028940	SPRING
00450	70902502	TACH CABLE
00500	20150107	FRONTPANEL COMPL.
00501	20230189	FUSE 6,25 A SLOW
00502	70903016	EE-PROM PRV. LABO 400 3016V02
00503	20140498	CPU

Spare-Part-List**Kendro Laboratory Products****75008159-01****LABOFUGE 400,120V/1N/50/60HZ,INC.ACC.**

Plant Osterode

from Serial-No.

to Serial-No. 258718

Index	Partno.	Text
00101	20360047	BOX SPANNER
00103	20054407	NUT
00105	20054544	SEAL
00108	20420025	NUT
00110	20054524	BOWL
00111	20054408	COVER
00120	20150026	CONTROL BOARD
00121	20460122	SCREW 3X6
00122	20055577	FOIL
00123	20054377	FRONTPANEL
00125	70902501	WIRING - DISPLAY
00127	20310302	BOLT
00128	20310301	CLAMP
00129	20290534	NOZZLE
00130	20180581	MAINS SWITCH
00201	70054369	HOUSING - TOP
00202	70054368	HOUSING
00203	70054367	LID
00204	20050543	WINDOW
00205	20054449	HANDLE
00210	20053064	HOOK
00211	20420151	NUT M10 X 1.25
00212	20480154	WASHER
00214	20052364	HINGE
00217	20420155	ACORN NUT
00220	20310311	GAS LID STAY
00221	20250085	WASHER 6,4
00230	70908552	LID LOCK
00231	20440306	SCREW M5X40
00232	70022477	FIXING CLAMP
00234	20310316	CLIP
00235	70055643	SHEET
00245	20290454	RUBBER BAND
00303	20054530	BASIS
00304	20300059	FOOT
00306	20420022	NUT-WASHER
00310	20054540	RUBBER MOUNT
00311	10104080	WASHER (SAND)

Spare-Part-List**Kendro Laboratory Products**

Plant Osterode

75008159-01**LABOFUGE 400,120V/1N/50/60HZ,INC.ACC.**

from Serial-No.

to Serial-No. 258718

Index	Partno.	Text
00312	20430220	SCREW M6X25
00313	20480251	SECURING WASHER
00314	20480193	WASHER
00315	20420070	NUT M 6
00316	20220530	ADHESIVE CLIP
00320	20023748	HINGE BRACKET
00324	20490287	PIN
00330	20210255	FAN
00331	20210225	PROTECTIVE GRILL
00340	20160686	RESISTOR 470R
00341	20220518	ELBOW
00342	20420125	CAGE NUT M 4
00344	70054543	BASISPLATE
00350	20901059	MAINS CABEL
00351	20220411	CORD GRIP
00354	20180138	CONNECTOR
00360	70054405	DEFLECTOR
00364	20220414	EDGE PROTECTION
00400	20210330	MOTOR
00410	20150028	SPEED CONTROL BOARD
00420	70054912	PLATE
00424	20140289	INSULATION SOCKET
00425	70054910	IMBALANCESWITCH - BASIS
00432	20250031	WASHER
00433	70900003	IMBALANCESWITCH
00443	20028940	SPRING
00450	70902502	TACH CABLE
00500	20150107	FRONTPANEL COMPL.
00501	20230189	FUSE 6,25 A SLOW
00502	70903004	EE-PROM
00503	70000477	PROCESSOR

Spare-Part-List**Kendro Laboratory Products**

Plant Osterode

75008159-02**LABOFUGE 400,120V/1N/50/60HZ,INC.ACC.**

from Serial-No. 258719

to Serial-No.

Index	Partno.	Text
00101	20360047	BOX SPANNER
00103	20054407	NUT
00105	20054544	SEAL
00108	20420025	NUT
00110	20054524	BOWL
00111	20054408	COVER
00120	20150026	CONTROL BOARD
00121	20460122	SCREW 3X6
00122	20055577	FOIL
00123	20054377	FRONTPANEL
00125	70902501	WIRING - DISPLAY
00127	20310302	BOLT
00128	20310301	CLAMP
00129	20290534	NOZZLE
00130	20180581	MAINS SWITCH
00201	70054369	HOUSING - TOP
00202	70054368	HOUSING
00203	70054367	LID
00204	20050543	WINDOW
00205	20054449	HANDLE
00210	20053064	HOOK
00211	20420151	NUT M10 X 1.25
00212	20480154	WASHER
00214	20052364	HINGE
00217	20420155	ACORN NUT
00220	20310311	GAS LID STAY
00221	20250085	WASHER 6,4
00230	70908552	LID LOCK
00231	20440306	SCREW M5X40
00232	70022477	FIXING CLAMP
00234	20310316	CLIP
00235	70055643	SHEET
00245	20290454	RUBBER BAND
00303	20054530	BASIS
00304	20300059	FOOT
00306	20420022	NUT-WASHER
00310	20054540	RUBBER MOUNT
00311	10104080	WASHER (SAND)

Spare-Part-List**Kendro Laboratory Products**

Plant Osterode

75008159-02**LABOFUGE 400,120V/1N/50/60HZ,INC.ACC.**

from Serial-No. 258719

to Serial-No.

Index	Partno.	Text
00312	20430220	SCREW M6X25
00313	20480251	SECURING WASHER
00314	20480193	WASHER
00315	20420070	NUT M 6
00316	20220530	ADHESIVE CLIP
00320	20023748	HINGE BRACKET
00324	20490287	PIN
00330	20210255	FAN
00331	20210225	PROTECTIVE GRILL
00340	20160686	RESISTOR 470R
00341	20220518	ELBOW
00342	20420125	CAGE NUT M 4
00344	70054543	BASISPLATE
00350	20901059	MAINS CABEL
00351	20220411	CORD GRIP
00354	20180138	CONNECTOR
00360	70054405	DEFLECTOR
00364	20220414	EDGE PROTECTION
00400	20210343	MOTOR
00410	20150028	SPEED CONTROL BOARD
00420	70054912	PLATE
00424	20140289	INSULATION SOCKET
00425	70054910	IMBALANCESWITCH - BASIS
00432	20250031	WASHER
00433	70900003	IMBALANCESWITCH
00443	20028940	SPRING
00450	70902502	TACH CABLE
00500	20150107	FRONTPANEL COMPL.
00501	20230189	FUSE 6,25 A SLOW
00502	70903015	EE-PROM
00503	20140498	CPU