


EXPLOSIONSZEICHNUNGEN
BREAKDOWNDRAWINGS



Spare-Part-List

Kendro Laboratory Products

Plant Osterode

75003236.01

BIOFUGE PICO 120V 50/60 HZ HERAEUS

from Serial-No.

to Serial-No.

Index	Partno.	Text
00003	20480242	SECURING WASHER
00004	20480244	LOCK WASHER
00006	20480251	SECURING WASHER
00030	20420069	HEXAGONAL NUT M3
00034	20460098	SCREW M 3 X 6
00037	20460094	SCREW M 4 X 8
00049	20440288	SCREW M 3 X 6
00055	20440285	SCREW
00056	20440298	SCREW M4X12
00058	20440307	SCREW M6X12
00060	20420070	NUT M 6
00072	20510285	SCREW M5 X 12
00094	20480108	FAN WASHER A 4,3
00100	20055918	HOUSING
00102	20056212	SEAL
00104	20360072	SPANNER 10 MM
00105	20290607	O-RING
00106	70056208	NUT
00110	70056282	LID
00111	20050543	WINDOW
00112	20058002	Handle
00114	20055924	NUT
00115	20350211	COLLAR
00116	70055921	AXIS
00117	20055922	SPRING (RIGHT)
00118	20055923	SPRING (LEFT)
00119	20055920	LID SEAL
00120	20901141	Latch 120 V JST
00121	20420164	Nut M6
00130	75003328	MICROLITRE ROTOR 24 X 1.5/2ML
00131	75003329	ROTOR LID AEROSOL TIGHT
00200	20055917	CHASSIS
00201	20058001	Overlay Biofuge pico
00202	20300101	Foot
00204	20190247	MAINSCABLE
00205	20180628	SOCKET
00207	20055928	SPACER
00208	20220416	CABLE-CLIP

Spare-Part-List**Kendro Laboratory Products**

Plant Osterode

75003236.01**BIOFUGE PICO 120V 50/60 HZ HERAEUS**

from Serial-No.

to Serial-No.

Index	Partno.	Text
00210	20220580	SPACER
00216	20180138	CONNECTOR
00220	20150293	Mainboard green display
00221	20230136	FUSE 4 AT
00222	70901072	NV-RAM
00223	20055932	INSULATION FOIL
00230	20210357	MOTOR
00231	20490376	BRACE
00232	20290567	O-RING
00233	20270076	SAFTY PAINT
00234	20480301	Sand paper washer
00235	20300074	ANTIVIBRATION MOUNT

Spare-Part-List

Kendro Laboratory Products

Plant Osterode

75003238.01

Biofuge pico 120V 50/60 Hz Sorvall

from Serial-No.

to Serial-No.

Index	Partno.	Text
00003	20480242	SECURING WASHER
00004	20480244	LOCK WASHER
00006	20480251	SECURING WASHER
00030	20420069	HEXAGONAL NUT M3
00034	20460098	SCREW M 3 X 6
00037	20460094	SCREW M 4 X 8
00049	20440288	SCREW M 3 X 6
00055	20440285	SCREW
00056	20440298	SCREW M4X12
00058	20440307	SCREW M6X12
00060	20420070	NUT M 6
00072	20510285	SCREW M5 X 12
00094	20480108	FAN WASHER A 4,3
00100	20055918	HOUSING
00102	20056212	SEAL
00103	75003327	FIXED ANGLE ROTOR 4 X PCR STRIP
00104	20360072	SPANNER 10 MM
00105	20290607	O-ring
00106	70056208	NUT
00110	70056282	LID
00111	20050543	WINDOW
00112	20058002	Handle
00114	20055924	NUT
00115	20350211	COLLAR
00116	70055921	AXIS
00117	20055922	SPRING (RIGHT)
00118	20055923	SPRING (LEFT)
00119	20055920	LID SEAL
00120	20901141	Latch 120 V JST
00121	20420164	Nut M6
00130	75003328	MICROLITRE ROTOR 24 X 1.5/2ML
00131	75003329	ROTOR LID AEROSOL TIGHT
00200	20055917	CHASSIS
00201	20058016	Overlay
00202	20300101	Foot
00204	20190247	MAINSCABLE
00205	20180628	SOCKET
00207	20055928	SPACER

Spare-Part-List**Kendro Laboratory Products**

Plant Osterode

75003238.01**Biofuge pico 120V 50/60 Hz Sorvall**

from Serial-No.

to Serial-No.

Index	Partno.	Text
00208	20220416	CABLE-CLIP
00210	20220580	SPACER
00216	20180138	CONNECTOR
00220	20150293	Mainboard
00221	20230136	FUSE 4 AT
00222	70901072	NV-RAM
00223	20055932	INSULATION FOIL
00230	20210357	MOTOR
00231	20490376	BRACE
00232	20290567	O-RING
00233	20270076	SAFTY PAINT
00234	20480301	Sand paper washer
00235	20300074	ANTIVIBRATION MOUNT