

ILLUSTRATED PARTS BREAKDOWN,
SPARES and ACCESSORIES
For
READACRIT and MP READACRIT Centrifuges




Figure 28. Illustrated Parts Breakdown,
READACRIT and MP READACRIT Centrifuges.

READACRIT 420591 & MP READACRIT 420593

#	DESCRIPTION	CATALOG NUMBER	REFERENCE NUMBER
N/A	Operator's manual	421451	
1	Motor, 110V	421449	42059101
2	Motor Brush assembly, 2 pk	421450, regular motor	42059102
2a	Motor Brush assembly, 2 pk	421474, "special kit" noted on motor	42059126
3	Timer w/ knob	421479	42059202
4	Knob only (for timer)	421480	42059203
5	Safety switch assembly	421454	42059106
13	Magnifier assembly	421465	42059117
14	Rubber channel	421466	42059118
15	Rubber foot	421204	42010113
17	Rubber grommet	421469	42059121
18	Rubber bumpers	421471	42059123
19	Wrench	421472	42059124
21	Balance tube, sedical tray	421482	42059302
<u>Discontinued parts</u>		<i>DQ'd 8/02</i>	
7	Cover assembly	421476	42059128
8	Catch assembly	421473	42059125
9	Head cover assembly, Gullwing	421461	42059113
10	Head cover, conversion kit	421462	42059114
11	Head cover, hardware kit	421463	42059115
12	Head assembly	421464	42059116
16	Bottom plate & motor bracket	421477	42059129
20	SEDI-CAL tray & cover assembly	421481	42059301
22	Two-place SEDI-CAL slide	421483	42059303
24	Latch	421459	42059111
6	Cordset assembly	no part # , not orderable	
N/A	Foot, old style	421467	42059119
N/A	Knob only, gullwing	421452	42059104
N/A	Timer Conversion kit	421475	42059127

*Part #s apply to Readacrit 420591 and MP Readacrit 420593