

Baxter / Heraeus Centrifuges

Clinifuge/Labofuge 200



PM Kit
OEM:70901069
P/N:38177



Lid NEW / Exchange
Site glass not included.
OEM:70020803
P/N:38526



Rubber Foot
OEM:20049932
P/N:38532

Haemofuge A

Haemofuge



Brush/Pair
P/N:CB-117



Motor Exchange
Rebuild Prices Available!
OEM:1414502
P/N:OCM-01E



Induction Motor Exchange
Rebuild Prices Available!
OEM:20210301
P/N:OCM-01T

Labofuge 300



New Motor Exchange
OEM:20210384
P/N:30150



Motor Exchange
Rebuild Prices Available!
OEM:20210384
P/N:OCM-01V



NEW Latch
Fits Models S/N
4003672 and below.
OEM:20901101
P/N:38173



New Latch
Fits Models S/N
2003673 and above.
OEM:20901131
P/N:20143



Gas Spring
OEM:20310428
P/N:38162



Lid Seal
OEM:20056506
P/N:38142



Hinge Bracket
OEM:20023748
P/N:38131