

# Becton Dickinson/Clay Adams

Becton Dickinson/  
Clay Adams

## Triac



**New Motor Exchange**  
With Gear  
S/N 145000 & Above  
OEM:421282  
**P/N:30122**




**NEW Motor Exchange**  
220V Motor Without Gear  
OEM:420601  
**P/N:30123**



**New Motor Exchange**  
Without Gear  
S/N 145000 & Below  
OEM:421282  
**P/N:30132**



**Brush, Pair**  
OEM:4202002  
**P/N:CB-003**



**Head Assy. w/ Cover**  
OEM:420202  
**P/N:420202**



**Tube Adapter**  
2ml 4-pack  
OEM:420203  
**P/N:420203**



**Tube Adapter**  
3ml 4-pack  
OEM:420204  
**P/N:420204**



**Motor Gasket**  
OEM:421264  
**P/N:421264**



**Timer**  
OEM:421265  
**P/N:421265T**



**NEW Lid Assembly**  
OEM:421268  
**P/N:421268**



**Cover Gasket**  
OEM:421269  
**P/N:421269**



**Mode Selector**  
OEM:421270  
**P/N:421270**



**Head Fastner Nut**  
OEM:421272  
**P/N:421272**



**Lower Head Cover**  
OEM:421273  
**P/N:421273**



**Wrench Spanner**  
OEM:421275  
**P/N:421275**



**Trunnion Cushion**  
(pk of 8)  
OEM:421276  
**P/N:421276**



**Tube Cushion**  
(pk of 8)  
OEM:421278  
**P/N:421278**



**Head Washer**  
OEM:421279  
**P/N:421279**



**Head Cover w/ Knob**  
OEM:421280  
**P/N:421280**



**Safety Switch**  
For S/N 145000 & Above  
OEM:421283  
**P/N:421283T**



**Tach Sensor**  
OEM:421284  
**P/N:421284**



**Rubber Feet**  
(pk of 4)  
OEM:421285  
**P/N:421285T**



**Transformer**  
OEM:421286  
**P/N:421286**



**PCB**  
For S/N 145000 & Above  
OEM:421287  
**P/N:421287**