

Brush Application List

Brush Prices Are Per Pair Unless Otherwise Noted

<u>Model</u>	<u>Notes</u>	<u>OEM Part #</u>	<u>Ozark Part #</u>
IEC			
2V, 2R, BPV, BPR		1724	CB-008
B-2, B-20, B35		2467	CB-010
B-20A	IEC Discontinued, 2 section type	44908	CB-012
B-20A	New Style, solid type	48264	CB-013
B-22		41775	CB-080
B-60	IEC Discontinued	42560	CB-034
C-Series	IEC Discontinued	1707	CB-024
CL, CH	Use P/N CB-006A	1780	CB-006
CL, CH, Centra CL2	New Style W/Shunt Wire	1780A	CB-006A
CM, C50		1731	CB-014
Centra-CLD		75202	CB-150
Centra-4		49203	CB-079
Centra-4B	Program 1.0 Motor W/White Caps	49374	CB-071
Centra-4B, CL4	Program 2.0 Motor W/Blk Caps	44989	CB-073
Centra-7	1/4x7/16x3/4 Old 2 section type	47263	CB-029
Centra-7	1/4x3/8x3/4	47889	CB-045
Centra-7	3/16x3/8x7/8 Most Units	48192	CB-035
Centra-8/8R, GP8/8R	IEC Motor	49801	CB-072
Centra GP8/8R	Baldor Motor	T01206S	CB-304
CU/CRU-5000	IEC Discontinued 2 section type	44908	CB-012
CU/CRU-5000	New Style Solid Type	48264	CB-013
C/CR-6000		41775	CB-080
H, HN-S, HN-SII	Use P/N CB-006A	1780	CB-006
HN-SII	New Style With Shunt Wire	1780A	CB-006A
HT		1860	CB-028 Discontinued
Equafuge		50757	CB-163
MB Microhematocrit	Specialty Motor	42055813	CB-001B
MB Microhematocrit	Ser. # up to AA2914		
	And Above AC6200	1840	CB-007
MB Microhematocrit	Jones Motor SerAA2914-AC6200	43653	CB-007A
MH			CB-006A
M-60	Round Motor	41822	CB-142 Discontinued
M-60	Square Motor	49846	CB-143
MP4, MP4R	Program 1.0 Motor W/White Caps	49374	CB-071
	Serial #s 24370600 & Below.		
	Refrig Serial #s 24380600 & Below		
MP4, MP4R	Program 2.0 Motor W/Blk Caps	44989	CB-073
	Serial #s 24371999 & Below.		
	Refrig Serial #s 24381999 & Below		
MP4, MP4R	Inverted Motor	50433	CB-162
	Serial #s 24371999 & Up.		
	Refrig Serial #s 24381999 & Up		
PR-1, PR-2	IEC Discontinued	1724H	CB-009
PR-6, PR-J		40069	CB-011
PR/DPR 6000	IEC Discontinued Old Style	44908	CB-012