



As an alternative to purchasing new units, we offer top quality reconditioned centrifuges.

Our reconditioned centrifuges look and operate like new, at a fraction of the price.

Some of the many reasons to consider a reconditioned centrifuge

1. Your Customer is important to us too!
2. Significant Savings--thoroughly reconditioned with warranty, providing maximum dollar value.
3. We will only provide instruments that have a reliable track history.
4. We provide service and support after the sale.
5. We will not sell you someone else's problems.
6. We take centrifuges in on trade.
7. We can solve your centrifuge problems and get parts on site in a timely fashion.
8. We help minimize down time.
9. We specialize in centrifuges, and have for over 25 years.

This is what you deserve in centrifuge sales and service!

Reconditioning Process

Our reconditioning process addresses all major components and subassemblies of the centrifuge, including the motor and latch assemblies. Touch panel overlays are replaced. Rotors are inspected and reconditioned.

Any mechanical timers and switches are replaced. Motor mounts are always replaced. These components are replaced with quality new and reconditioned components from our stock. Because we sell parts to customers worldwide you can be assured of top quality.

The centrifuge is then inspected for any cosmetic defects. Paint is matched to original manufacturers color by our in-house paint technicians who restore the units to "like new" appearance. The finished unit passes on to our quality assurance technicians. A complete alignment of the electronics is performed assuring proper acceleration, deceleration, speed detection, timing and imbalance detection.

Finally, our refurbished units are shipped utilizing a combination of our Sealed Air Instapak foam, custom built pallets and heavy duty cardboard boxes.

We ship any centrifuge over 100 pounds via trucklines individually palletized and packaged to assure safe arrival at your site.

At Ozark Biomedical, we are experts in assuring your unit gets to you in perfect condition.

