

Ozark Biomedical

2012 Catalog of Centrifuges, Motors & Parts

1001 Commerce Place - Beebe, Arkansas 72012

Toll Free: 800-457-7576 Fax: 501-882-2122

Visit us online at www.ozarkbiomedical.com

ISO 9001:2008 Registered Company



Your Centrifuge Specialty Store
Solving Centrifuge Problems Since 1978

We are an international company with products in all 50 states and 22 countries. Since its inception Ozark Biomedical has been driven by the customers demand for one company to turn to with all their centrifuge needs. Our mission is to be that one company!

What can we do for you?

**New or Reconditioned Centrifuges, Repair Depot for Motors, Latches, Circuit Boards.
Remanufactured Motors, Spindles, Latches.**

New and Repaired Circuit Boards, New and Used Parts, Technical Support.

Our quality reconditioned motors, spindles, latches and circuit boards provide savings often exceeding 50% below OEM prices. Our motor brushes are a fraction of the cost of OEM brushes.

This adds up to significant savings over the life of your instrument.

As specialists in this field we are focused on helping you diagnose your centrifuge problems and providing you with the products you need, efficiently and at significant cost savings to your facility. Our good working relationships with the many centrifuge manufacturers assures you we can find a solution to your need, regardless of the severity.

We carry parts for most centrifuges including

Abbott, Baxter, Becton Dickinson, Beckman, Bio-Dynamics, Brinkman, Clay Adams, Dade, Drucker, Global Focus, Hamilton Bell, Hermle, IEC, Jouan, Kendro, Mistral, Sorvall

The next time you have a centrifuge problem you only have one call to make.

Toll Free: 1-800-457-7576 - Fax: 866-340-8496

The Savings will be significant!




**Front Cover taken by Tim Ernst:
"Keefe Falls along the Ozark Highlands
Trail"**

Visit his website at cloudland.net



As an alternative to purchasing new, we offer top quality reconditioned centrifuges.
Our reconditioned centrifuges look and operate like new at a fraction of the price.

Consider a Reconditioned Centrifuge from



Latches

Circuit Boards

Motors



Brushes

Switches
Motor Mounts
Lid Stays
Timers
Overlays
Fuses
Tachometers
Bearings
Armatures
New Parts
Used Parts
Hard-to-Find Parts
Wiring Harnesses
Relays
Solenoids



Spindles & Gyros

OZARK CENTRIFUGE MOTORS



With a 6 month Limited Warranty
our Replacement Motors offer the best warranty in the business!
Our Quality Assurance program pinpoints failures and
eliminates the cause, providing failure rates of
less than 2%

Visit Us at www.ozarkbiomedical.com



**= Quality Products You Can Count On!
Accept No Substitute!**



These motors don't belong in a landfill

We are environmentally friendly. We have recycled over
20,000 centrifuge motors since 2000.

Before you throw another motor away, realize there is a better way!

Visit Us at www.ozarkbiomedical.com



Centrifuges, Motors & Replacement Parts

Our motors meet or exceed the specifications of new motors.

The remanufacturing process includes all components necessary to bring the motor back to original specifications.

Our process includes new brushes, bearings, and any necessary motor housing repair or shaft replacement. A rewind armature with new commutator is used if needed.

All rebuilt armatures are balanced as well as dipped and baked to assure maximum performance. All motors are bench tested and most are instrument tested under a load.

All motors and spindles come with a 6-month full replacement warranty.

We sell motors and armatures 2 ways:

Exchanges are shipped the day ordered. Upon receiving the product we ask that you return your exchange in the same box within 15 days. All paper work and return shipping label are enclosed for your convenience. Your help returning Exchange Items on time is greatly appreciated!

Rebuild: You send your product to us for repair and receive the remanufactured product in return. Our turn-a-round is usually 1 day. Call for RGA# before sending any product to us! We recommend buying exchange motors when time is of utmost importance. It is not cost effective to request a rebuild over an exchange if you are going to ship priority both ways. With the exchange you can have your motor the next morning by 10:30AM.



Abbott Centrifuges

3531

Abbott/Ametek



Lid Lock Assy
120V Solenoid
OEM:20003621
P/N:38122



Front Lid Seal
OEM:20020712
P/N:38123



Rear Lid Seal
OEM:20020713
P/N:38124



Lid Sight Glass
OEM:20050543
P/N:38134



Hinge Kit
OEM:20037565
P/N:38133



Motor Exchange
Rebuild Price Available!
OEM:20210290
P/N:OCM-20B



Rubber Foot
OEM:20049932
P/N:38532



Outer Housing
OEM:20020726
P/N:38272

Tachometer & Accessories

Speed Products



Superlube 3oz.
OEM:38551
P/N:38551



Reflective Tabs
25 Dots per Sheet
OEM:51002
P/N:51002



Reflective Tape
30"x.5" Strip
OEM:51004
P/N:51004

Did you know?

Ozark Biomedical is a ISO9001:2008 Registered Company.

**It's just part of our continued commitment to excellence
and customer satisfaction.**

Baxter / Heraeus Centrifuges

Medifuge 3611/3632

Biofuge Primo



Motor Exchange
Rebuild Price Available!
OEM:210276
P/N:OCM-01M



Lid Seal
Fits Non-Refrig. Units
OEM:20056506
P/N:38142



Lid Hook
OEM:20056527
P/N:38524

Biofuge 13 / Biofuge 15



NEW Latch Exchange
220V Solenoid. Confirm Solenoid
Voltage Before Ordering!
OEM:70003700
P/N:20113



Latch Exchange
220V Solenoid. Confirm Solenoid
Voltage Before Ordering!
OEM:70003700
P/N:20114



Lid Lock Assy
120V Solenoid.
OEM:20003621
P/N:38122



Biofuge 15 Motor Rebuild
Must Call for RMA#
OEM:2120277
P/N:OCM-01P



Biofuge 13 Motor Exchange
Rebuild Price Available!
OEM:20210290
P/N:OCM-20B



Front Lid Seal
OEM:20020712
P/N:38123



Rear Lid Seal
OEM:20020713
P/N:38124



Motor Mounts
Sold Individually
OEM:20020732
P/N:38125



Hinge Kit
OEM:20037565
P/N:38133



Rubber Foot
OEM:20049932
P/N:38532



Lid Sight Glass
OEM:20050543
P/N:38134

Baxter / Heraeus Centrifuges

Biofuge Pico



NEW Latch

Fits Models S/N 40017506 & Above
OEM:20901101
P/N:38173



Motor Exchange

Rebuild Price Available!
OEM:20210357
P/N:OCM-01AA

Clinifuge/Labofuge 200



New Motor Exchange

OEM:20210275
P/N:30130



Motor Exchange

OEM:20210275
P/N:OCM-01Q



Roll Pin for Motor Shaft

OEM:RP-01Q
P/N:38120



Lid Lock Assy.

OEM:20003621
P/N:38122



Lid Gasket

OEM:20020717
P/N:38126



Front Panel

OEM:20022598
P/N:38127



Antivibration Mount

OEM:20022874
P/N:38128



Key Foil

OEM:20022901
P/N:38129



Hinge Kit

OEM:20037565
P/N:38133



Lid Sight Glass

OEM:20050543
P/N:38134



Rotor Nut

OEM:20053358
P/N:38138



Power Switch

OEM:20180581
P/N:38153



Start Capacitor

OEM:20540344
P/N:38164



New Main PCB

Sold Exchange
OEM:20150178
P/N:38168



Main PCB Exchange

For PCB's Ending in 17801
OEM:20150178
P/N:38169



Main PCB Exchange

For PCB's Ending in 07601
OEM:20150178
P/N:38170

Baxter / Heraeus Centrifuges

Clinifuge/Labofuge 200



PM Kit
OEM:70901069
P/N:38177



Lid NEW / Exchange
Site glass not included.
OEM:70020803
P/N:38526



Rubber Foot
OEM:20049932
P/N:38532

Haemofuge A

Haemofuge



Brush/Pair
P/N:CB-117



Motor Exchange
Rebuild Prices Available!
OEM:1414502
P/N:OCM-01E



Induction Motor Exchange
Rebuild Prices Available!
OEM:20210301
P/N:OCM-01T

Labofuge 300



New Motor Exchange
OEM:20210384
P/N:30150



Motor Exchange
Rebuild Prices Available!
OEM:20210384
P/N:OCM-01V



NEW Latch
Fits Models S/N
4003672 and below.
OEM:20901101
P/N:38173



New Latch
Fits Models S/N
2003673 and above.
OEM:20901131
P/N:20143



Gas Spring
OEM:20310428
P/N:38162



Lid Seal
OEM:20056506
P/N:38142



Hinge Bracket
OEM:20023748
P/N:38131

Baxter / Heraeus Centrifuges

Labofuge 400



Motor Exchange
Rebuild Price Available!
OEM:20210343
P/N:OCM-01X



Latch Exchange
For S/N 221056 & below
OEM:50079829
P/N:20111



Solenoid coil
120V Solenoid for P/N 20111
OEM:5001004
P/N:20112



Latch Exchange
Fits Models S/N
221056 & above w/
P/N 70908552 on Latch
OEM:70908552
P/N:20139



Hinge
OEM:20052364
P/N:38135



Lid Hook
OEM:20052842
P/N:38137



Lid Gasket
OEM:20054544
P/N:38140



Front Panel Overlay
OEM:20055577
P/N:38141



Display PCB Exchange
OEM:20150026
P/N:38145



Main PCB Exchange
For Non-Refrig. Units
OEM:20150107
P/N:38533



Gas Lid Stay
Fits Non-Refrig. Units.
OEM:20310311
P/N:38159



Gas Lid Stay
Fits Refrig. Units
OEM:20310318
P/N:38160



Front Panel Assy.
OEM:20054377
P/N:38171



Main PCB Exchange
w/ White Connectors
for Labofuge 400/R.
OEM:20151092
P/N:38580

Labofuge A



Lid Stay
OEM:20310311
P/N:38159



Brush, Pair
OEM:001754
P/N:CB-052



Motor Exchange
Rebuild Prices Available!
OEM:210181
P/N:OCM-01A

Baxter / Heraeus Centrifuges

Labofuge B



Brush, Pair
OEM:02228
P/N:CB-053



Motor Exchange
Rebuild Price Available!
OEM:20210173
P/N:OCM-01D

Omnifuge RT



Brush/Pair
OEM:003854
P/N:CB-102



Motor Exchange
OEM:210212
P/N:OCM-01K


Megafuge 2630



Latch Exchange
OEM:70902521
P/N:20115



Tach Pick-up
OEM:20001678
P/N:38121



Lid Hook
OEM:20052842
P/N:38137



Hall Effect Vane Sensor
OEM:CYHME301
P/N:38254



Brush Pair
OEM:002790
P/N:CB-094



Motor, Exchange.
Rebuild Price Available!
OEM:210258
P/N:OCM-01C

Medifuge



Motor Assy NEW
Sold Exchange
OEM:20210205
P/N:30121



Motor Exchange
Rebuild Price Available!
OEM:20210205
P/N:OCM-01B



Brush, Pair
OEM:001193
P/N:CB-049

Baxter / Heraeus Centrifuges

Megafuge 1.0 / 2.0



Lid Seal
For Refrig. Units
OEM:20290516
P/N:38172



Collet Chuck
OEM:70003100
P/N:38175



Tach Cable
OEM:70902502
P/N:38179



New Imbalance Switch
Fits 220V Megafuge 1.0
4.5" long w/ 2-pin connector
OEM:70902503
P/N:38180



Main PCB Exchange
Fits Non-Refrig. Models
3491 S/N 208897 and below
3492 S/N 210905 and below
3493 S/N 212041 and below
P/N:38149



Main PCB Exchange
with white connectors
OEM:20150090R
P/N:38528



NV Ram
For P/N 38149
OEM:70903400
P/N:38184



NV Ram
For P/N 38534 & 38528
Rev. V03
OEM:70903010
P/N:38183



PM Kit
Fits Non-Refrig Units
OEM:70902656
P/N:38182



Imbalance Switch
12.5" long w/ 4-pin connector
OEM:70062556
P/N:38510



Speed Nut Kit
OEM:082604
P/N:38185



Foil Overlay
Fits Non-Refrig. Units
OEM:20055572
P/N:38187



Rotor Lid Pins
For 3360 Rotor
OEM:20510224
P/N:38520



Latch Exchange
For S/N 221056 & below
OEM:50079829
P/N:20111



Solenoid coil
120V Solenoid for P/N 20111
OEM:5001004
P/N:20112



NEW Latch Exchange
220V Solenoid. Confirm Solenoid
Voltage Before Ordering!
OEM:70003700
P/N:20113



Latch Exchange
220V Solenoid. Confirm Solenoid
Voltage Before Ordering!
OEM:70003700
P/N:20114



Motor Exchange
Rebuild Price Available!
5-wire Motor
OEM:20210332
P/N:OCM-01N



Motor Exchange
Rebuild Price Available!
7-wire Motor
OEM:20210317
P/N:OCM-01N-1

Beckman Centrifuges

Allegra 6



Bowl Gasket
OEM:349745
P/N:10025



Control PCB Exchange
OEM:358946
P/N:10027



Safety Lock Exchange
OEM:363414
P/N:20118



Thermistor Assy
OEM:349787
P/N:10033



Display PCB Exchange
OEM:357190
P/N:10043



Motor Boot
OEM:358009
P/N:10044



Striker
OEM:358020
P/N:10045



Speed Limiter Kit
OEM:364239
P/N:10065



Timer 30 min.
OEM:927457
P/N:10088



Fan, Chasis
OEM:961189
P/N:10090



Brake Resistor
13 Ohm, 100w
OEM:927526
P/N:10114



Speed/Temp
OEM:961314
P/N:10122



Lower Mount Spacer
OEM:MP-1104
P/N:10170



Latch Assy. Exchange
Non-Refrig. Models
OEM:360548
P/N:20116



Latch Assy. Exchange
Refrig. Models
OEM:360548
P/N:20117



Hinge Kit
OEM:360744
P/N:10173



Brush, Pair
Solid Type
OEM:356527
P/N:CB-081



Brush, kit
Solid Type
OEM:356527
P/N:CB-081A



Brush, Pair
Autostop
OEM:364241
P/N:CB-155



Brush, kit
Autostop Type
OEM:364241
P/N:CB-155A



Motor Exchange
Rebuild Price Available!
OEM:367369
P/N:OCM-02R

Beckman Centrifuges

Allegra X-22



Tach PCB
Rectangular Exchange
OEM:361349
P/N:10049



Tach Pick-Up
Square
OEM:368253
P/N:10130



Motor Mount
Each
OEM:961803
P/N:10119



Motor Boot
OEM:364206
P/N:10064



Power Switch
OEM:392236
P/N:10148



Power PCB Exchange
For 120V Units
OEM:365738
P/N:10129



New Motor Exchange
OEM:392223
P/N:30153



Motor Exchange
Rebuild Price Available!
OEM:392223
P/N:OCM-02X

Allegra X-12



Motor NEW Exchange
OEM:392357
P/N:30152



Motor Exchange
Rebuild Price Available!
OEM:392357
P/N:OCM-02Z



Motor Rebuild Only
OEM:368250
P/N:OCM-02Y

Allegra 25/25R

Accuspin



Belt Drive
OEM:885688T
P/N:10086



Spindle Assy
OEM:343030
P/N:20514



Brush, Pair
For Elinico Motor
OEM:349544
P/N:CB-043



Brush, Pair
For Carter Motor
OEM:349869
P/N:CB-044



Carter Motor Exchange
Rebuild Price Available!
OEM:345517
P/N:OCM-02C

Beckman Centrifuges

J6 Series



Magnetic Disk
OEM:343892
P/N:10021



Spindle Pulley
OEM:369312
P/N:10075




Motor, Exchange
For J6/B/BP
Rebuild Price Available!
OEM:368610
P/N:OCM-02



Motor, Exchange
For J6-M
Rebuild Price Available!
OEM:343897
P/N:OCM-02K



Drive Belt
3/16" thick
OEM:870961
P/N:10105



Cover Handle
OEM:344589
P/N:10141



Spindle
OEM:364937
P/N:20508



Handle Assy.
OEM:344667
P/N:10142



Door Closed Switch
OEM:885582
P/N:10099



Brush, Pair
with Indicator Wire
Thick Brush
OEM:364934
P/N:CB-041



Brush, Pair
with Indicator Wire
Thin Brush
OEM:349454
P/N:CB-144

Microfuge 11 / 12



Timer Assy.
OEM:360859
P/N:10060



New Motor Exchange
OEM:344451
P/N:30174



Motor Exchange
Rebuild Price Available!
OEM:344451
P/N:OCM-02F



Brush Pair
OEM:343662
P/N:CB-046

Beckman Centrifuges

GP

Beckman



Latch Spacer
OEM:349797
P/N:10034



Bowl Mount Stand-Off
OEM:349799
P/N:10035



Motor Mount Loop
OEM:355821
P/N:10036



Motor, Exchange
Rebuild Price Available!
OEM:355969
P/N:OCM-02A



Bowl Gasket
OEM:349745
P/N:10025



Vibration Mount Spring
OEM:349756
P/N:10026



Tach Sensor
OEM:349785T
P/N:10032



Speed/Temp/Timer Knob
OEM:974672
P/N:10030



Thermistor Assy
OEM:349787
P/N:10033



Hinge Kit
OEM:360744
P/N:10173



Motor Grommet
OEM:349786
P/N:10160



Brake Switch
OEM:927462
P/N:10039



Motor Boot
OEM:358009
P/N:10044



Brake Switch
OEM:927462
P/N:10039



Vibration Mount Spacer
OEM:362274
P/N:10050



Latch Assy. Exchange
For Non-Refrig. Units
OEM:360548
P/N:20116



Latch Assy. Exchange
For Refrig. Units
OEM:360548
P/N:20117



Brush, Pair
Solid Type
OEM:356527
P/N:CB-081



Brush, kit
Solid Type
OEM:356527
P/N:CB-081A



Brush, Pair
Autostop
OEM:364241
P/N:CB-155



Brush, kit
Autostop Type
OEM:364241
P/N:CB-155A



PCB Assy Exchange
OEM:358946
P/N:10024



Display PCB Exchange
OEM:357190
P/N:10029



Vibration Mount Kit
OEM:N/A
P/N:10146

Beckman Centrifuges

GP



Vibration Mount Hardware

OEM:N/A

P/N:10085



Lid Open Switch

OEM:927461

P/N:10110



Bridge Rectifier

OEM:927473

P/N:10115



Timer 30 minute

OEM:927457

P/N:10088



Speed/Temp

OEM:961314

P/N:10122



Imbalance Switch

OEM:927901

P/N:10111



115V Chassis Fan

OEM:961189

P/N:10134



Brake Resistor

13 Ohm 100W


OEM:927526
P/N:10114



Speed Limiter Kit

OEM:364239

P/N:10065



Roll Pin

OEM:964556

P/N:10128



Doorlock Harness

OEM:349784T

P/N:10031

J2-21 Series



Motor Exchange

Rebuild Price Available!

OEM:341742

P/N:OCM-02P



Drive Belt

3/16" thick

OEM:878821

P/N:10098



Spindle

OEM:364934

P/N:20502



Door Solenoid

OEM:362771

P/N:10147



Brush, Pair

with Indicator Wire

Thick Brush

OEM:364934

P/N:CB-041



Brush, Pair

with Indicator Wire

Thin Brush

OEM:349454

P/N:CB-144

Beckman Centrifuges

TJ-6 Series



Meter RPM
OEM:339209
P/N:10006



Striker
OEM:339214
P/N:10007



Motor Bellow
OEM:339224
P/N:10008



Imb. / Lid Switch
OEM:340874
P/N:10010




Brake/Speed Knob
Push-on Type
OEM:341298
P/N:10011



Tach Pick-Up
OEM:200-24D
P/N:10014



Speed Control
OEM:341757
P/N:10015



Spindle pulley
OEM:341814
P/N:10017



Cushions (24pk)
For 339282 Carrier
OEM:342602
P/N:10018



Cushion (18pk)
For 339279 Carrier
OEM:342603
P/N:10019



Diode
OEM:343159
P/N:10020



Latch Knob
OEM:878519
P/N:10038



TJ 6 Timer Kit
Timer Kit
OEM:362276
P/N:10063



TH-4 Rotor Insert
OEM:399266
P/N:10087



Brake Resistor
2 Ohm 100w
OEM:878534
P/N:10089



Stop LED
OEM:878583
P/N:10097



Rubber Foot
OEM:358798
P/N:10100



Start Switch
OEM:878609
P/N:10092



Rubber Shock Mount
OEM:878520
P/N:10094



Bridge Rectifier
OEM:878526
P/N:10095



Power Switch
OEM:878629
P/N:10096



Fan Blade
OEM:885023
P/N:10104



Relay
OEM:878815
P/N:10106



Drive Belt
OEM:878513
P/N:10107

Beckman Centrifuges

TJ-6 Series



Fuse Holder
OEM:BFH-1
P/N:10118



Hinge Kit
OEM:360744
P/N:10173



Zero Speed Indicator PCB
Sold Exchange
OEM:340858
P/N:10156



Timer Overlay
OEM:61855-0
P/N:10157



Bowl Gasket
OEM:339228
P/N:10185



Carter Motor Exchange
OEM:345372
P/N:OCM-02B



Brush Pair
For Elinco Motor
OEM:349544
P/N:CB-021



Brush Pair
For Carter Motor
OEM:340536
P/N:CB-020

Microfuge 22R & R



New Motor Exchange
OEM:365812
P/N:30159



Timer 30 min.
OEM:927457
P/N:10088



Brush/Pair
Short Brush
OEM:356883
P/N:CB-096



Brush/Pair
Long Brush
OEM:356881
P/N:CB-100

Ozark Biomedical strives to provide high-quality parts and services to our customers...and others have noticed.

Ozark Biomedical is a registered ISO9001:2008 Company!

Beckman Centrifuges

Rotor Inserts



Aerosolve Tube Racks

Set of 2

Includes Cushions.
24 Place, 14x110mm Tubes
OEM:359483

P/N:10175



Modular-Disk Adapter

Set of 2

12ml 16mm x 120mm
OEM:359471

P/N:10164



Modular-Disk Adapter

Set of 2

10ml 14mm x 110mm
OEM:359470

P/N:10162



Modular-Disk Adapter

Set of 2

Includes Cushions.
30 Place, 13x110mm Tubes
OEM:359478

P/N:10184



JS4.2 SM Blood Bag Cup

Set of 2

OEM:356993

P/N:10169



PTS-2000 Cannister Cover

Set of 2

OEM:362160

P/N:10172



Modular-Disk Adapter

Set of 4

10ml 14mm x 110mm
OEM:359470

P/N:10163



Modular-Disk Adapter

Set of 4

12ml 16mm x 120mm
OEM:359471

P/N:10165



Bail For Modular-Disk Adapters

OEM:349817

P/N:10183

Beckman Centrifuges

Spinchron DLX



Door Seal
OEM:365439
P/N:10002



Power PCB
OEM:392505
P/N:10003



NEW Power PCB Exchange
OEM:392505
P/N:10161



Display PCB Exchange
OEM:365776
P/N:10069



Motor Boot
OEM:358009
P/N:10044



Striker
OEM:358020
P/N:10045



Display Cable
OEM:365404
P/N:10067



Front Panel Spacer
OEM:392036
P/N:10041



Tach Pick-up
OEM:369514
P/N:10074



Tach Pick-up
OEM:369515
P/N:10073



Cooling Fan
OEM:961189
P/N:10090



GFCI
OEM:970423
P/N:10108



Brake Resistor
100 ohm 50 watt
OEM:974782
P/N:10121



Vibration Damper Kit
OEM:367899
P/N:10153



Grommet
OEM:349786
P/N:10160



Ring mount
6.5"x.75"
OEM:367218
P/N:10042



Imb. Switch
OEM:927901
P/N:10111



Touch Pad w/ PCB
OEM:365776
P/N:10070



Touch Pad Exchange
OEM:365735
P/N:10068



Circuit Breaker 15A
OEM:893388
P/N:10177



Safety Lock Exchange
OEM:365408
P/N:20119



New Latch Exchange
OEM:365706
P/N:20122



Latch Exchange
OEM:365706
P/N:20123



Motor Exchange
Rebuild Price Available!
OEM:365434
P/N:OCM-02Q



Rubber Foot
OEM:358798
P/N: 10100

Beckman Centrifuges

Spinchron 15 / Allegra 21



Display PCB Exchange

For 220V Non-Refrig. Units

OEM:329563

P/N:10001



Display PCB Exchange

For 120V Non-Refrig. Units

OEM:365745

P/N:10084



Power PCB Exchange

For 120V Units

OEM:365738

P/N:10129



Power PCB Exchange

For 220V Units

OEM:361357

P/N:10071



Lid Spring

Left 2.26"

OEM:363764

P/N:10055



Lid Spring

Right 2.26"

OEM:363765

P/N:10056



Tapered Sleeve

OEM:361328

P/N:10051



Motor Capseal

OEM:361330

P/N:10053



Touch Pad Exchange

OEM:365735

P/N:10068



Rotor ID Pickup

Fits the 120V Units

OEM:365752

P/N:10077



Tach Cable

Fits the 120V Units

OEM:365754

P/N:10078



Motor Exchange

Rebuild Price Available!

OEM:361348

P/N:OCM-021



115V Chassis Fan

OEM:961189

P/N:10090



220V Chassis Fan

OEM:961838

P/N:10103



120V Power Inlet

OEM:964643

P/N:10112



Tach PCB Exchange

Rectangular Style

OEM:361349

P/N:10049



Motor Mount Each

OEM:961803

P/N:10119



Rotor Tie Down

OEM:361367

P/N:10135



Motor Boot

OEM:364206

P/N:10064



New 120V Latch Exchange

OEM:365706

P/N:20122



120V Latch Exchange

OEM:365706

P/N:20123



220V Latch Exchange

OEM:961845

P/N:20127

Beckman Centrifuges

Spinchron



PCB Assy Exchange
OEM:358946
P/N:10023



Bowl Gasket
OEM:349745
P/N:10025



Vibration Mount Spring
OEM:349756
P/N:10026



Power Switch
OEM:927904
P/N:10123



Knob for Speed/Temp/Timer
OEM:974672
P/N:10030



Thermistor Assy
OEM:349787
P/N:10033



Motor Mount Loop
OEM:355821
P/N:10036



Hinge Kit
OEM:360744
P/N:10173



Brake Switch
OEM:927462
P/N:10039




Display PCB Exchange
OEM:357190
P/N:10043



Motor Boot
OEM:358009
P/N:10044



Striker
OEM:358020
P/N:10045



Door Molded Window
OEM:358753
P/N:10046



Vibration Mount Spacer
OEM:362274
P/N:10050



Doorlock Harness
OEM:359328
P/N:10052



Tach Sensor
OEM:358949
P/N:10057



Speed Limiter Kit
OEM:364239
P/N:10065



Rotor Nut For GH-3.8
OEM:366642
P/N:10081



Vibration Mount Hardware
OEM:N/A
P/N:10085



Speed/Temp
OEM:961314
P/N:10122



Timer 30 minute
OEM:927457
P/N:10088



115V Chassis Fan
OEM:961189
P/N:10090



Lid Open Switch
OEM:927461
P/N:10110



Imbalance Switch
OEM:927901
P/N:10111



Brake Resistor
13 Ohm 100W
OEM:927526
P/N:10114



Bridge Rectifier
OEM:927473
P/N:10115



Receptacle power
OEM:927867
P/N:10116




Safety Lock Exchange
OEM:363414
P/N:20118

Beckman Centrifuges

Spinchron

Beckman



Roll Pin
OEM:964556
P/N:10128



Cooling Fan Cover
OEM:349729
P/N:10132



Lid Door Filter
OEM:361996
P/N:10138



O-Ring Rotor Nut
OEM:927600
P/N:10182



Vibration Mount Kit
OEM:N/A
P/N:10146



Knob for Timer
OEM:974672
P/N:10150



Motor Grommet
OEM:349786
P/N:10160



Motor Exchange
Rebuild Price Available!
OEM:359325
P/N:OCM-02G



Brush, Pair
Autostop brush
OEM:364241
P/N:CB-155



Brush/Pair Kit
Contains grommets, brush caps,
locknuts, 2 pair of CB-155
OEM:364241
P/N:CB-155A



Brush, Pair
OEM:356527
P/N:CB-081



Brush/Pair Kit
Contains grommets, brush caps,
locknuts, 2 pair of CB-081
OEM:356527
P/N:CB-081



Latch Assy. Exchange
For Non-Refrig. Units
OEM:360549
P/N:20120



Latch Assy. Exchange
For Refrig. Units
OEM:360549
P/N:20121

Becton Dickinson/Clay Adams

QBC 4207



Brush, Pair
For S/N 184999 & Below
OEM:42445709
P/N:CB-083



Brush, Pair
For S/N 185000 & Above
OEM:42445713
P/N:CB-085



Motor Rebuild
Must Call for RGA#
OEM:42445502
P/N:OCM-03L



Cover Knob
OEM:465829
P/N:465829

Agglut Viewer



Bulbs, 2-Pk
12v 5w
OEM:421484
P/N:421484

Marsters Incubator



Power Switch
OEM:421814
P/N:421814



Lamp Assy
For Model 5380
OEM:421817
P/N:421817



Thermometer
Right Angle
OEM:421818
P/N:421818




Pilot Light
OEM:421820
P/N:421820



Heating Element
OEM:421821
P/N:421821




Relay
OEM:421823
P/N:421823



Feet, 4-pk
OEM:421825
P/N:421825



Inserts, 5-pk
OEM:421826
P/N:421826



Temp-Trim Top Plate
OEM:421828
P/N:421828



Temp Trim PC Board
OEM:421829
P/N:421829N



Temp Trim Thermistor
OEM:421830
P/N:421830



Temp Trim Pot
OEM:421831
P/N:421831

Becton Dickinson/Clay Adams

Becton Dickinson/
Clay Adams

Nutator Mixer



New Motor
OEM:421504
P/N:30169



On/Off Switch
OEM:421505
P/N:421505

Safetyhead



Brush Pair
For Jones & Black Universal Motors
OEM:42055102
P/N:CB-001

Sero-Fuge



12 Place Rotor
OEM:420545
P/N:420545



Timer
OEM:421365
P/N:421365T



Hold Switch
OEM:421366
P/N:421366T



New Spindle Exchange
OEM:421367
P/N:421367



Spindle Exchange
OEM:421367
P/N:421367R



Rubber Channel
OEM:421368
P/N:421368T



Rubber Strip Liner
OEM:421371
P/N:421371



Rubber Feet 4-pk
OEM:421388
P/N:421388T



Motor Exchange
Rebuild Price Available!
OEM:421375
P/N:OCM-03G

Becton Dickinson/Clay Adams

Sero-Fuge II



12 Place Rotor
OEM:420545
P/N:420545



Timer with Knob
OEM:421379
P/N:421379T



Spindle Exchange
OEM:421367
P/N:421367R



New Spindle Exchange
OEM:421367
P/N:421367



Safety Switch
OEM:421382
P/N:421382T



Rubber Channel
OEM:421368
P/N:421368T



Rubber Stip Liner
OEM:421371
P/N:421371



Rubber Feet 4-pk
OEM:421388
P/N:421388T



Motor Exchange
Rebuild Price Available!
OEM:421375
P/N:OCM-03G



New Motor Exchange
OEM:421375
P/N:30127



Catch Block
OEM:421395
P/N:421395



Latch & Hardware
New Style
OEM:421459
P/N:421459



Brush Assy
OEM:421376
P/N:421376



Governor Assy.
OEM:421377
P/N:421377



Governor Brush, Pair
OEM:42054120
P/N:CB-070



Timer Knob
OEM:421380
P/N:421380




Hi/Lo Switch
OEM:421391
P/N:421391



Base Seal
OEM:421393
P/N:421393T



Cover Assy New style
OEM:421385
P/N:421385T



Spacer Pk of 4
OEM:421394
P/N:421394

Becton Dickinson/Clay Adams

Compact / Compact II

Becton Dickinson/
Clay Adams



Motor Exchange
Rebuild Price Available!
Fits Compact / Analytical
OEM:42015101
P/N:OCM-03P



Motor Exchange
Rebuild Price Available!
Fits Compact II
OEM:421319
P/N:OCM-03T



New Motor Exchange
Fits Compact II
OEM:421319
P/N:30118



6-place Rotor Head
OEM:421320
P/N:421320



Rubber Foot
OEM:421136
P/N:421136T



Drive for Motor
OEM:421260
P/N:421260



Head Screw Assy
OEM:421261
P/N:421261



30 Min. Timer Kit
OEM:421323
P/N:421323T



Capacitor
6 MFD
OEM:421324
P/N:421324



Seal Ring
OEM:421321
P/N:421321



Knob
OEM:421322
P/N:421322



Cover Assy
OEM:421333
P/N:421333



Safety Switch
OEM:421326
P/N:421326



Cover Latch Assy
OEM:421329
P/N:421329

Becton Dickinson/Clay Adams

Sero-Fuge 2000 Series



Latch Exchange
OEM:421359
P/N:20132



New Latch Exchange
OEM:421359
P/N:20141



New Motor Exchange
OEM:421351
P/N:30125



Motor Exchange
Rebuild Price Available!
OEM:421351
P/N:OCM-03R




Fuse
4 amp time delay
OEM:421388
P/N:421338



Latch Hook Kit
OEM:421340
P/N:421340



Power Entry Module
OEM:421341
P/N:421341



Lid Seal Gasket
OEM:421343
P/N:421343



New PCB Exchange
OEM:421344
P/N:421344



PCB Exchange
OEM:421344
P/N:421344R



Keypad
For Sero-Fuge 2001
OEM:421345
P/N:421345



Keypad
For Sero-Fuge 2002
OEM:421354
P/N:421354



Lid Switch Cable
OEM:321358
P/N:421358



Hinge Set Kit
OEM:422953
P/N:422953



Lidseal Retrofit Kit
OEM:429614
P/N:429614



Rubber Foot
OEM:421136
P/N:421136T

Becton Dickinson/Clay Adams

Becton Dickinson/
Clay Adams

Triac



New Motor Exchange
With Gear
S/N 145000 & Above
OEM:421282
P/N:30122




NEW Motor Exchange
220V Motor Without Gear
OEM:420601
P/N:30123



New Motor Exchange
Without Gear
S/N 145000 & Below
OEM:421282
P/N:30132



Brush, Pair
OEM:4202002
P/N:CB-003



Head Assy. w/ Cover
OEM:420202
P/N:420202



Tube Adapter
2ml 4-pack
OEM:420203
P/N:420203



Tube Adapter
3ml 4-pack
OEM:420204
P/N:420204



Motor Gasket
OEM:421264
P/N:421264



Timer
OEM:421265
P/N:421265T



NEW Lid Assembly
OEM:421268
P/N:421268



Cover Gasket
OEM:421269
P/N:421269



Mode Selector
OEM:421270
P/N:421270



Head Fastner Nut
OEM:421272
P/N:421272



Lower Head Cover
OEM:421273
P/N:421273



Wrench Spanner
OEM:421275
P/N:421275



Trunnion Cushion
(pk of 8)
OEM:421276
P/N:421276



Tube Cushion
(pk of 8)
OEM:421278
P/N:421278



Head Washer
OEM:421279
P/N:421279



Head Cover w/ Knob
OEM:421280
P/N:421280



Safety Switch
For S/N 145000 & Above
OEM:421283
P/N:421283T



Tach Sensor
OEM:421284
P/N:421284



Rubber Feet
(pk of 4)
OEM:421285
P/N:421285T




Transformer
OEM:421286
P/N:421286



PCB
For S/N 145000 & Above
OEM:421287
P/N:421287

Becton Dickinson/Clay Adams

MHCT / Autocrit



Head Adapter Nut
OEM:30894
P/N:38333



Safety Locknut
OEM:421403
P/N:421403



Timer & Switch
OEM:421404
P/N:421404



Timer Knob
OEM:421405
P/N:421405



Safety Switch
OEM:421406
P/N:421406



Latch Button Hardware
OEM:421408
P/N:421408



Head Cover
OEM:421410
P/N:421410



Head Liner Rubber
(pk of 3)
OEM:421411
P/N:421411



Shock Mounts Kit
OEM:421412
P/N:421412



Rubber Feet
(pk of 4)
OEM:421413
P/N:421413



Spanner Wrench
OEM:421414
P/N:421414



Lid Safety Switch
OEM:421422
P/N:421422T



Scale Plate
OEM:421444
P/N:421444



New Motor Exchange
OEM:421401
P/N:30177



Brush, Pair
For Black Universal Motor
OEM:42055102
P/N:CB-001



Brush, Pair
For Jones & Specialty Motor
OEM:42055813
P/N:CB-001B

Becton Dickinson/Clay Adams

MHCT II / Autocrit Ultra 3



New Motor Exchange
OEM:421401
P/N:30145



Brush, Pair
For Jones Motor
OEM:421429
P/N:421429



Brush, Pair
Fits RPI Motor ONLY
OEM:NA
P/N:CB-001A



Brush, Pair
For Jones & Specialty Motor
OEM:42055813
P/N:CB-001B



Notched Rubber Liner
OEM:421447
P/N:421447



Lid Safety Switch
OEM:421422
P/N:421422T



Latch Solenoid Assy.
OEM:421424
P/N:421424



Lid Gasket
OEM:421425
P/N:421425



Lamp Assy
OEM:421427
P/N:421427



Cover Spring Assy.
OEM:421428
P/N:421428



Head Adapter Nut
OEM:30894
P/N:38333



Magnifier
OEM:421431
P/N:421431



Head Cover
OEM:421438
P/N:421438



Rubber Foot
OEM:421136
P/N:421136T



Guide Pin
OEM:421446
P/N:421446

Dynac



New Motor Exchange
OEM:421192
P/N:30165



Hemogard Spacer .75"
(pk of 4)
OEM:420252
P/N:420252



Cushions-green
(pk of 12)
OEM:420943
P/N:420943



0901 Rubber Adapters
(pk of 12)
OEM:420944
P/N:420944



Timer & Knob
OEM:421194
P/N:421194



Timer Knob
OEM:421195
P/N:421195



Switch Brake
OEM:421196
P/N:421196



Latch Set
OEM:421198
P/N:421198

Becton Dickinson/Clay Adams

Dynac



Speed Control
OEM:421199
P/N:421199



Speed Control Knob
OEM:421200
P/N:421200



Lamp Assy
OEM:421427
P/N:421427



Head Screw & Retainer
OEM:421202
P/N:421202



Rubber Bumper
(pk Of 3)
OEM:421203
P/N:421203



Resistor
15 ohm
OEM:421205
P/N:421205U



Brush, Pair.
Fits Jones Motor
OEM:42000102
P/N:CB-002



Motor Exchange
Rebuild Price Available!
OEM:421192
P/N:OCM-03B

Dynac II / III



New High Torque Motor
Exchange For Dynac II
OEM:421242
P/N:30124



Low Torque Motor
Exchange For Dynac II
Rebuild Price Available!
OEM:421206
P/N:OCM-03C



New High Torque Motor
Exchange for Dynac III
OEM:421242
P/N:30168



Brush, Pair
Fits Jones Motor
OEM:42055813
P/N:CB-002



Cover w/ Hinge Assy
OEM:421216
P/N:421216T



Fuse Holder
OEM:421229
P/N:421229T



Interrupter Assy
OEM:421234
P/N:421234



Latch Assy
OEM:421217
P/N:421217



Vibration Isolator
(pk of 3)
OEM:421226
P/N:421226



0902 Shields
(pk of 6)
OEM:420902
P/N:420902U



**Rubber Cushion &
Leather Disc for 50ml Shields**
OEM:420941
P/N:420941



Safety Switch
OEM:421222
P/N:421222T

Becton Dickinson/Clay Adams

Dynac II / III

Becton Dickinson/
Clay Adams



Motor Gasket
OEM:421210
P/N:421210



Timer Assy
OEM:421212
P/N:421212



Power Switch
OEM:421214
P/N:421214T



Brake Switch
OEM:421215
P/N:421215



Solenoid Cover
OEM:421223
P/N:421223



Brake Resistor
25 Ohm/30W
OEM:421224
P/N:421224



New Tachometer
OEM:421225
P/N:421225



New Speed Control
Sold Exchange
OEM:421219
P/N:421219



Transformer
OEM:421232
P/N:421232



K2 Relay
OEM:421233
P/N:421233



Brake Resistor
For Dynac II, 50 ohms
OEM:421230
P/N:421230U



PC Board Assy Exchange
For Dynac II
OEM:421231
P/N:421231R



New Tach PCB Exchange
For Dynac II
OEM:421236
P/N:421236



New Tach PCB Exchange
For Dynac III
OEM:421243
P/N:421243



Front Panel Overlay
For Dynac II
OEM:421209
P/N:421209



New Front Label
For Dynac III
OEM:421251
P/N:421251T



Solenoid Assy
OEM:421248
P/N:421248



Safety Switch Assy
OEM:421254
P/N:421254



NEW Rev. B Main PCB
For Dynac III, Exchange
OEM:421250
P/N:421250



Rev. B Main PCB
For Dynac III, Exchange
OEM:421250
P/N:421250R



Brake Resistor
For Dynac III
OEM:421252
P/N:421252

Becton Dickinson/Clay Adams

Readacrit



New Motor
OEM:421449
P/N:30131



Brush, Pair
For Black or Silver Motor
w/small Brush Caps
OEM:42059102
P/N:CB-005



Brush, Pair
For Silver Motor
w/ Large Brush Cap
OEM:42059126
P/N:CB-089



Brush Kit
OEM:421450
P/N:421450



Safety Switch
New Style
OEM:421454
P/N:421454T



Latch & Hardware
New Style
OEM:421459
P/N:421459



Ball Catch
OEM:421460
P/N:421460



New Gulwing Assy.
OEM:421461
P/N:421461T



Head Cover Hardware
OEM:421463
P/N:421463



New Head Assy.
OEM:421464
P/N:421464



Rubber Foot
OEM:421136
P/N:421136T



Catch block
OEM:421473
P/N:421473



Rubber Channel
OEM:421466
P/N:421466



Rotor Cushions
(pk of 12)
OEM:421471
P/N:421471



Wrench
OEM:421472
P/N:421472



Timer Knob
OEM:421480
P/N:421480



Timer Conversion Kit
OEM:421475
P/N:421475



Lid Cover Assy
OEM:421476
P/N:421476



Timer
OEM:421479
P/N:421479T



Sedi-Cal Slide
OEM:421483
P/N:421483



Balance Tube
OEM:421482
P/N:421482



Safety Switch
Old-Style
OEM:42059105
P/N:42059105

Biodynamics

Select-a-Fuge 630



Brush, Pair
S/N 70999 & Below
OEM:919807
P/N:CB-026



Brush, Pair
S/N 71000 & Above
OEM:919808
P/N:CB-027



Motor Exchange
Rebuild Price Available!
OEM:1-61703
P/N:OCM-12A

Eppendorf /Brinkmann

5402



Brush, Pair
OEM:2236527-4
P/N:CB-078



Motor Exchange
Rebuild Price Available!
OEM:22-37519-9
P/N:OCM-11F



Timer
OEM:927457
P/N:10088



Brush, Pair
OEM:2236521-5
P/N:CB-075



Motor Exchange
Rebuild Price Available!
OEM:22-37-240-8
P/N:OCM-11A

5413



Brush, Pair
OEM:2236521-5
P/N:CB-075



Motor Exchange
Rebuild Price Available!
OEM:22-36-511-8
P/N:OCM-11B



Motor Exchange
Rebuild Price Available!
OEM:NA
P/N:OCM-11L

5804 / 5810

Eppendorf / Brinkmann

5415



Lock Microswitch
OEM:2237096-1
P/N:2237096-1



Microswitch
OEM:2237161-4
P/N:2237161-4



Power Switch
OEM:2237216-5
P/N:2237216-5



Lid Lock Microswitch
OEM:2237508-3
P/N:2237508-3



Brush, Pair
For S/N 3716 & Below
OEM:2236524-0
P/N:CB-076



Brush, Pair
For S/N 3717 & Above
OEM:2236523-1
P/N:CB-077



Motor Exchange
Bottom Tach
Rebuild Price Available!
OEM:2237245-9
P/N:OCM-11C



Motor Mount
OEM:2237250-5
P/N:2237250-5



Opti-Coupler
OEM:2237252-1
P/N:2237252-1

5415C



Code Disk
OEM:2237714-1
P/N:2237714-1



Opti-Coupler
OEM:2237716-1
P/N:2237716-1



Fan
OEM:22377213
P/N:22377213



Side Tach
OEM:2237767-1
P/N:2237767-1



New Motor Exchange
w/Extended Tach
OEM:022375181
P/N:30156



Motor Exchange
Side-Tach, 4 Wires
For S/N 25377 & Up
OEM:022375181
P/N:OCM-11D



Motor Exchange
Bottom Tach, 4 Wires
S/N 25376 & Down
OEM:02237710-8
P/N:OCM-11E



Brush Pair
OEM:2236527-4
P/N:CB-078

Baxter / Dade

Immufuge / Immufuge II



Motor Rebuild
Must call for RGA #
OEM:NA
P/N:OCM-05



Motor Exchange
Rebuild Price Available!
OEM:043020-100
P/N:OCM-05A



New 12-Place Rotor
Sold Exchange
OEM:B5055-2
P/N:38522



12-Place Rotor
Sold Exchange
OEM:B5055-2
P/N:38219



Latch Exchange
OEM:S45010-001
P/N:20102



Brake Stop Plastic
OEM:020018-001
P/N:38189



Brake Lever Assy
w/ Brake
OEM:020067-001
P/N:38190



Solenoid Brake
OEM:020072-100
P/N:38191



Lid Spring
(pair)
OEM:043022-100
P/N:38192



Lid Catch
OEM:043094-001
P/N:38194



Lid Switch & Boot
OEM:S43092-001
P/N:38195



Grommet, Solenoid
OEM:13029
P/N:38197



Solid State Relay
OEM:NA
P/N:38224



SSR PCB Exchange
For S/N 64799 and Below
OEM:S43069-001R
P/N:38225



Rotor Hub Pin
OEM:S14194-100
P/N:38230



Panel Label Assy
OEM:S43004-001
P/N:38236



Spindle Assy Metal
OEM:S45013-001
P/N:38238



5 Volt Power Supply
OEM:S45062-001
P/N:38239



Spring Clip Kit
OEM:B5055-9
P/N:38256



SSR PCB Exchange
For S/N 64800 & Above
OEM:S43069-001
P/N:38500



On/Off Switch
OEM:683S0103
P/N:38216



O-Ring Each
OEM:NA
P/N:38535



Rubber Foot
OEM:421136
P/N:421136T



Main PCB Exchange
OEM:S43064-001
P/N:38223

Baxter / Dade

DAC II



Motor Exchange
OEM:S23283-100
P/N:OCM-05D



Magnet Brush, Pair
OEM:S45099-001
P/N:CB-147



Motor Brush, Pair
OEM:S45098-001
P/N:CB-148



Brush Holder Kit
OEM:S45165-001
P/N:38214



Decant Limit Pot
OEM:18548
P/N:38199



Spin Speed Limit Pot
OEM:18549
P/N:38200



Relay K1-K2
OEM:27536
P/N:38202



60V SSR for SSR2
OEM:27549
P/N:38204



Latch Upgrade Kit
OEM:NA
P/N:38205



Saline Flow Switch
OEM:28092
P/N:38206



Fuse Holder
OEM:28527
P/N:38207



Saline Detector Body
OEM:NA
P/N:38208



OEM Pump Tubing
OEM:6402-18
P/N:38209



Pump Tubing
OEM:B5054-23
P/N:38228




Pump Head Assy
OEM:S23050-100
P/N:38232



Tube Holder Assy
OEM:B5054-10
P/N:38210




Transformer T1
OEM:S29519
P/N:38235



Tubing Kit
OEM:300181JS
P/N:38217



Spring For Lid
OEM:S45041-100
P/N:38222



Fill Tube
OEM:B5054-20
P/N:38226



Main PCB Exchange
OEM:S45122-001
P/N:38211



Darlington Transistor
OEM:07F8524
P/N:38213

Drucker

611T / 613 / 614



New Motor
Fits 611T, 613B, 614B/V
OEM:7735050
P/N:30101

708E / 708T



Motor Exchange
Fits 708E/EL/T
Rebuild Price Available!
OEM:7767153
P/N:OCM-06H



Brush, Pair
Fits 708E
OEM:42000102
P/N:CB-002

755V



New Motor Exchange
Fits 755V, 853VES, &
Plasmafuge 6/12
OEM:7735016
P/N:30183



Motor Exchange
Fits 755V, 853VES, &
Plasmafuge 6/12
Rebuild Price Available!
OEM:7735016
P/N:OCM-06I

809T



Brush, Pair
Fits 803, 806, 809,
809T, 809R, 818R
OEM:7735003
P/N:CB-030

818R



Brush, Pair
Fits 803, 806, 809,
809T, 809R, 818R
OEM:7735003
P/N:CB-030

L-708




Motor Rebuild
Must Call for RGA#
OEM:7767053
P/N:OCM-06A



Brush, Pair
Fits L708, L708E,
411, 720E
OEM:7735023
P/N:CB-018

Fisher


Accuspin 1



Lid Hook
OEM:20052842
P/N:38137



Latch Assy Exchange
OEM:50079829
P/N:20111



Main PCB Exchange
Fits Accuspin 1/R
OEM:20151092
P/N:38528



Gas Lid Stay
For Non-Refrig. Units
OEM:20052842
P/N:38137



Gas Lid Stay
For Refrig. Units
OEM:20052842
P/N:38158



Motor Exchange
Rebuild Price Available!
OEM:20210332
P/N:OCM-01N

Accuspin 400



Latch Assy. Exchange
OEM:70908552
P/N:20108



Solenoid coil
OEM:5001004
P/N:20112



Lid Hook
OEM:20056527
P/N:38524



Motor Exchange
Rebuild Price Available!
OEM:20210385
P/N:OCM-07R

Marathon 3000



New Motor Exchange
OEM:51035OR
P/N:30151



Motor Exchange
Rebuild Price Available!
OEM:51035OR
P/N:OCM-07Q



Power PCB Exchange
OEM:51403OR-R
P/N:38309



Motor Mount
OEM:51053CR
P/N:38429



Crowned Shaft adapter
Fits 243 Rotor
OEM:50969OF
P/N:38424



Collar
OEM:486210R
P/N:38506



Gas Shock
Fits Marathon 3000R
OEM:51062OR
P/N:38405

Fisher

Centrifric 225



Motor Exchange
Rebuild Price Available!
OEM:SPN38218
P/N:OCM-07B



Brush, Pair
OEM:SPN39329
P/N:CB-106



Thumb Screw
OEM:38816
P/N:38491



Timer, 30 Min
OEM:SPN107126R
P/N:38489



Timer, 60 Min
OEM:SPN107126R
P/N:38490



Speed Control Assy.
Sold Exchange
OEM:SPN88255
P/N:38492

Centrifric 228



Motor Exchange
Rebuild Price Available!
OEM:50782
P/N:OCM-07J



Timer Knob
OEM:106869
P/N:38583



15-Min. Timer
OEM:106868
P/N:38584



60-Min. Timer
OEM:70867
P/N:38488

13K



Motor Exchange
Rebuild Price Available!
OEM:HL259007302
P/N:OCM-07F



Brush, Pair
OEM:SPN83142
P/N:CB-119

Fisher

Marathon 3200



Motor Exchange
Rebuild Price Available!
 OEM:50900OF
P/N:OCM-08MM



Power PCB Exchange
 OEM:51403OR-R
P/N:38309



Crowned Shaft adapter
 Fits 243 Rotor
 OEM:50969OF
P/N:38424



Collar
 OEM:486210R
P/N:38506



Gas Shock
 Fits Marathon 3200
 OEM:50907OR
P/N:38484



Motor Boot
 OEM:50902AR
P/N:38403



Latch Assy. Exchange
 OEM:50929R
P/N:20134



Motor Mount
 OEM:50911OR
P/N:38439



Fan
 OEM:50915OR
P/N:38455

235A / 235B / 235C

335



Motor Exchange
Rebuild Price Available!
 OEM:SPN70743
P/N:OCM-07A



Brush, Pair
 OEM:43425
P/N:CB-039



Motor Exchange
Rebuild Price Available!
 OEM:NA
P/N:OCM-07I



Brush, Pair
 OEM:43425
P/N:CB-039

59A

6K



Motor Exchange
Rebuild Price Available!
 OEM:NA
P/N:OCM-07H



Brush, Pair
 OEM:SPN06649
P/N:CB-038



Motor Exchange
Rebuild Price Available!
 OEM:HL243008103
P/N:OCM-07E



Brush, Pair
 OEM:SPN83788
P/N:CB-121



Ozark Biomedical is built around the needs of our customers. We offer new and rebuilt motors, latches, and PCB's for many of today's laboratory centrifuges. In many cases, you will save nearly 50% on parts versus buying from the OEM. All while still meeting strict quality demands of ISO 9001:2008. That is the savings you want while ensuring the quality you deserve.

Visit www.ozarkbiomedical.com for all your centrifuge parts & service needs.

Ozark Biomedical also offers quality New, Reconditioned, and Certified centrifuges at significant savings. All our centrifuges are backed by our industry-recognized support.

We offer new centrifuges from Clay Adams, Drucker, Hettich, Thermo Fisher, to name a few. These models are brand new, direct from OEM with warranty coverage.

Our reconditioned centrifuges are hand-picked, reliable units that are meticulously rebuilt from the chasis up. These units then undergo a stringent Quality Control process that meet or exceed OEM specifications.

Ozark Biomedical's line of recertified centrifuges are gently used units that have been inspected for wear, serviced, and tested to meet OEM standards. These units offer the highest degree of savings for end users. There is plenty of life left in these centrifuges.

Visit www.centrifugestore.com to view our centrifuge lines available for purchase. Use the online contact form to let us assist you in finding a particular centrifuge.



Rupp Bowman / Global Focus

Silencer 2100 / 2200 / 2300



Motor Exchange
Rebuild Price Available!
Fits 2100 Units
OEM:SIL00719
P/N:OCM-18F



Motor Exchange
Rebuild Price Available!
Fits 2200 /2300 Units
OEM:SIL00756
P/N:OCM-18G



Brush, Pair
Fits 2100 / 2200 / 2300
OEM:00033
P/N:CB-159



Cable and Motor J3-J2
OEM:SIL90200
P/N:35022



Latch Exchange
Fits 2100 /2200 / 2300
OEM:SIL00671R
P/N:20135



Magnet
Fits 2200 / 2300
OEM:SIL00682
P/N:35009



Drive/Brake Relay
Fits 2100 / 2200
OEM:SIL01102
P/N:35014



Cable Assy
Fits 2200 Temp & Imbalance
OEM:Roc0001
P/N:35019



Tach Pickup
Fits 2100 / 2200 / 2300
OEM:SIL00838
P/N:35013



New Power PCB
Fits 2100 / 2200 / 2300
OEM:SIL00829
P/N:35012



Display PCB Exchange
Fits 2100
OEM:SIL00807
P/N:35010



Brush Spring
Fits 2200 / 2300
OEM:SIL0095
P/N:35008



Brush Cap
Fits 2200 / 2300
OEM:SIL0092
P/N:35007



Imbalance Switch
Fits 2200
OEM:SIL00818R
P/N:35011



Display PCB Exchange
Fits 2200 / 2300
OEM:SIL00808
P/N:35017

Hamilton Bell

V6000 / 6500

Hamilton Bell



New Motor Exchange
OEM:V7001
P/N:30116



Motor Exchange
OEM:V7001
P/N:OCM-16B



Detent Fastener
OEM:V6055
P/N:41002



Rotor O-Ring
OEM:V6062
P/N:41003



Bowl Gasket
OEM:V6065
P/N:41004



Timer & Knob
OEM:V6070N
P/N:41005



Timer Screw Each
OEM:V6071
P/N:41006



Timer Knob
OEM:V6075
P/N:41007



Capacitor
OEM:V6105
P/N:41008



Cover Assy.
OEM:V6025
P/N:41009



Motor Screw Cap
OEM:V6115
P/N:41010



Power Inlet
OEM:V6125
P/N:41011



Micro Switch
OEM:V6130
P/N:41012



On/Off Switch
OEM:V6140
P/N:41014



On/Off Indicator
OEM:V6145
P/N:41015



Shock Mount Each
OEM:V6150
P/N:41016



Nylon Washer
OEM:V6150
P/N:41017



Rubber Cushions
OEM:V6170
P/N:41018



15ml Shield Each
OEM:V6175
P/N:41019



10ml Shield Each
OEM:V6180
P/N:41020



Shaft Interlock Rod
OEM:V6045
P/N:41021



Interlock Spring
OEM:V6048
P/N:41022



Adapter Cushions
pk of 6
OEM:V6150
P/N:41024



Rubber Foot
OEM:V6135T
P/N:421136T

Hermle

Z230 / Z230A / Z232 / Z232 MK2



Motor Exchange
Fits Z230A / Z232
Rebuild Price Available!
OEM:Co230A-M
P/N:OCM-17K



Motor Exchange
Fits Z230
Rebuild Price Available!
OEM:CO230-15M2
P/N:OCM-17B



Brush Pair
Fits Z230/A & Z232
OEM:Co230-15B
P/N:CB-110



Brush Pair
Fits Z230AMK2
& Z232 MK2
OEM:NA
P/N:CB-166

Z230M / Z230H / Z230HA / Z229 / Z320



Motor Exchange
Fits Z230M
Rebuild Price Available!
OEM:CO230-E
P/N:OCM-17C



Motor Exchange
Fits Z230M / Z230H / Z229
Rebuild Price Available!
OEM:244.00.34.02
P/N:OCM-17C-1



Motor Exchange
Fits Z320
Rebuild Price Available!
OEM:Co320-E
P/N:OCM-17D



Brush Pair
OEM:Co230-G
P/N:CB-109

Z231M



Motor Exchange
Rebuild Price Available!
OEM:Co231M
P/N:OCM-17G



Brush Pair
OEM:Co230-G
P/N:CB-109

Z230MR



Motor Rebuild Only
Call for latest pricing!
OEM:Co230HA-M
P/N:OCM-17J



Brush Pair
OEM:Co230-G
P/N:CB-109

Z300



Motor Exchange
Rebuild Price Available!
OEM:Co300-E
P/N:OCM-17T

ZK / Z364



Motor Exchange
Rebuild Price Available!
OEM:Co300-E
P/N:OCM-17T



Brush Pair
OEM:Co364-G
P/N:CB-104

Z323K



Motor Exchange
Rebuild Price Available!
OEM:CO 323-E
P/N:OCM-17V

Thermo Fisher / IEC

Centra CL2

Thermo Fisher /
IEC



New Armature Exchange
OEM:50594OF
P/N:30108



Armature Exchange
Rebuild Price Available!
Short Shaft
OEM:48411OR
P/N:OCMA-08LL



Armature Exchange
Rebuild Price Available!
Long Shaft
OEM:50594OF
P/N:OCMA-08LL-1



Brush, Pair
New Style w/ Shunt Wire
OEM:1780AF
P/N:CB-006A



Motor Housing Exchange
OEM:42470AF
P/N:38391



Spring Washer
OEM:41778OR
P/N:38348



New Interlock PCB
Sold Exchange
OEM:44566OR
P/N:CB-38371



Clear Window
OEM:48116OR
P/N:38380



Cover Gasket
OEM:T01434
P/N:38538



New Timer/Display PCB
Sold Exchange
Fits S/N 42606000 & Up
OEM:44471OF
P/N:38396



Solenoid Assy.
OEM:50634BR
P/N:38420



Membrane Panel
OEM:50788AR
P/N:38426



Rotor Magnet
OEM:50344AR
P/N:38468



Magnetic Rotor Kit
OEM:50344OF
P/N:38469



Suction Foot Each
OEM:50148OR
P/N:38529

Thermo Fisher / IEC

Centra CL3



Motor Exchange
Rebuild Price Available!
OEM:50900OF
P/N:OCM-08MM



Latch Assy. Exchange
OEM:50594
P/N:20134



Motor Boot
OEM:50594
P/N:38403



Solenoid
OEM:50634CR
P/N:38421



Motor Shaft Seal
OEM:50943OF
P/N:38423



Crown Adapter
For 243 Rotor
OEM:50969OF
P/N:38424



Gas Shock
OEM:50907OR
P/N:38484



Gasket Kit
OEM:50901OF
P/N:38477



Resistor
50 Ohm / 50 Watt
OEM:60430AR
P/N:38476

Centra CL4



Motor Exchange
Rebuild Price Available!
OEM:43310
P/N:OCM-08X



Brush Pair
OEM:44989
P/N:CB-073



Hinge
OEM:49825AR
P/N:38461

Centra 4 / 4R



Motor Rebuild Only
Must Call for RGA #
OEM:49291
P/N:OCM-08



Brush Pair
OEM:49203
P/N:CB-079



Motor Mount Each
OEM:961803
P/N:10119



Latch Exchange
OEM:40469U
P/N:20128

Thermo Fisher / IEC

Centra CL



New Armature Exchange
Include Lower Bearings & Fan
OEM:42469
P/N:30172



New Armature Exchange
No Bearings
OEM:1788
P/N:30173



Brush, Pair
OEM:1780AF
P/N:CB-006A



Plunger
OEM:49658OR
P/N:38460



Rubber Feet
4-pk
OEM:421388
P/N:421388



Switch & Transformer
OEM:36536OR
P/N:38338

Centra 4B



Motor Exchange
Rebuild Price Available!
OEM:49711
P/N:OCM-08H



Brush Pair
OEM:49374
P/N:CB-071



Brush Pair
OEM:44989
P/N:CB-073



Motor Mount Each
OEM:961803
P/N:10119



New Latch
OEM:50405AR
P/N:38416


Centra 8 / 8R



Motor Exchange
Rebuild Price Available!
OEM:49763
P/N:OCM-08A



Brush Pair
OEM:498010F
P/N:CB-072



Switch Button Control
OEM:466538BR
P/N:38356



Hinge, Each
OEM:49825AR
P/N:38461

Thermo Fisher / IEC

Centra 7 / 7R



Motor Exchange
Rebuild Price Available!
OEM:48107
P/N:OCM-08B



Brush Pair
3/16x3/8x7/8
OEM:48192
P/N:CB-035



Brush Pair
1/4x3/8x3/4
OEM:47889
P/N:CB-045



Brush Pair
1/4x7/16x3/4
OEM:47263
P/N:CB-029



Stop Switch
4-prong
OEM:47115AR
P/N:38357



Run Switch
OEM:47116AR
P/N:38358



Power Switch
OEM:47112A
P/N:38414



Motor Boot
OEM:48860
P/N:38375



***Used* Motor Mount**
Pk of 4
OEM:47088
P/N:38413



New Motor Mount
Sold Each
OEM:45924
P/N:38545



Relay
OEM:61081AR
P/N:38440

Centra B / B+



Motor Exchange
Rebuild Price Available!
OEM:44609OR
P/N:OCM-08R



Rubber Foot
OEM:421136
P/N:421136T



Latch Catch
OEM:50405BR
P/N:38417



New Latch
OEM:50405AR
P/N:38416

Thermo Fisher / IEC

Centra W / Cellwasher II Plus



Motor Exchange
Rebuild Price Available!

6-pin GE Motor
w/ Extended Shaft
OEM:49324

P/N:OCM-08V



Motor Exchange
Rebuild Price Available!

5-pin GE Motor
w/ Extended Shaft
OEM:49990

P/N:OCM-08W



Motor Exchange
Rebuild Price Available!

5-pin Ashland Electric Motor
w/ Extended Shaft
OEM:50792A

P/N:OCM-08W-1



Motor Exchange
Rebuild Price Available!

5-pin GE Motor
w/o Clutch Bearing
OEM:49990

P/N:OCM-08Y



Motor Exchange
Rebuild Price Available!

5-pin Ashland Electric or Aeromotors
w/o Clutch Bearing
OEM:50304

P/N:OCM-08Y-1



Latch Assy. Exchange

OEM:42474R

P/N:20133



Guard Bowl

OEM:40255OR

P/N:38303



Decant Plunger

OEM:44919OF-S

P/N:38400



Dowel

OEM:44903OR

P/N:38370



Control PCB Exchange

Eprom Ver. 3.6
OEM:44443OR-R

P/N:38386



Control PCB Exchange

Eprom Ver. 3.8
OEM:44443OR-R

P/N:38511



Power Outlet

w/o Fuse Drawer
OEM:49945OR

P/N:38464



Pump Tubing

OEM:65451OR

P/N:38449



Tubing Kit

OEM:50019OF-S

P/N:38465



Motor Mount

OEM:49354OR

P/N:38459



Fan Blade

OEM:49388OR

P/N:38458



O-Ring Each

OEM:41199OR

P/N:38514



Grommet Each

OEM:49348AR

P/N:38475



Plunger Each

OEM:49348BR

P/N:38406



Brake Resistor

200 Ohm

OEM:60432OR

P/N:38517

Thermo Fisher / IEC

Centra W / Cellwasher II Plus



Pivot Ring
OEM:44906OF
P/N:38369



Rotor Stand
OEM:44945OF
P/N:38411



Tube Holder Pivot
OEM:44921OR
P/N:38402



Tube Holder
OEM:44924OR
P/N:38410

CRU5000 / CU 5000



Motor Exchange
Rebuild Price Available!
S/N 23472584 & below
OEM:45995
P/N:OCM-08G



Motor Exchange
Rebuild Price Available!
S/N 23472584 & above
OEM:65032
P/N:OCM-08G-1



Brush Pair
Fits S/N 23451407 & Down
OEM:48264
P/N:CB-013



Brush Pair
Fits S/N 23451407 & Up
OEM:41775
P/N:CB-080



***Used* Motor Mount**
pk of 4
OEM:47088
P/N:38413

DPR/PR 6000



Motor Exchange
Rebuild Price Available!
Fits PR-6000
OEM:44639
P/N:OCM-08D



Motor Exchange
Rebuild Price Available!
Fits DPR-6000
OEM:46687OF
P/N:OCM-08J



Brush Pair
OEM:48264
P/N:CB-013



Flexible Shaft Exchange
OEM:47051R
P/N:20516

Thermo Fisher / IEC

Equafuge 1L

PR-6 / PR-J



Motor Exchange
Rebuild Price Available!
OEM:50766
P/N:OCM-08EE



Brush Pair
OEM:50757
P/N:CB-163



Motor Exchange
Rebuild Price Available!
OEM:NA
P/N:OCM-08P



Brush Pair
OEM:40069
P/N:CB-011

GP6



New Motor Exchange
OEM:51035OR
P/N:30151



Motor Exchange
Rebuild Price Available!
OEM:51035OR
P/N:OCM-07Q



Solenoid
OEM:50634CR
P/N:38421



Motor Isolator
OEM:51053CR
P/N:38429

UV



Motor Exchange
Rebuild Price Available!
OEM:1939
P/N:OCM-08F



Brush Pair
OEM:1946
P/N:CB-015



Brush, Each
Auto Transformer
OEM:1927
P/N:CB-074

Medispin



New Motor w/ Line Filter
OEM:50339OF
P/N:30106



Upgrade Kit
OEM:50173OF
P/N:38310



Foot (pk of 4)
OEM:50148OR
P/N:38529

Thermo Fisher / IEC

Size 2 Model K



Motor Exchange
Rebuild Price Available!
OEM:1985
P/N:OCM-08E



Brush Pair
OEM:1913
P/N:CB-023



Rubber Tip for Tach
OEM:058830R
P/N:38318



Brush
Auto Transformer
OEM:1785
P/N:CB-146

Spinette



New Motor Exchange
OEM:46620
P/N:30105



Motor Exchange
Rebuild Price Available!
OEM:46620
P/N:OCM-08I



Lid Gasket
OEM:50897OR
P/N:38428

HNS II



Motor Exchange
Rebuild Price Available!
OEM:49955OF
P/N:OCM-08K



New Armature Exchange
Carter Armature w/ Bearings
OEM:43687OR
P/N:30103



Armature Exchange
Rebuild Price Available!
GE Armature
2.165" Diameter (55mm)
OEM:41101
P/N:OCMA-08K



Armature Exchange
Rebuild Price Available!
Carter Armature
2.06" Diameter (52.5mm)
OEM:43687OR
P/N:OCM-08K-1



Brush Pair
OEM:1780AF
P/N:CB-006A



Brush Holders
OEM:09037OF
P/N:38324



Motor Cap
OEM:42249
P/N:38349



Spring Washer
OEM:41778OR
P/N:38348



3-way Switch
OEM:9943
P/N:38326



Brush Holder Spring
OEM:36624
P/N:38328



Armature Key
OEM:32385OR
P/N:38337



Handle Spring
OEM:45489ER
P/N:38362

Thermo Fisher / IEC

HNS II



Leg
OEM:42324OF
P/N:38389



Speed Control Pot.
OEM:61621OR
P/N:38447



Speed Control PCB
Sold Exchange
OEM:45004OR
P/N:38518



Brake Switch
OEM:9942
P/N:421196

GP8 /R



New Motor Exchange
IEC Motor
OEM:65475OF-S
P/N:30102



Motor Exchange
Rebuild Price Available!
IEC Motor
OEM:65475OF-S
P/N:OCM-08BB



Motor Exchange
Rebuild Price Available!
Baldor Motor
OEM:65527OF-S
P/N:OCM-08BB-1



Brush Pair
Fits IEC Motor
OEM:498010F
P/N:CB-072



Brush Pair
Fits Baldor Motor
OEM:T01206S
P/N:CB-304



New Latch Exchange
OEM:43181AR
P/N:20129



Latch Exchange
OEM:43181R
P/N:20130



O-Ring Gasket
Fits Bucket 3218
OEM:478240R
P/N:38376



Window Plug
OEM:48116OR
P/N:38380



Logic PCB Exchange
OEM:44467R
P/N:38388



5-pk of 1A fuses
OEM:846-2007
P/N:38404



Motor Mount
OEM:51053CR
P/N:38429



Brake Resistor
10 Ohm
OEM:60430GR
P/N:38431



Hinge
OEM:65479OR
P/N:38453



Tack Magnet
OEM:65520OR
P/N:38454



Brake/Run Relay
OEM:8860147
P/N:38457

Thermo Fisher / IEC

GP8/8R



Power Outlet
w/o Fuse Drawer
OEM:49945OR
P/N:38464



Membrane Panel
Fits GP8 Only
OEM:65465OR
P/N:38479



Logic PCB Exchange
For GP8R Only
OEM:44465OR
P/N:38485



Circuit Breaker
OEM:62718AR
P/N:38504



Collar for Shaft
OEM:486210R
P/N:38506

MP4 / R



Motor Exchange
Rebuild Price Available!
Program 1.0 Motor
w/ White Brush Caps
S/N 24370600 & below.
Refrig. Units S/N 24380600 & below.
OEM:43310
P/N:OCM-08X



Motor Exchange
Rebuild Price Available!
Program 2.0 Motor
w/ Black Brush Caps
S/N 24371999 & below.
Refrig. Units S/N 24381999 & below.
OEM:44990OF
P/N:OCM-08Z



Motor Exchange
Rebuild Price Available!
Inverted Motor
S/N 24371999 & Up.
Refrig. Units S/N 24381999 & Up.
OEM:50407OR
P/N:OCM-08CC



New Motor Exchange
Inverted Motor
S/N 24371999 & Up.
Refrig. Units S/N 24381999 & Up.
OEM:44996OF
P/N:30115



Brush Pair
Fits OCM-08X
OEM:49374
P/N:CB-071



Brush Pair
Fits OCM-08Z
OEM:44989
P/N:CB-073



Brush Pair
Fits OCM-08CC
OEM:50433
P/N:CB-162



New Latch Exchange
OEM:43181AR
P/N:20129



Latch Exchange
OEM:43181R
P/N:20130



Main PCB Exchange
Fits S/N 24370724 & Above
OEM:44469OR-R
P/N:38470



Main PCB Exchange
Fits S/N 24381500 & Below
OEM:44998OF-R
P/N:38478



Circuit Breaker
OEM:62717B
P/N:38448

Thermo Fisher / IEC

MP4 / R



Motor Mount
OEM:961803
P/N:10119



Tach Sensor Only
OEM:NA
P/N:38314



Shaft Adapter
New Style
OEM:48399OR
P/N:38367



Motor Bellow
S/N 24385000 & Up
OEM:50347OR
P/N:38474



Microswitch
OEM:62716T
P/N:38394



Hinge, Each
OEM:49825AR
P/N:38461



Power Outlet
w/o Fuse Drawer
OEM:499450R
P/N:38464



Fuse Drawer
OEM:43173BR
P/N:38531



New Rotor Nut
OEM:656370F
P/N:38543

Why is ISO9001:2008 Certification so important?

It shows a commitment to Quality Assurance at every level of production. This saves time and money for the company, and it's customers too.

This is why Ozark Biomedical is proud to be a ISO9001:2008 Certified Company.

Thermo Fisher / IEC

MB / MMB



Motor Exchange
MB Motor
OEM:44890OF
P/N:OCM-08C



Motor Exchange
MMB Motor
OEM:44890OF-S
P/N:OCM-08C-1



Brush Pair
For S/N AA2914 & Below
& S/N AC6200 & Above
OEM:1840
P/N:CB-007



Brush Pair
For Jones Motors
S/N AA2914-AC6200
OEM:43653
P/N:CB-007A

Brush Pair
For Specialty Motor
OEM:42055813
P/N:CB-001B

Old-Style Front Panel
OEM:49788T
P/N:38306

New-Style Front Panel
OEM:49788OR
P/N:38307



Timer Kit
OEM:768OF
P/N:38313



Lid Gasket
OEM:01511OR
P/N:38311



Cushion
OEM:01505OR
P/N:38312



Spanner Wrench
OEM:01883OR
P/N:38316



Rotor Cover
OEM:7912OF
P/N:38319



Head Adapter
OEM:7915T
P/N:38320



Head Cover Nut Kit
OEM:7913T
P/N:38322




Mount Spacers
OEM:30882OR
P/N:38329



Head Adapter Nut
OEM:30894
P/N:38333



Latch
OEM:30878OF
P/N:38330



Screen
OEM:30889OR
P/N:38334



Vent Cover Screen
OEM:30900OR
P/N:38335



Motor Mount Kit
OEM:30992
P/N:38336



Spring Plunger
OEM:65455OR
P/N:38450



Rubber Foot
OEM:421136
P/N:421136T



Brake Switch
OEM:9942
P/N:421196



Timer Knob
OEM:974672
P/N:10150

Thermo / Jouan

312 / 412 Series



Motor Exchange

Rebuild Price Available!

312 Series Only
OEM:25509559

P/N:OCM-09C



Motor Exchange

Rebuild Price Available!

412 Series Only
OEM:51245484

P/N:OCM-09H



Brush Pair

6mmx16mmx32mm
Fits OCM-09H / 09I
OEM:89000214

P/N:CB-134



Brush Pair

6mmx10mmx37mm
Fits OCM-09C
OEM:89000056

P/N:CB-135



Control PCB Exchange

312 / 412 Series
OEM:51245483

P/N:36039



New Display PCB

Sold Exchange
312 / 412 Series
OEM:85280072

P/N:36048



Motor Mount

For C412 Series
OEM:26007026

P/N:36012



Relay

312 / 412 Series
OEM:NA

P/N:36080



Timer Assy.

312 / 412 Series
OEM:86000498

P/N:36052



Temp. Sensor

312 / 412 Series
OEM:86001274

P/N:36053



Imbalance Sensor

312 / 412 Series
OEM:85240244


P/N:36070



Rotor Nut Wrench

312 / 412 Series
OEM:25055501


P/N:36087



Lid Cap

312 / 412 Series
OEM:3174805

P/N:36013



Rotor Nut

312 / 412 Series
OEM:25292004

P/N:36004



Lid Seal

For C412 Series
OEM:25437049

P/N:36005



Motor Boot

312 / 412 Series
OEM:25439072

P/N:36008



Door Switch

Lever Style
OEM:3171945

P/N:36009



Lid Lock Grommet

312 / 412 Series
OEM:26554028

P/N:36024

Thermo / Jouan

312 / 412 Series



Knob
312 / 412 Series
OEM:26091006
P/N:36014



Knob Cap
312 / 412 Series
OEM:26136005
P/N:36015



Power Switch
312 / 412 Series
OEM:86001434
P/N:36017




EMI Filter
312 / 412 Series
OEM:26365012
P/N:36020



Motor Grommet
312 / 412 Series
OEM:26554023
P/N:36023



Brake Resistor
312 / 412 Series
OEM:26720062
P/N:36028



Plastic Viewport
312 / 412 Series
OEM:85240275
P/N:36025



Rubber Foot, Each
312 / 412 Series
OEM:26566011
P/N:36026



Speed / Temp Pot.
312 / 412 Series
OEM:26649085
P/N:36027



Gas Shock
412 & 412R Series
OEM:26962008
P/N:36032



Gas Shock
Non-Refrig. 312 Series
OEM:85240275
P/N:36025



Gas Shock
Refrig. 312 Series
OEM:26566011
P/N:36026



Keypad Overlay
Non-Refrig. 312 / 412 Series
OEM:85210059
P/N:36043



Latch Solenoid
312 / 412 Series
OEM:85130085
P/N:36041



Striker
312 / 412 Series
OEM:85240195
P/N:36044

422 Series



Bridge Rectifier
422 Series
OEM:3172591
P/N:36029



Eprom Version J
C422 Series
OEM:85290016
P/N:36050



Relay
422 Series
OEM:86000472
P/N:36051



Power Rack Exchange
422 Series
OEM:25153152R
P/N:36066

Thermo / Jouan

422 Series



Motor Exchange

Rebuild Price Available!

422 Series

OEM:25509567

P/N:OCM-09I



Latch Solenoid

422 Series

OEM:85130085

P/N:36041



Brush Pair

6mmx16mmx32mm

Fits OCM-09I

OEM:89000214

P/N:CB-134



Striker

422 Series

OEM:85240195

P/N:36044



Brush Pair

6mmx10mmx37mm

For Older, Narrow Brush Motor

OEM:89000056

P/N:CB-135




Power Switch

422 Series

OEM:86001434

P/N:36017



Plastic Viewport

422 Series

OEM:85240275

P/N:36025



Rubber Foot, Each

422 Series

OEM:26566011

P/N:36026



Thyister

422 Series

OEM:26896007

P/N:36030



Temp. Sensor

422 Series

OEM:86001274

P/N:36053



Lid Seal

422 Series

OEM:25437049

P/N:36005



Motor Boot

422 Series

OEM:25439072

P/N:36008




Touch Pad

422 Series

OEM:25176008

P/N:36072



Rotor Nut

422 Series

OEM:25292004

P/N:36004



Door Switch

Lever Style

OEM:3171945

P/N:36009



Lid Lock Grommet

422 Series

OEM:26554028

P/N:36024




Motor Mount

422 Series

OEM:26007026

P/N:36012



Lid Cap

422 Series

OEM:3174805

P/N:36013



Power PCB Exchange

422 Series

OEM:26005001

P/N:36011

Thermo / Jouan

422 Series



Motor Rebuild Only
For K422 Series
OEM:25509607
P/N:OCM-09F



Brush Pair
OEM:26160032
P/N:CB-099

4232



Motor Exchange
Rebuild Price Available!
OEM:89000511
P/N:OCM-09R



Brush Pair
OEM:86000637
P/N:CB-133

B412/4235



Motor Exchange
Rebuild Price Available!
OEM:11201499
P/N:OCM-09P



Brush Pair
OEM:86000056
P/N:CB-133

C411/C300



Motor Exchange
Rebuild Price Available!
OEM:25509545
P/N:OCM-09A



Brush Pair
OEM:86000637
P/N:CB-133

C4I



Motor Exchange
Rebuild Price Available!
OEM:89002888
P/N:OCM-09FF

4237



Motor Exchange
Rebuild Price Available!
OEM:11201527
P/N:OCM-09EE

B311 / BR311



Motor Exchange
Rebuild Price Available!
120V B/BR311
OEM:25509608
P/N:OCM-09B



Brush Pair
Before 06/1988 .2515"
OEM:25160027
P/N:CB-063A



Brush Pair
After 06/1988 .2300"
OEM:26160031
P/N:CB-097

Thermo / Jouan

B4I



Motor Exchange
Rebuild Price Available!
OEM:89000878
P/N:OCM-09BB



New Main PCB
OEM:89002923
P/N:36088



Gas Shock
OEM:26962012
P/N:36034

C3I



New Motor Exchange
OEM:11250202
P/N:30154



Motor Exchange
Rebuild Price Available!
OEM:11250202
P/N:OCM-09AA



Gas Shock
OEM:26962012
P/N:36034



New Main PCB
OEM:89002923
P/N:36088



Lock Assy.
OEM:89000848
P/N:36064

4217 / 4218



New Motor Exchange
OEM:89000629
P/N:30138



Motor Exchange
Rebuild Price Available!
OEM:89000629
P/N:OCM-09T



Brush Pair
OEM:86000637
P/N:CB-133

MR18.12



Motor Exchange
Rebuild Price Available!
OEM:89000064
P/N:OCM-09Z-1



Brush Pair
OEM:26160031
P/N:CB-097



Temp. Sensor
312 / 412 Series
OEM:86001274
P/N:36053

Sanyo

Mistral 1000 / R



Motor Exchange
Rebuild Price Available!
OEM:26091006
P/N:OCM-15



Brush Pair
Soldering Required
OEM:61600
P/N:CB-105



Brush Pair
No Soldering Required
OEM:NA
P/N:CB-105A

Mistral 2000 / R



Motor Exchange
Rebuild Price Available!
Mistral 2000
OEM:66500-816
P/N:OCM-15A



Motor Exchange
Rebuild Price Available!
Mistral 2000R
OEM:66500-816
P/N:OCM-15A-1



Brush Pair
Non-Autostop
OEM:34500-015
P/N:CB-113



Brush Pair
Auto-Stop Brush
OEM:NA
P/N:CB-113A



Lid Screw, Each
OEM:34500-015
P/N:40001



Rotor Bolt
OEM:34500-015
P/N:40002



Fuse Holder
OEM:34500-015
P/N:40005



Lock Microswitch
OEM:34500-015
P/N:40006



Solenoid
OEM:34500-015
P/N:40011



Tension Spring
OEM:78100-087
P/N:40013



Imbalance Switch
OEM:61200-197
P/N:40016



Triac Kit
OEM:61200-262T
P/N:40018



Lid Hinge
OEM:81900-474U
P/N:40037



Tension Spring
OEM:78100-087
P/N:40031



Left Torsion Spring
OEM:78100-112
P/N:40035

Sanyo

Mistral 2000 / R



Right Hinge Stop
OEM:81900-475U
P/N:40038



Left Hinge Stop
OEM:81900-476U
P/N:40039



Plunger
OEM:86120-094
P/N:40040



Lid Hook
OEM:87900-031U
P/N:40041



SCR Kit
OEM:NA
P/N:40043



Plunger Head
OEM:72020-208
P/N:40044

Mistral 3000E / I



Motor Exchange
Rebuild Price Available!
50 hz motor
OEM:66500-862
P/N:OCM-15B



Motor Exchange
60 hz motor
OEM:66500-862
P/N:OCM-15B-1



CPU Board Exchange
For Mistral 3000E
OEM:61080-187R
P/N:40015



CPU Board Exchange
For Mistral 3000I
OEM:61080-186R
P/N:40022



Display PCB Exchange
For Mistral 3000 I
OEM:61080178R
P/N:40014



Fan
OEM:32100-020
P/N:10090



Inverter Assy. Exchange
OEM:43060-004R
P/N:40017



Microswitch
OEM:36020-027
P/N:40007



Gaiter
OEM:71500-039
P/N:40023



Drive Belt
OEM:74030-002T
P/N:40029



Compression Spring
OEM:78100-050
P/N:40033

Savant

HSC 10K



Motor Exchange

Rebuild Price Available!

OEM:NA

P/N:OCM-26



Brush Pair

OEM:MOO-0010-01

P/N:CB-165

Separation Technology

Stat 60



Motor Exchange Motor Mount, Each

Rebuild Price Available!

For Stat 60

OEM:250-100

P/N:OCM-27



For Stat 60

OEM:730-101

P/N:730-101

Thermo Shandon

Cytospin Series



Motor Exchange

Rebuild Price Available!

For Cytospin 2 / 3

OEM:P07402

P/N:OCM-22A



Transistor

OEM:P09407

P/N:38481

Sorvall

Capspin / Microspin

Econospin



Motor Exchange
Rebuild Price Available!
OEM:2045
P/N:OCM-10S



Motor Exchange
Rebuild Price Available!
OEM:E850
P/N:OCM-10T



Brush Pair
OEM:H0276
P/N:CB-090



Motor Exchange
Rebuild Price Available!
OEM:H-1232
P/N:OCM-10G



Brush Pair
OEM:H0334
P/N:CB-091

Evolution

Microspin



Motor & Gyro
Sold Exchange
Rebuild Price Available!
OEM:72840
P/N:30155



Motor Exchange
Rebuild Price Available!
OEM:73448
P/N:30157



Motor Exchange
Rebuild Price Available!
OEM:E850
P/N:OCM-10T



Brush Pair
OEM:H0276
P/N:CB-090

Omnispin

RMC-14



Motor Exchange
Rebuild Price Available!
OEM:H-6000
P/N:OCM-10Z



Brush Pair
OEM:01785OR
P/N:CB-145



Motor Exchange
Rebuild Price Available!
OEM:12802
P/N:OCM-10CC



Brush Pair
OEM:12980
P/N:CB-164

MC12C / MC12V



Motor Exchange
Rebuild Price Available!
OEM:H-6000
P/N:OCM-10V



Motor Exchange
Rebuild Price Available!
OEM:NA
P/N:OCM-10V-1



Brush Pair
OEM:12934
P/N:CB-156

Sorvall

RC-2 / SS-3 / SS-4



New Motor Exchange
For RC-2, SS-4, SS-3
OEM:12353
P/N:30142



Motor Exchange
For RC-2, 6
prong connector
Rebuild Price Available!
OEM:12353
P/N:OCM-10L



Motor Exchange
For SS-3
Rebuild Price Available!
OEM:25245
P/N:OCM-10X



Brush Pair
OEM:61809
P/N:CB-032

RC-2B



Motor Exchange
Rebuild Price Available!
OEM:50222
P/N:OCM-10U



Brush Pair
OEM:12284
P/N:CB-036



Spline Drive Exchange
OEM:20466R
P/N:20510



New Speed Control
Sold Exchange
OEM:12217
P/N:38498

RC-3



Motor Rebuild Only
OEM:12221
P/N:OCM-10J



Brush Pair
OEM:21613
P/N:CB-019



Gyro Exchange
OEM:12710
P/N:20503



Rotor Adapter
OEM:44091
P/N:38581

RC-3B / 3B Plus



Motor Exchange
Rebuild Price Available!
OEM:55344
P/N:OCM-10K



Motor Exchange
Rebuild Price Available!
OEM:55242
P/N:OCM-10K-1



Brush Pair
OEM:55176
P/N:CB-042



Brush Pair
Included Brush Caps
OEM:NA
P/N:CB-042A



Gyro Exchange
OEM:55241R
P/N:20518



Lower Coupling
For Spline Drive
OEM:20462
P/N:38055



Bottom Shaft
For Gyro
OEM:NA
P/N:20521



Power PCB Exchange
OEM:55241R
P/N:38063

Sorvall

RC-3B / 3B Plus



Gas Lid Stay
OEM:66388
P/N:38069



Rotor Nut
OEM:20462
P/N:38044



Panel Light
OEM:NA
P/N:38053



Indicator PCB Exchange
OEM:55089
P/N:38065



Magnetic Pickup
OEM:55129
P/N:38066



Temp Sensor RTD
OEM:20462
P/N:38055



Modulator PC Kit
OEM:55304
P/N:38074



Gyro Insulation
OEM:55213
P/N:38073



Leveler
OEM:50227
P/N:38079



Microswitch
OEM:62716T
P/N:38394



Power Heat Sink Exchange
OEM:55186R
P/N:38502



Panel Indicator Lamp
OEM:NA
P/N:38539



Timer, 30 Min.
OEM:66388
P/N:38039



HLR-6 Rotor & Buckets
OEM:11350
P/N:38574

RC-3B Plus Induction Motor



Motor Rebuild
Must call for RGA #
OEM:56430
P/N:OCM-10EE



Optical Pickup
OEM:55635
P/N:38077



Gyro Boot
OEM:74196
P/N:38041



HLR-6 Rotor & Buckets
OEM:11350
P/N:38574

Sorvall

RC-3C / 3C Plus



Motor Rebuild
Must call for RGA #
OEM:55242
P/N:OCM-10K-1



Brush Pair
OEM:55176
P/N:CB-042



Brush Pair
Includes Brush Caps
OEM:NA
P/N:CB-042A



Optical Pickup
OEM:55635
P/N:38077



Gyro Exchange
OEM:55241R
P/N:20518



Lower Coupling
For Spline Drive
OEM:20462
P/N:38055



Optical Pickup Disk
OEM:50735
P/N:38062



Imb. Switch
OEM:55201T
P/N:38072



Swivel Caster
OEM:62086
P/N:38109



Gas Lid Stay
OEM:66388
P/N:38069



HLR-6 Rotor & Buckets
OEM:11350
P/N:38574

RC-5B



New Motor Exchange
OEM:50222
P/N:30139



New Motor Exchange
OEM:20210326
P/N:30143



Motor Exchange
Rebuild Price Available!
Spline Drive
OEM:50222
P/N:OCM-10P



Motor Exchange
Rebuild Price Available!
Quiet Drive
OEM:20210326
P/N:OCM-10P-1

Sorvall

RC-5B



Brush Pair
.250" or 6.3mm Thick
OEM:12284
P/N:CB-036



Lower Coupling Assy.
OEM:12345
P/N:38047



Rubber Insert
OEM:67262
P/N:20520



Speed Nut
OEM:66189T
P/N:38087



Lower Coupling
For Spline Drive
OEM:20462
P/N:38055



Flex Bolt Mount
OEM:50408
P/N:38060



Motor Baffle
OEM:55BMB1
P/N:38068



Rubber Stem Bushing
OEM:60669
P/N:38110



Spline Drive Exchange
OEM:20466R
P/N:20510



Quiet Drive Exchange
OEM:12817R
P/N:20507

RC-5B Plus



Motor Exchange
Rebuild Price Available!
OEM:74510
P/N:OCM-10R-1



Brush Pair
OEM:12284
P/N:CB-036



Brush Pair
.375" or 9.5mm Thick
OEM:74196
P/N:CB-036-1

RC-5C



Motor Exchange
Rebuild Price Available!
OEM:12991
P/N:OCM-10R



Brush Pair
.250" or 6.3mm Thick
OEM:12284
P/N:CB-036



Brush Pair
.375" or 9.5mm Thick
OEM:74196
P/N:CB-036-1



Optical Pickup
OEM:55635
P/N:38077

Sorvall

RC-5C



Lower Coupling Assy.
OEM:12345
P/N:38047



Rubber Insert
OEM:67262
P/N:20520



Speed Nut
OEM:66189T
P/N:38087

RC-5C Plus



Motor Exchange
Rebuild Price Available!
OEM:74510
P/N:OCM-10R-1



Brush Pair
.250" or 6.3mm Thick
OEM:12284
P/N:CB-036



Brush Pair
.375" or 9.5mm Thick
OEM:74196
P/N:CB-036-1



Optical Pickup
OEM:55635
P/N:38077

Technospin



Motor Exchange
Rebuild Price Available!
OEM:H-6000
P/N:OCM-10N



Brush Pair
OEM:H5877
P/N:CB-087



Motor Exchange
Rebuild Price Available!
OEM:07405
P/N:OCM-10Y



Brush Pair
OEM:07487
P/N:CB-168



Drive Belt
OEM:878513T
P/N:10107

RT7

Sorvall

GLC Series

Sorvall



Motor Exchange
For GLC 1 / 2
Rebuild Price Available!
OEM:06175
P/N:OCM-10



Motor Exchange
For GLC 2B, GE Motor
Rebuild Price Available!
OEM:06369
P/N:OCM-10B



Motor Exchange
For GLC 2B / 3
Rebuild Price Available!
OEM:06175
P/N:OCM-10A



Motor Exchange
For GLC 2B / 3
Rebuild Price Available!
OEM:66301
P/N:OCM-10C




Motor Exchange
For GLC 4
Rebuild Price Available!
OEM:12735
P/N:OCM-10D



New Motor Exchange
For GLC 4
OEM:06020
P/N:30146



Brush Pair
For Grey or Black
Dumore Motor
OEM:12293
P/N:CB-016



Brush Pair
For Howard Motor
OEM:12279
P/N:CB-017



Lid Seal
For GLC 4
OEM:06543
P/N:38037




Timer, 30 Min.
OEM:07513
P/N:38039



Speed Control
For GLC 2 Series
OEM:12245
P/N:38045



PCB Exchange
For GLC 2B
OEM:12246N
P/N:38046



Relay
OEM:62542
P/N:38081



Brake Switch
For GLC 2B
OEM:66158
P/N:38085



Lower Coupling
For GLC 4
OEM:12245
P/N:38045



Knob
For GLC 2 / 3 Series
OEM:62127
P/N:38112



New Techmeter
For GLC 2
OEM:06164
P/N:38541



Hinge Pin
OEM:06186
P/N:38035

Sorvall

T / RT-6000



Motor Exchange
Rebuild Price Available!
OEM:07887
P/N:OCM-10E



Brush Pair
OEM:12293
P/N:CB-016



Drive Belt
OEM:878513T
P/N:10107



30 Min. Timer
OEM:07513
P/N:38039



Right Latch Exchange
OEM:07525R
P/N:20136



Left Latch Exchange
OEM:07526R
P/N:20137



Woodruff Key
OEM:64877
P/N:38113



Spindle Boot
OEM:07638
P/N:38041



Spindle Exchange
OEM:07889R
P/N:20512

Cellwasher II



Motor Exchange
GE Motor
Rebuild Price Available!
OEM:12975
P/N:OCM-10M



Motor Exchange
EAD Motor
Rebuild Price Available!
OEM:12975
P/N:OCM-10M-1



Front Display Exchange
OEM:04591R
P/N:38003



Latch Exchange
OEM:12792R
P/N:20131



Pump Head Assy.
OEM:67632
P/N:38070



DA-12 Exchange
OEM:04285U
P/N:38104



Collector Ring
OEM:12981
P/N:38052



Tubing Kit
OEM:12977T
P/N:38005

Sorvall

Cellwasher II

Sorvall




Latch Exchange
OEM:12792R
P/N:20131



Bowl Guard
OEM:04565
P/N:38001



Decant Coil
OEM:04214
P/N:38007



Drive Ring
OEM:04277T
P/N:38008



DA-12 Clip Kit
OEM:04330
P/N:38009



Rear Wet Guard
OEM:04517
P/N:38010



Pump Tubing
OEM:04214
P/N:38011



Retainer Ring
OEM:04538
P/N:38012



Nozzle Assy.
OEM:04552
P/N:38014



GE Capacitors
OEM:04568T
P/N:38015



Power Switch
OEM:04570
P/N:38017



Corecom Assy.
OEM:04574
P/N:38019



Latch Plate Magnet
OEM:67653
P/N:38097



Stainless Steel Hinge
OEM:04698T
P/N:38026



Front Panel Overlay
OEM:04371
P/N:38018



**Latch Plate
w/o magnet**
OEM:04864U
P/N:38028



Drain Tray Seal
OEM:04872
P/N:38029



Drain Tray
OEM:04874
P/N:38030



Sealing Gasket
OEM:04975
P/N:38033



Needle Value
OEM:12925T
P/N:38050



Collector Ring
OEM:12981
P/N:38052



Distributor
OEM:43260
P/N:38059




Spring Clip
OEM:43231
P/N:38058



Bottle Kit
For Drain Tray Overflow
OEM:JP-63526-40T
P/N:38071

Sorvall

Cellwasher II



Motor Mount
OEM:90983
P/N:38080



5V Regulator
OEM:67078
P/N:38090



Capacitor
5600uF, 35V
OEM:67628
P/N:38095



Transformer
OEM:67635
P/N:38096



Flow Switch Assy.
OEM:04620
P/N:38024



Rubber Foot
OEM:68568
P/N:38099



O-Ring for Retainer
OEM:69722
P/N:38100



Fuse Kit
6.25A Slo-Blo (Qty. 5)
OEM:NA
P/N:38101



Saline Shield
OEM:04941
P/N:38116



Latch Solenoid
OEM:NA
P/N:38118



Cooling Fan
OEM:93006
P/N:38486



Lid Assy. Exchange
OEM:04660
P/N:38521



Pump Motor
OEM:04578
P/N:38020



New Bowl Assy. Exchange
OEM:12796
P/N:38049



Bowl Assy. Exchange
OEM:12796R
P/N:38002

Ozark Biomedical is a ISO9001:2008 Registered Company.

**Quality Parts and Services, along with Customer Satisfaction,
are Our Top Priorities.**



Centrifuges, Motors & Replacement Parts

Brushes are a crucial element in any motor.

Proper size, grade of carbon and tension on the commutator are all important elements.

If the grade of carbon is too hard, premature failure of the commutator is likely.

If the grade of carbon is too soft you will be changing brushes frequently.

Spring tension is equally important for the same reasons.

Motor design and current draw determine if a brush should have a shunt wire or not.

Our Brushes Are Manufactured From OEM Specifications!

Manufacturers change part numbers periodically.

Our catalog is constantly evolving and has become a leading source of gathered information.

Please consult the service manual and original manufacturer if any questions arise concerning brush application. Always feel free to call us with any questions.

Our Brushes are sold in pairs! All pricing in this catalog are for brush pairs!

Special Brush Discounts:

10% off orders over \$100.00

15% off orders over \$200.00

20% off orders over \$300.00

25% off orders over \$500.00

We carry Brush Holders and Brush Caps for most centrifuge motors!

For Brushes Not Listed Please Inquire!

We Will Gladly Add Brushes We Do Not Have to Our Stock!

Tech Tip!

Did you know that high leakage current in many centrifuges is due to a buildup of carbon dust in and around the motor! Blow those motors out when you change brushes.

Vist Us at www.ozarkbiomedical.com

Brush Application List

Brush Prices Are Per Pair Unless Otherwise Noted

<u>Model</u>	<u>Notes</u>	<u>OEM Part #</u>	<u>Ozark Part #</u>
Abbott			
TDX(LN9527-01)	MSI&Carter Motor		CB-003
	Baldor Motor		CB-304
	Bodine Motor		CB-026
Baxter			
Biofuge A	Model 1302	70001786	CB-050
Biofuge B		70001195	CB-051
Biofuge 17/17R		70002779	CB-101
Cryofuge 6000,			
Cryofuge 8500		70006914	CB-103
Cryofuge 8000		70006698	CB-158
Haemofuge	Model 1216		CB-117
Labofuge A	Model 2502	70001754	CB-052
Labofuge B	Model 2650	70002228	CB-053
Medifuge	Model 1215	70001193	CB-049
Megafuge	Model 2630	70002790	CB-094
Omnifuge		70003854	CB-102
Varifuge RF		70001811	CB-093
Beckman			
Accuspin	Carter Motor	349869	CB-044
	Elinco Motor	349544	CB-043
GP/GS-6/Spinchron Allegra 6/R	Solid Type Brushes Only.	356527	CB-081
GP/GS-6/Spinchron Allegra 6/R	Kit-Includes caps, grommets, locknuts & 1 pair brushes.		CB-081A
GP/GS-6/Spinchron Allegra 6/R	Autostop Brushes Only.	363442	CB-155
GP/GS-6/Spinchron Allegra 6/R	Kit-includes caps, grommets, locknuts & 2 pair brushes.	363443	CB-155A
J Drives	Includes one with Indicator wire.	349456	CB-041
	Longer thin brush for Beckman.	349454	CB-144
	Rebuilt motor #329417		
L340,L350,L540	2-Section Brush	333384	CB-033
Microfuge 11 & 12		343662	CB-046
Microfuge B	Dayton Motor	343156	CB-046A
Microfuge B	Howard Motor	343157	CB-046B
Microfuge E	Short Brush	356883	CB-096
Microfuge E	Long Brush	356881	CB-100
Microfuge Lite		365635	CB-173
TJ-6/6R	Carter Motor	340536	CB-020
TJ-6/6R	Elinco Motor	349544	CB-021
Microfuge 152			CB-084
Bio-Dynamics			
630 Select-A-Fuge	Ser. #70999 & Below	919807	CB-026
	Ser. #71000 & Up	919808	CB-027
640 Basic	Model 2008	106403	CB-027B
646 Hematocrit	Model 2011	106474	CB-027C
Hematocrit	Hettich Motor		CB-108
Universal #24	Model 1200	106433	CB-27A

Brush Application List

Brush Prices Are Per Pair Unless Otherwise Noted

<u>Model</u>	<u>Notes</u>	<u>OEM Part #</u>	<u>Ozark Part #</u>
--------------	--------------	-------------------	---------------------

Brinkmann / Eppendorf

3200	Silver Motor	2236520-7	CB-054
	Black Motor	2236521-5	CB-055
5412,5413,5414		2236521-5	CB-075
5415	Ser. # 3716 & Below	2236524-0	CB-076
	Ser. # 3717 & Up	2236523-1	CB-077
5415,5402	Replaces CB-077	2236527-4	CB-078

Clay Adams (Becton Dickinson)

Autocrit, MCHT	Black Universal Motor	42055102	CB-001
Autocrit II, MCHT II	Black Howard Motor	42055802	CB-001A
Autocrit II, MCHT II	RPI Replacement Motor	IEB024	CB-001C
Autocrit Ultra-3	Specialty Motor	42055813	CB-001B
	Jones Motor	42055813	CB-001B
Dynac, Dynac II	Jones Motor	42000102	CB-002
Physicians Jr.		0121600001	CB-004
QBC	Ser. # 184999 & Below	42445709	CB-083
	Ser. # 185000 & Up	42445713	CB-085
	Governor Brush	42445710	CB-131
Micro 2000		86000637	CB-207
Readacrit/MP	Blk or Silver Mtr small caps	42059102	CB-005
	Silver Motor W/large caps	42059126	CB-089
Safeguard, Safety Head	Jones Motor	42000102	CB-002
	Short Grey Motor		CB-002A
Serofuge II	Governor Brush	42054120	CB-070
Triacs		42020002	CB-003

Cobe

C-3 Centrifuge System	Blood Pump	500128200	CB-149
-----------------------	------------	-----------	--------

Dade

C7-M		014127-001	CB-047
Dac I	Magnet Brush	014533-100	CB-048
	Motor Brush	015532	CB-048A
Dac II	Magnet Brush	S45099-001	CB-147
	Motor Brush	S45098-001	CB-148

Drucker

L708, L708E, 411, 720E		7735023	CB-018
803, 806, 809, 809T		7735003	CB-030
809R, 818R			
896C2, MPC2001R		7735013	CB-031
MPC1500, MPC3000			
708E			CB-002

Brush Application List

Brush Prices Are Per Pair Unless Otherwise Noted

<u>Model</u>	<u>Notes</u>	<u>OEM Part #</u>	<u>Ozark Part #</u>
Fisher			
59A		SPN06649	CB-038
Centrific 225/50		SPN39329	CB-106
235A, 235B, 235C, 335		43425	CB-039
6K		SPN83788	CB-121
10K, 10KR		HL320140	CB-123
13K		SPN83142	CB-119
13KM		SPN83094	CB-120
21K, 21KR, 21KBR		HL320177	CB-122
Micro A, Micro H		SPN43425	CB-039A
Fresenius Dialysis			
Deareation Pump	Hard Brush/Soldering Required	670104	CB-201
Flow Pump/Blood Pump	Soft Brush/Soldering Required	643762	CB-202
Hamilton Bell			
Jr. Conical(1600)			CB-082
Hawksley			
Micro-Centrifuge		01980	CB-095
Hermle / National Labnet			
Z206A		Co200-B	CB-208
Z229, Z230M, Z230MA		Co230-G	CB-109
Z230MR, Z230H, Z230HA			
Z231M, Z320			
Z230AMK, Z232MK2	verify model# on back of unit Z230AMK2 says Z230A on front of unit but takes this brush	C2300-B	CB-166
Z230, Z232, Z230A		Co230-15B	CB-110
Z230MR-2, Z180M	Z180M with Silver Motor	C1236-20B2	CB-170
Z180M	Z180M w/ Black Carter Motor	C1236-20B	CB-171
Z360K, Z380, ZK380		Co360-G	CB-111
Z364, ZK364		Co364-G	CB-104
Z510, ZK510		Co510-G	CB-112
Hettich			
Rotanta		H5877	CB-140
Rotixa		H5780	CB-145
Universal		H0334	CB-141
EBA-20			CB-203

Brush Application List

Brush Prices Are Per Pair Unless Otherwise Noted

<u>Model</u>	<u>Notes</u>	<u>OEM Part #</u>	<u>Ozark Part #</u>
IEC			
2V, 2R, BPV, BPR		1724	CB-008
B-2, B-20, B35		2467	CB-010
B-20A	IEC Discontinued, 2 section type	44908	CB-012
B-20A	New Style, solid type	48264	CB-013
B-22		41775	CB-080
B-60	IEC Discontinued	42560	CB-034
C-Series	IEC Discontinued	1707	CB-024
CL, CH	Use P/N CB-006A	1780	CB-006
CL, CH, Centra CL2	New Style W/Shunt Wire	1780A	CB-006A
CM, C50		1731	CB-014
Centra-CLD		75202	CB-150
Centra-4		49203	CB-079
Centra-4B	Program 1.0 Motor W/White Caps	49374	CB-071
Centra-4B, CL4	Program 2.0 Motor W/Blk Caps	44989	CB-073
Centra-7	1/4x7/16x3/4 Old 2 section type	47263	CB-029
Centra-7	1/4x3/8x3/4	47889	CB-045
Centra-7	3/16x3/8x7/8 Most Units	48192	CB-035
Centra-8/8R, GP8/8R	IEC Motor	49801	CB-072
Centra GP8/8R	Baldor Motor	T01206S	CB-304
CU/CRU-5000	IEC Discontinued 2 section type	44908	CB-012
CU/CRU-5000	New Style Solid Type	48264	CB-013
C/CR-6000		41775	CB-080
H, HN-S, HN-SII	Use P/N CB-006A	1780	CB-006
HN-SII	New Style With Shunt Wire	1780A	CB-006A
HT		1860	CB-028 Discontinued
Equafuge		50757	CB-163
MB Microhematocrit	Specialty Motor	42055813	CB-001B
MB Microhematocrit	Ser. # up to AA2914		
	And Above AC6200	1840	CB-007
MB Microhematocrit	Jones Motor SerAA2914-AC6200	43653	CB-007A
MH			CB-006A
M-60	Round Motor	41822	CB-142 Discontinued
M-60	Square Motor	49846	CB-143
MP4, MP4R	Program 1.0 Motor W/White Caps	49374	CB-071
	Serial #s 24370600 & Below.		
	Refrig Serial #s 24380600 & Below		
MP4, MP4R	Program 2.0 Motor W/Blk Caps	44989	CB-073
	Serial #s 24371999 & Below.		
	Refrig Serial #s 24381999 & Below		
MP4, MP4R	Inverted Motor	50433	CB-162
	Serial #s 24371999 & Up.		
	Refrig Serial #s 24381999 & Up		
PR-1, PR-2	IEC Discontinued	1724H	CB-009
PR-6, PR-J		40069	CB-011
PR/DPR 6000	IEC Discontinued Old Style	44908	CB-012

Brush Application List

Brush Prices Are Per Pair Unless Otherwise Noted

<u>Model</u>	<u>Notes</u>	<u>OEM Part #</u>	<u>Ozark Part #</u>
IEC			
PR/DPR 6000	New Style Solid Type	48264	CB-013
PR-7000		41775	CB-080
SBV, SBR, U	(Hand Brake)	1708	CB-025
Size 2/ Model K	Autotransformer Brush	1785	CB-146
Size 2/Model K, U, KB, EXD		1913	CB-023
UV	Autotransformer Brush	1927	CB-074
UV,CS		1946	CB-015
JOUAN			
4203,4204 Micro Fugette	Use 86000637	3860	CB-132
A12, A13, A14, ALC 4206, 4214,4218 , 4235, 2403/04		86000637	CB-133
BBVV		26160045	CB-063B
B/BR311, V511	Before 06/1988 .2515"	25160027	CB-063A
B/BR311	After 06/1988 .2360"	26160031	CB-097
C/CR312, C/CR412, 422	Up-graded ALC Motor	89000214	CB-134
G/GR412/422	(6mmx16mmx32mm)		
C/CR312, C/CR412, 422	Old Style Motor	89000056	CB-135
G/GR412/422	(6mmx10mmx37mm)		
C/CR-3000, C/CR/G/		89000056	CB-135
GR/GT411			
Hemac		26160058	CB-137 Discontinued
Hemac		29993144	CB-138 Discontinued
K110SX, LR522		26160056	CB-129
KR422		26160032	CB-099
M/MR 1411	Before 06/1988	26160029	CB-098
M/MR 1411	After 06/1988	26160031	CB-097
M/MR 1822			CB-097
SR12/ 2022		26160030	CB-139 Discontinued
Lab Line			
Max Q 5000 Shaker		370-356-00	CB-205
3595 Shaker	W/Indiana General Motor	370-156-00	CB-151
3595 Shaker	W/Applied Motor	370-156-01	CB-152
3595, 3525, 3526, 3529	W/Bodine Motor	370-258-00	CB-153
3590, 3597			
Lourdes			
20G. R-30 Clinifuge	Floor Model	1811	CB-040
MLA			
Tytem	200R	70005	CB-154
Mistral/ MSE/Sanyo			
1000	Soldering Required	61600	CB-105
	If motor has brush caps		
	Order CB-113		
1000	New Style no soldering		CB-105A
2000		34500-015	CB-113
2000(Auto-Stop)			CB-113A
GF4, GT2, GT4			CB-088
Micro-Centar		34500-014	CB-161

Brush Application List

Brush Prices Are Per Pair Unless Otherwise Noted

<u>Model</u>	<u>Notes</u>	<u>OEM Part #</u>	<u>Ozark Part #</u>
Precision			
Universal	Serial #s Starting W/10	223245	CB-056
Universal	Serial #s Starting W/11	223504	CB-057
Hi-Speed	Cat# 67310 & Up		
	Cat# 67390	25A18A-2	CB-058
Rupp Bowman/Global Focus			
H103NA, H103NB		00030	CB-124
H103NARD			
H103FRS		00050	CB-126
H103FRL		00050	CB-172
H2O		00029	CB-125
S103NA		00032	CB-127
S103ND		00058	CB-128
2100, 2200/R, 2300/R		00033	CB-159
Sakura			
Cytotec	Model 4325,4332	A0-93-0010	CB-160
Savant			
HSC10K/KA/HSC15R		MOO-0010-01	CB-165
SFR11K/SFR13K		137-1018-00	CB-169
Shandon			
Cytospin 1		PO1455-67	CB-067
Cytospin 2 & 3	Soldering Required	PO9351	CB-068
Sorvall			
Angle Head			CB-107 Discontinued
Capspin, Easyspin, Microspin		H0276	CB-090
Econospin		H0334	CB-091
GLC, GLC2,	Howard Motor	12279	CB-017
GLC2B, GLC3			
GLC, GLC2, GLC2B,	Grey Or Black Dumore Motor	12293	CB-016
GLC3, GLC4			
MC-12V, MC-12C		12934	CB-156
Omnispin		H5780	CB-145
RC-2		61809	CB-032
RC-2B, 5,5B, 5C	.250" or 6.3mm Thick	12284	CB-036
RC-5C(New Style Motor)	.375" or 9.5mm Thick	74196	CB-036-1
RC-3	Serial #s 70400 & Above	21613	CB-019
RC-3B, 3C Plus	2 Section With Indicator Wire	55176	CB-042
	CB-042 with Brush Caps		CB-042A
RMC-14	Refrigerated Table Top	12980	CB-164
SS-3, SS-4		61809	CB-032
T-6000, T-6000D,		12293	CB-016
RT-6000			
Technospin		H5877	CB-087
RT-7	Includes Brush Caps	07487	CB-168

Brush Application List

Brush Prices Are Per Pair Unless Otherwise Noted

<u>Model</u>	<u>Notes</u>	<u>OEM Part #</u>	<u>Ozark Part #</u>
Stryker			
840 Series Cast Cutter		870-1-13	CB-064
810 Autopsy Saw			
Electo-surgical Unit	Model 1301-115		CB-157
Vulcan			
Microspin 24	With Rae Motor		CB-114
Microspin 24	With Redman Motor		CB-118
Quantum,			CB-204
Microspin 24	With MET Motor *These fit holders 5/16" x 5/16", if holders measure 1/4" x 3/8" use CB-018*		



Still don't see a brush for your Centrifuge?
Call 800-457-7576 to speak with an Application Specialist that
can help locate a brush meeting your
specifications!

**Centrifuge problems wearing you out?
Help is a phone call away!**

800-457-7576



**Picture by Tim Ernst: “Black Bear, Ozark National Forest, Arkansas”
Visit his website at cloudland.net**



As an alternative to purchasing new units, we offer top quality reconditioned centrifuges.

Our reconditioned centrifuges look and operate like new, at a fraction of the price.

Some of the many reasons to consider a reconditioned centrifuge

1. Your Customer is important to us too!
2. Significant Savings--thoroughly reconditioned with warranty, providing maximum dollar value.
3. We will only provide instruments that have a reliable track history.
4. We provide service and support after the sale.
5. We will not sell you someone else's problems.
6. We take centrifuges in on trade.
7. We can solve your centrifuge problems and get parts on site in a timely fashion.
8. We help minimize down time.
9. We specialize in centrifuges, and have for over 25 years.

This is what you deserve in centrifuge sales and service!

Reconditioning Process

Our reconditioning process addresses all major components and subassemblies of the centrifuge, including the motor and latch assemblies. Touch panel overlays are replaced. Rotors are inspected and reconditioned.

Any mechanical timers and switches are replaced. Motor mounts are always replaced. These components are replaced with quality new and reconditioned components from our stock. Because we sell parts to customers worldwide you can be assured of top quality.

The centrifuge is then inspected for any cosmetic defects. Paint is matched to original manufacturers color by our in-house paint technicians who restore the units to "like new" appearance. The finished unit passes on to our quality assurance technicians. A complete alignment of the electronics is performed assuring proper acceleration, deceleration, speed detection, timing and imbalance detection.

Finally, our refurbished units are shipped utilizing a combination of our Sealed Air Instapak foam, custom built pallets and heavy duty cardboard boxes.

We ship any centrifuge over 100 pounds via trucklines individually palletized and packaged to assure safe arrival at your site.

At Ozark Biomedical, we are experts in assuring your unit gets to you in perfect condition.



To Place An Order:

Call 1-800-457-7576 or Fax orders to 501-882-3986

Order Online at www.ozarkbiomedical.com

We have no minimum orders.

Orders placed prior to 4:00 PM Central Standard Time will be shipped the same day!

Terms:

Net terms available following credit approval or COD. We accept Visa, Mastercard and American Express.

Past due invoices are subject to a late charge of 1.5% per month

No allowance for shortages, overcharge or damage if not reported within 10 days.

All returned parts are subject to a 30% restocking fee

Shipping and handling is prepaid and added to the invoice. Minimum charge is \$7.00

Parts sold on exchange basis not returned within 15 days or which are deemed non-repairable are subject to a core charge. The core charge is the difference between our price and full list price.

Your help returning these items is greatly appreciated.

This catalog supersedes all previous editions.

All product names used in this document are trademarks or registered trademarks of their respective holders.

Special Discounts

10% off brush orders over \$100.00

15% off brush orders over \$200.00

20% off brush orders over \$300.00

25% off brush orders over \$500.00

Brush prices are per pair!

Please read Carefully:

Ozark Biomedical's liability is limited to the cost of the defective part. The buyer assumes all risk and liability resulting from the use of these products, whether individually or in combination with other products.

Prices subject to change without notice!



1001 Commerce Place Beebe, Arkansas 72012

Prsrt Std
U.S. Postage
PAID
Permit No. 9
Beebe, AR