

BD Clay Adams Autocrit Ultra 3 Microhematocrit Centrifuge



BD Clay Adams Autocrit Ultra 3 Microhematocrit Centrifuge, Model 420575, BD Diagnostics is a reliable and dependable addition to the [BD Centrifuges](#) family of products. Combining top-notch and uncompromising quality with an affordable price, the **BD Clay Adams Autocrit Ultra 3 Microhematocrit Centrifuge, Model 420575, BD Diagnostics 420575/ 15176-419** can fulfill your laboratory needs while still offering a great value for the money.

Product Category Description:

Designed for hematology applications. Features direct-reading spiral scale built into the head. Makes specimen handling and reading more efficient.

Electronic reset timer automatically returns to selected interval after each cycle. Automatic locking latch prevents lid from opening when power is on. Automatic braking at the end of each cycle prevents abrupt stops for a gradual, consistent deceleration. Transparent carrying tray holds capillary tubes from sample collection to readout. Shock absorbance system minimizes vibration.

Description:	CENTRIFUGE AUTOCRIT III
Maximum Speed:	11,700 rpm
Maximum RCF:	13,700 g
Capacity:	24 x 75 mm capillary tubes
Operating Temperature:	Ambient
Timer:	15 min.
Dimensions:	23H x 25W x 31D cm
Weight:	10.2 kg (22 lbs.)
Electrical:	120V, 60Hz, 216W