

Vibration Mount Kit Allegra 6

The original mounts for this unit have been discontinued by OEM.



The replacement kit from OEM includes a new motor mounting plate and mounts. It is inferior in Ozarks opinion due to premature failures during gross imbalances. Also technician will have to completely remove and disassemble mounting plate to install new one. (If a customer already has the new motor mounting plate and mounts and wants to continue using the metal stand-offs, P/N 10153 can be sold as direct replacement parts.)



P/N 367899

Our kit is the same as the original Spinchron (GS6) motor mount kit with the exception it includes 3 additional spacers for the Allegra 6 mounting plate difference. Mounts can be installed one at a time without removing mounting plate. If later model plate and mount is found in unit these mounts will work without the spacer.



P/N 10158

Spacer